# ALP-LQ (Alkaline phosphatase)

p-Nitrophenylphosphate. kinetic. Liquid. DGKC


## Quantitative determination of alkaline phosphatase (ALP)

### IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Alkaline phosphatase (ALP) catalyses the hydrolysis of p-nitrophenyl phosphate at pH 10.4, liberating p-nitrophenol and phosphate, according to the following reaction:

The rate of p-nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of alkaline phosphatase present in the sample<sup>1,2</sup>.

### CLINICAL SIGNIFICANCE

Alkaline phosphatase is an enzyme present in almost all weaves of the organism, being particularly high in bone, liver, placenta, intestine and kidney. Both increases and decreases of plasma ALP are of importance clinically.

Causes of increased plasma ALP: Paget's disease of bone, obstructive liver disease, hepatitis, hepatotoxicity caused by drugs or osteomalacia.

Causes of decreased plasma ALP: Cretinism and vitamin C deficiency<sup>1,5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

<b>R 1</b> Buffer	Diethanolamine (DEA) pH 10.4 Magnesium chloride	1 mmol/L 0.5 mmol/L
<b>R 2</b> Substrate	p-Nitrophenylphosphate (pNPP)	10 mmol/L

### PREPARATION

All the reagents are ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm  $\geq$  1.50.

### ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

### SAMPLES

Serum or heparinized plasma<sup>1</sup>. Use unhemolyzed serum, separated from the clot as soon as possible. Stability: 3 days at 2-8°C.

### Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1.00	1.22	2.64
30°C	0.82	1.00	1.33
37°C	0.61	0.75	1.00

### REFERENCE VALUES<sup>1</sup>

	25°C	30°C	37°C
Children (1-14 years)	< 400 U/L	< 480 U/L	< 645 U/L
Adults	60 - 170 U/L	73 - 207 U/L	98 - 279 U/L

Factors affecting ALP activities in a normal population include exercise, periods of repaid growth in children and pregnancy.

These values are for orientation purpose; each laboratory should establish its own reference range.

## APPLICATION SPINTECH 240

Item Name ALP <u>DATA INFORMATION</u> Units U/L Decimals 0 <u>ANALYSIS</u> Type RATE W.Length 1 405 W.Length 1 700 Method DGKC <u>CORR</u> SLOPE INTER 1.000 x + 0		<u>CALIBRATION</u> TYPE Linear  STANDARD #1 * #4 #2 #5 #3 #6 <u>NORMAL RANGE (37°C)</u> SERUM MALE LOW HIGH FEMALE URINE	
Item Name ALP <u>ASPIRATION</u> KIND Single <input checked="" type="checkbox"/> Double  VOLUME SAMPLE 5 µL REAGENT 1 240 µL REAGENT 2 60 µL  Third Mix <input checked="" type="checkbox"/> OFF ON R1 Blank <input checked="" type="checkbox"/> Water R1-B		<u>DATA PROCESS</u> <u>ABSORBANCE LIMIT</u> READ LOW -3.000 START END HIGH 3.000 MAIN 36 48 SUB  <u>ENDPOINT LIMIT</u> LINEAR CHECK (%) 90  <u>FACTOR</u> Blank Correction 1.000  <u>PROZONE CHECK</u> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	
<u>MONITOR</u> 0 LEVEL POINT 1 SPAN 3.000			

### QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of 0,6845 U/L to *linearity limit* of 1200 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

### Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	174	443	175	434
SD	0,72	1,56	6,88	11,93
CV (%)	0,41	0,35	3,93	2,75

**Sensitivity:** 1 U/L = 0,0003 ΔA/min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>:0,99938.

Regression equation: y = 1,025x - 1,105.

The results of the performance characteristics depend on the analyzer used.

### BIBLIOGRAPHY

1. Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
2. Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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- Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
- Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: TK41233	Cont.	R 1: 10 x 25 mL
		R 2: 10 x 7 mL

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