
**Quantitative determination of iron IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

The iron is dissociated from transferring-iron complex in weakly acid medium. Liberated iron is reduced into the bivalent form by means of ascorbic acid. Ferrous ions give with FerroZine a coloured complex:

 The intensity of the color formed is proportional to the iron concentration in the sample<sup>1,2</sup>.

**CLINICAL SIGNIFICANCE**

The iron is the component of a great number of enzymes. The myoglobin, muscular protein, contains iron, as well as the liver.

Iron is necessary for the hemoglobin production, molecule that transports oxygen inside red globules. Their deficit in the last causes the ferropenic anemia. High levels of iron are found in hemochromatosis, cirrhosis, hepatitis and in increased transferrin levels.

 The variation day to day is quite marked in healthy people<sup>1,5,6</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b>	Buffer	Acetate pH 4.9	100 mmol/L
<b>R 2</b>	Reductant	Ascorbic acid	99.7%
<b>R 3</b>	Color	FerroZine	40 mmol/L

**PREPARATION**

Dissolve (→) the contents of one tube of reductant (R2) in one bottle of R 1 Buffer. Cap and mix gently to dissolve contents. Stability: 3 months at 2-8°C or 1 month at 15-25°C. R3 is ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 570 nm ≥ 0.020.

**ADDITIONAL EQUIPMENT**

- Autoanalyzer Spintech 240.
- General laboratory equipment (Note 1).

**SAMPLES**

 Serum or heparinized plasma.  
 Free of hemolysis and separated from cells as rapidly as possible.  
 Stability of the sample: 2-8°C for 7 days<sup>1</sup>.

**REFERENCE VALUES<sup>5</sup>**

 Male 65 - 175 µg/dL ≅ 11.6 - 31.3 µmol/L (Note 4)  
 Female 40 - 150 µg/dL ≅ 7.16 - 26.85 µmol/L (Note 4)

These values are for orientation purpose; each laboratory should establish its own reference range.

**Conversion factor:** µg/dL x 0.179 = µmol/L.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**APPLICATION SPINTECH 240**

Item Name FE		CALIBRATION	
<b>DATA INFORMATION</b>		TYPE	Linear
Units	µg/dL		
Decimals	0		
<b>ANALYSIS</b>		STANDARD	
Type	END	#1	* #4
		#2	#5
		#3	#6
W.Length 1	570	<b>NORMAL RANGE</b>	
Method	Ferrozine	SERUM	MALE HIGH
			FEMALE
<b>CORR</b>			
SLOPE	INTER		
1.000 x +	0		
Item Name FE		<b>DATA PROCESS</b>	
<b>ASPIRATION</b>		<b>ABSORBANCE LIMIT</b>	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	LOW -3.000
		START END	HIGH 3.000
	<b>VOLUME</b>	MAIN	50 52
SAMPLE	25 µL	SUB	28 29
REAGENT 1	225 µL		
REAGENT 2	20 µL		
		<b>ENDPOINT LIMIT 3</b>	
		<b>LINEAR CHECK (%)</b>	
		<b>FACTOR</b>	
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	Blank Correction 1.000	
R1 Blank	Water <input checked="" type="checkbox"/> R1-B		
<b>MONITOR</b>		<b>PROZONE CHECK</b>	
0 LEVEL POINT	1	START END	LIMIT (%)
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

 Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than **40 days**.

**PERFORMANCE CHARACTERISTICS**
**Measuring range:** From detection limit of 0,850 µg/dL to linearity limit of 1000 µg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)	
Mean (µg/dL)	113	250
SD	0,89	0,72
CV (%)	0,79	0,29

	Inter-assay (n=20)	
	111	249
	3,51	6,29
	3,17	2,52

**Sensitivity:** 1 µg/dL = 0,00104 A.

**Accuracy:** Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents.

The results obtained using 50 samples were the following:

 Correlation coefficient (r)<sup>2</sup>: 0,9934.

Regression equation: y = 1,0243x - 3,877.

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. It is recommended to use disposable material. If glassware is used the material should be soaking for 6 h in diluted HCl (20% v/v) and then thoroughly rinsed with distilled water and dried before use.
2. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
3. Use clean disposable pipette tips for its dispensation.
4. Strongly method dependent.

**BIBLIOGRAPHY**

1. Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.
2. Itano M M D. Cap Serum Iron Survey 1978 (70): 516-522.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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**BIBLIOGRAPHIE**

- Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.
- Itano M M D. Cap Serum Iron Survey 1978 (70): 516-522.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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