Protein in urine and CSF

Pyrogallol red. Colorimetric

Quantitative determination of total urinary and CSF protein IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Proteins react in acid solution with pyrogallol red and molybdate to form a colored complex. The intensity of the color formed is proportional to the protein concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

In healthy persons, the urine contains no protein or only a trace amount of protein; normally the glomeruli prevent passage of protein from the blood to the glomerular filtrate. Glomerular injury causes increased permeability to plasma proteins, resulting in proteinuria, which refers to the presence of protein in the urine.

A persistent finding of proteinuria is the single most important indication of renal disease.

Elevated concentration of protein in cerebro-spinal fluid (CSF) can be caused by infections and intracranial pressure^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Pyrogallol red	50 µmol/L
	Sodium molybdate	0.04 mmol/L

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 600 nm \geq 0.30.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

- Urine 24 h: Stability 8 days at 2-8°C.
- Cerebrospinal fluid (CSF): Stable 4 days at 2-8°C.

REFERENCE VALUES⁵

Urine:	< 100 mg/24 h (< 150 mg/24 h in pregnancy)
CSF:	Children 300 -1000 mg/L
	Adults 150 - 450 mg/L

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

A calibration using a urine protein calibrator is needed (Ref.: 1002450). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

APPLICATION SPINTECH 240

Item Name PTOR <u>DATA INFORMATION</u> Units mg/L Decimals 0 <u>ANALYSIS</u> Type END W.Length 1 600 Method Colorimetric CORR SLOPE INTER 1.000 x + 0		<u>CALIBRATION</u> TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 <u>NORMAL RANGE</u> SERUM ADULT LOW HIGH CHILDREN URINE	
Item Name PTOR <u>ASPIRATION</u> KIND <input checked="" type="checkbox"/> Single Double VOLUME SAMPLE 5 µL REAGENT 1 250 µL Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B <u>MONITOR</u> 0 LEVEL POINT 1 SPAN 3.000		<u>DATA PROCESS</u> <u>READ</u> START END MAIN 50 52 SUB ENDPOINT LIMIT 3 LINEAR CHECK (%) <u>FACTOR</u> Blank Correction 1.000 <u>PROZONE CHECK</u> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit 9,44 mg/L to linearity limit of 4000 mg/L. If the concentration is greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n= 20)			Inter-assay (n= 20)		
	Mean (mg/L)	SD	CV (%)	216	499	1018
	220	3,7	2,28	18,3	26,1	166,1
	536	4,0	0,75	7,35	5,22	16,43
	1014	5,2	0,51			

Sensitivity: 1mg/L = 0,00026 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,9338

Regression equation: y = 0,4294x - 5,4159

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Orsonneau JL et al. An improved Pyrogallol Red-Molybdate Method for Determining Total Urinary Protein. Clin Chem 1989 (35):2233-2236.
2. Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1316-1324 and 418.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAFIA

1. Orsonneau JL et al. An improved Pyrogallol Red-Molybdate Method for Determining Total Urinary Protein. Clin Chem 1989; 35:2233-2236.
2. Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1316-1324 and 418.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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