**Quantitative determination of total protein IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Proteins give an intensive violet-blue complex with copper salts in an alkaline medium. Iodide is included as an antioxidant. The intensity of the color formed is proportional to the total protein concentration in the sample<sup>1,4</sup>.

**CLINICAL SIGNIFICANCE**

The proteins are macromolecular organic compounds, widely distributed in the organism. They act like structural and transport elements. The proteins of the serum are divided into two fractions, albumin and globulins. The determination of total proteins is useful in the detection of:

- High protein levels caused by hemoconcentration like in the dehydrations or increase in the concentration of specific proteins.
- Low protein level caused by hemodilution by an impaired synthesis or loss (as by hemorrhage) or excessive protein catabolism<sup>4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R</b>	Sodium potassium tartrate	15 mmol/L
	Sodium iodide	100 mmol/L
Biuret	Potassium iodide	5 mmol/L
	Copper (II) sulphate	19 mmol/L

**PRECAUTIONS**

H314-Causes severe skin burns and eye damage. H412-Harmful to aquatic life with long lasting effects.

Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

The reagent is ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 546 nm  $\geq$  0.22.

**ADDITIONAL EQUIPMENT**

- Autoanalyzer Spintech 240.
- General laboratory equipment.

**SAMPLES**

Serum or heparinized plasma<sup>1</sup>:

Stability of the sample: 1 month at refrigerator (2-8°C).

**REFERENCE VALUES<sup>1</sup>**

Adults: 6.6 – 8.3 g/dL

Newborn: 5.2 – 9.1 g/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**APPLICATION SPINTECH 240**

Item Name TP <u>DATA INFORMATION</u> Units g/dL Decimals 2 <u>ANALYSIS</u> Type END W.Length 1 546 Method Biuret <u>CORR</u> SLOPE INTER 1.000 x + 0		<u>CALIBRATION</u> TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 <u>NORMAL RANGE</u> LOW HIGH SERUM MALE FEMALE URINE	
Item Name TP <u>ASPIRATION</u> KIND <input checked="" type="checkbox"/> Single <input type="checkbox"/> Double VOLUME SAMPLE 7 $\mu$ L REAGENT 1 280 $\mu$ L Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B <u>MONITOR</u> 0 LEVEL POINT 1 SPAN 3.000		<u>DATA PROCESS</u> <u>READ</u> START END MAIN 50 52 SUB ENDPOINT LIMIT 3 LINEAR CHECK (%) <u>FACTOR</u> Blank Correction 1.000 <u>PROZONE CHECK</u> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 0,007 g/dL to linearity limit of 14 g/dL. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (g/dL)	SD	Mean (g/dL)	SD
Mean (g/dL)	6,53	4,89	6,77	5,08
SD	0,01	0,01	0,07	0,05
CV (%)	0,21	0,24	1,05	0,94

**Sensitivity:** 1 g/dL = 0,0825 A.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,97002

Regression equation y= 0,954x +0,511.

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

**BIBLIOGRAPHY**

1. Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1316-1324 and 418.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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**APPLICATION AU SPINTECH 240**

<p>Item Name TP</p> <p><b>DATA INFORMATION</b></p> <p>Units g/dL Decimals 2</p> <p><b>ANALYSIS</b></p> <p>Type END</p> <p>W.Length 1 546</p> <p>Method Biuret</p> <p><b>CORR</b></p> <p>SLOPE INTER 1.000 x + 0</p>	<p><b>CALIBRATION</b></p> <p>TYPE Linear</p> <p>STANDARD</p> <table style="width: 100%;"> <tr> <td>#1</td> <td>*</td> <td>#4</td> </tr> <tr> <td>#2</td> <td></td> <td>#5</td> </tr> <tr> <td>#3</td> <td></td> <td>#6</td> </tr> </table> <p><b>NORMAL RANGE</b></p> <table style="width: 100%;"> <tr> <td></td> <td>LOW</td> <td>HIGH</td> </tr> <tr> <td>SERUM</td> <td>MALE</td> <td>FEMALE</td> </tr> <tr> <td>URINE</td> <td></td> <td></td> </tr> </table>	#1	*	#4	#2		#5	#3		#6		LOW	HIGH	SERUM	MALE	FEMALE	URINE								
#1	*	#4																							
#2		#5																							
#3		#6																							
	LOW	HIGH																							
SERUM	MALE	FEMALE																							
URINE																									
<p>Item Name TP</p> <p><b>ASPIRATION</b></p> <p>KIND <input checked="" type="checkbox"/> Single <input type="checkbox"/> Double</p> <p>VOLUME</p> <p>SAMPLE 7 µL REAGENT 1 280 µL</p> <p>Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B</p> <p><b>MONITOR</b></p> <p>0 LEVEL POINT 1 SPAN 3.000</p>	<p><b>DATA PROCESS</b></p> <p><b>ABSORBANCE LIMIT</b></p> <table style="width: 100%;"> <tr> <td>READ</td> <td>LOW</td> <td>-3.000</td> </tr> <tr> <td>START</td> <td>END</td> <td>HIGH 3.000</td> </tr> <tr> <td>MAIN 50</td> <td>52</td> <td></td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table> <p>ENDPOINT LIMIT 3 LINEAR CHECK (%)</p> <p><b>FACTOR</b></p> <p>Blank Correction 1.000</p> <p><b>PROZONE CHECK</b></p> <table style="width: 100%;"> <tr> <td>START</td> <td>END</td> <td>LIMIT (%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td><input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High</td> </tr> <tr> <td>THIRD</td> <td></td> <td><input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High</td> </tr> </table>	READ	LOW	-3.000	START	END	HIGH 3.000	MAIN 50	52		SUB			START	END	LIMIT (%)	FIRST			SECOND		<input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High	THIRD		<input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High
READ	LOW	-3.000																							
START	END	HIGH 3.000																							
MAIN 50	52																								
SUB																									
START	END	LIMIT (%)																							
FIRST																									
SECOND		<input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High																							
THIRD		<input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High																							

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