# HDL Cholesterol D

Direct. Enzymatic colorimetric

## Quantitative determination of HDL cholesterol

### IVD



Store at 2-8°C

### PRINCIPLE OF THE METHOD


 Directly determination of serum HDLc (high-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation of the sample<sup>3,5</sup>

The assay takes place in two steps.

- 1° Elimination of lipoprotein no-HDL

- 2° Measurement of HDLc

The intensity of the color formed is proportional to the HDLc concentration in the sample.

### CLINICAL SIGNIFICANCE

HDL particles are high-density lipoproteins that transport cholesterol from the body tissues to the liver. Since HDL can remove cholesterol from the arteries and carry it back to the liver for their excretion, HDL is known as "good cholesterol" because high levels are thought to lower the risk of heart disease and coronary artery disease.

 A low HDL cholesterol levels, is considered a greater heart disease risk<sup>1,2,4</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

R 1	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 6.6	100 mM
	N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline (HDAOS)	0,7 mM
	Cholesterol Esterase	≥ 800 U/L
	Cholesterol oxidase	≥ 500 U/L
	Catalase	≥ 300 U/L
	Ascorbic oxidase	≥ 3000 U/L
R 2	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 7.0	1,1 mM
	4 - Aminoantipyrine (4-AA)	100 mM
	Peroxidase	≥ 3500 U/L

### PREPARATION

- R 1 and R 2: Are ready to use.

### STORAGE AND STABILITY (Note 1)

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze the reagents.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.

### ADDITIONAL EQUIPMENT

 - Autoanalyzer Spintech 240.  
 - General laboratory equipment.

### SAMPLES

Serum, heparinized plasma or EDTA plasma. If any samples show precipitates, centrifuge before using.

Stability of the sample: 6 days at 2-8°C and 1 year when stored at -70°C.

### REFERENCE VALUES

	Men	Women
Low risk	> 50 mg/dL	> 60 mg/dL
Normal risk	35 - 50 mg/dL	45 - 60 mg/dL
High risk	< 35 mg/dL	< 45 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

**Conversion factor:** mg/dL x 0,0259= mmol/L.

### QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## APPLICATION SPINTECH 240

Item Name HDL <u>DATA INFORMATION</u> Units mg/dL Decimals 1 <u>ANALYSIS</u> Type END W.Length 1 570 Method DIRECT <u>CORR</u> SLOPE INTER 1.000 x + 0		<u>CALIBRATION</u> TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 <u>NORMAL RANGE</u> LOW HIGH SERUM MALE FEMALE	
Item Name HDL <u>ASPIRATION</u> KIND Single <input checked="" type="checkbox"/> Double VOLUME SAMPLE 3 µL REAGENT 1 225 µL REAGENT 2 70 µL Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B		<u>DATA PROCESS</u> <u>READ</u> START END MAIN 50 51 SUB 28 29 ENDPOINT LIMIT 3 LINEAR CHECK (%) <u>FACTOR</u> Blank Correction 1.000 <u>PROZONE CHECK</u> START END LIMIT (%) FIRST <input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High SECOND <input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High THIRD	
<u>MONITOR</u> 0 LEVEL POINT 1 SPAN 3.000			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 5,0 mg/dL to linearity limit of 151 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

### Precision:

Mean (mg/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
	28,0	76,1	27,5	75,3
SD	0,25	0,81	1,26	2,04
CV (%)	0,89	1,06	4,60	2,71

**Sensitivity:** 1mg/L = 0,001399 (A).

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

 Correlation coefficient (r)<sup>2</sup>: 0,938.

Regression equation: y= 0,9825x - 1,41606.

The results of the performance characteristics depend on the analyzer used.

### NOTES

1. The reagent 2 presents yellowish coloration due to the peroxidase, but it does not affect its functionality.

2. SPINREACT has instruction sheets for several automatic analyzers.

Instructions for many of them are available on request.

### BIBLIOGRAPHY

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219

### PACKAGING

Ref:TK1001096	Cont.	R1: 10 x 24 mL, R2: 10 x 8 mL
Ref:TK1001097		R1: 2 x 24 mL, R2: 2 x 8 mL

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Item Name HDL		CALIBRATION	
DATA INFORMATION		TYPE	
Units	mg/dL	Linear	
Decimals	1		
ANALYSIS		STANDARD	
Type	END	#1	* #4
W.Length 1	570	#2	#5
		#3	#6
Method DIRECT		NORMAL RANGE	
		SERUM	MALE LOW HIGH FEMALE
CORR SLOPE			
1.000 x +	INTER 0		
Item Name HDL		ASPIRATION	
DATA PROCESS		ABSORBANCE LIMIT	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	LOW -3,000 HIGH 3,000
VOLUME		START END	
SAMPLE	3 µL	MAIN 50 51	
REAGENT 1	225 µL	SUB 28 29	
REAGENT 2	70 µL	ENDPOINT LIMIT 3 LINEAR CHECK (%)	
Third Mix		FACTOR	
R1 Blank	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	Blank Correction 1,000	
	Water <input type="checkbox"/> R1-B <input checked="" type="checkbox"/>		
MONITOR		PROZONE CHECK	
0 LEVEL POINT	1	FIRST	START END LIMIT (%)
SPAN	3,000	SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

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- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 – 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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### BIBLIOGRAPHIE

1. National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
2. Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
3. Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
4. Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
5. Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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	Homens	Mulheres
Risco menor	> 50 mg/dL	> 60 mg/dL
Risco normal	35 - 50 mg/dL	45 - 60 mg/dL
Risco elevado	< 35 mg/dL	< 45 mg/dL

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### BIBLIOGRAFIA

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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