Uric acid-LQ


Uricase -POD. Liquid

Quantitative determination of uric acid IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Uric acid is oxidized by uricase to allantoin and hydrogen peroxide (2H₂O₂), which under the influence of POD, 4-aminophenazone (4-AP) and 2-4 Dichlorophenol sulfonate (DCPS) forms a red quinoneimine compound:

The intensity of the red color formed is proportional to the uric acid concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Uric acid and its salts are end products of the purine metabolism. With progressive renal insufficiency, there is retention in blood of urea, creatinine and uric acid. Elevate uric acid level may be indicative of renal insufficiency and is commonly associated with gout^{1,5,6}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Phosphate pH 7.4	50 mmol/L
Buffer	2-4 Dichlorophenol sulfonate (DCPS)	4 mmol/L
R 2	Uricase	60 U/L
Enzymes	Peroxidase (POD)	660 U/L
	Ascorbate oxidase	200 U/L
	4 - Aminophenazone (4-AP)	1 mmol/L

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 520 nm ≥ 0.16.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

- Serum or plasma¹: Stability 3-5 days at 2-8°C or 6 months at -20°C.
- Urine (24 h)¹: Stability 4 days at 15-25°C, pH >8. Dilute sample 1/50 in distilled water. Mix. Multiply results by 50 (dilution factor); If urine is cloudy; warm the specimen to 60°C for 10 min to dissolve precipitated urates and uric acid. Do not refrigerate.

REFERENCE VALUES¹

Serum or plasma:
 Women 2.5 - 6.8 mg/dL ≅ 149 - 405 μmol/L
 Men 3.6 - 7.7 mg/dL ≅ 214 - 458 μmol/L
 Urine: 250 - 750 mg/24 h ≅ 1.49 - 4.5 mmol/24 h
 These values are for orientation purpose; each laboratory should establish its own reference range.

Conversion factor: mg/dL x 59.5= μmol/L.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210). If control values are found outside the defined range, check the instrument, reagents and technique for problems. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINTECH 240

Item Name URIC			
DATA INFORMATION			
Units	mg/dL		
Decimals	2		
ANALYSIS			
Type	END		
W.Length 1	505		
Method	Uricase-PAP.Liquid		
CORR			
SLOPE	INTER		
1.000 x +	0		
Item Name URIC			
ASPIRATION			
KIND	Single	✓ Double	
VOLUME			
SAMPLE	7	μL	
REAGENT 1	150	μL	
REAGENT 2	150	μL	
Third Mix	✓ OFF	ON	
R1 Blank	Water	✓ R1-B	
MONITOR			
0 LEVEL POINT	1		
SPAN	3.000		
CALIBRATION			
TYPE	Linear		
STANDARD			
#1	*	#4	
#2		#5	
#3		#6	
NORMAL RANGE			
	LOW	HIGH	
SERUM	MALE		
	FEMALE		
URINE			
DATA PROCESS			
READ			
START	END	LOW	-3.000
MAIN 50	52	HIGH	3.000
SUB			
ENDPOINT LIMIT 3			
LINEAR CHECK (%)			
FACTOR			
Blank Correction	1.000		
PROZONE CHECK			
START	END	LIMIT (%)	
FIRST		✓ Low High	
SECOND		✓ Low High	
THIRD		✓ Low High	

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,01647 mg/dL to linearity limit of 40 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/L)	SD	CV (%)	
Mean (mg/L)	4,46	10,37	4,71	11,02
SD	0,02	0,05	0,06	0,15
CV (%)	0,46	0,44	1,20	1,37

Sensitivity: 1 mg/dL = 0,0323 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99734.

Regression equation: y=0,816x + 0,319.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Schultz A. Uric acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1261-1266 and 418.
2. Fossati P et al. Clin Chem 1980;26:227-231.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995

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