
**Quantitative determination of urea**
**IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

 Urea in the sample is hydrolyzed enzymatically into ammonia (NH<sub>4</sub><sup>+</sup>) and carbon dioxide (CO<sub>2</sub>).

 Ammonia ions formed reacts with α-ketoglutarate in a reaction catalysed by glutamate dehydrogenase (GLDH) with simultaneous oxidation of NADH to NAD<sup>+</sup>:

 The decrease in concentration of NADH, is proportional to urea concentration in the sample<sup>1</sup>.

**CLINICAL SIGNIFICANCE**

Urea is the final result of the metabolism of proteins; it is formed in the liver from their destruction.

 It can appear the urea elevated in blood (uremia) in: diets with excess of proteins, renal diseases, heart failure, gastrointestinal hemorrhage, dehydration or renal obstruction<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b>	TRIS pH 7.8	80 mmol/L
Buffer	α-Ketoglutarate	6 mmol/L
	Urease	75000 U/L
<b>R 2</b>	GLDH	60000 U/L
Enzymes	NADH	0.32 mmol/L

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1.00.

**ADDITIONAL EQUIPMENT**

- Autoanalyzer Spintech 240.
- General laboratory equipment<sup>(Note 1)</sup>.

**SAMPLES**

 - Serum or heparinized plasma<sup>1</sup>: Do not use ammonium salts or fluoride as anticoagulants.

 - Urine<sup>1</sup>: Dilute sample 1/50 in distilled water. Mix. Multiply the results by 50 (dilution factor). Preserve urine samples at pH < 4.

Urea is stable at 2-8°C for 5 days.

**REFERENCE VALUES<sup>4,5</sup>**

Serum or plasma:

15-45 mg/dL ≅ 2.5-7.5 mmol/L

Urine:

26 – 43 g/24 h ≅ 428-714 mmol/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

**Conversion factor:** mg/dL x 0.1665 = mmol/L.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**APPLICATION SPINTECH 240**

Item Name UREA			
<b>DATA INFORMATION</b>			
Units	mg/dL		
Decimals	0		
<b>ANALYSIS</b>			
Type	END		
W.Length 1	340		
Method	Urease UV		
<b>CORR</b>			
SLOPE	INTER		
1.000 x +	0		
Item Name UREA			
<b>ASPIRATION</b>			
KIND	Single	✓ Double	
<b>VOLUME</b>			
SAMPLE	3	μL	
REAGENT 1	240	μL	
REAGENT 2	60	μL	
Third Mix	✓ OFF	ON	
R1 Blank	Water	✓ R1-B	
<b>MONITOR</b>			
0 LEVEL POINT	1		
SPAN	3.000		
		<b>DATA PROCESS</b>	
		<b>ABSORBANCE LIMIT</b>	
READ	START	END	LOW HIGH
MAIN	39	40	-3.000 3.000
SUB	35	36	
		<b>ENDPOINT LIMIT 3</b>	
		<b>LINEAR CHECK (%)</b>	
		<b>FACTOR</b>	
		Blank Correction 1.000	
		<b>PROZONE CHECK</b>	
		START	END LIMIT (%)
FIRST			✓ Low High
SECOND			✓ Low High
THIRD			✓ Low High

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

**PERFORMANCE CHARACTERISTICS**
**Measuring range:** From detection limit 0,743 mg/dL to linearity limit 400 mg/dL.

If the concentration is greater than linearity limit dilute 1/2 the sample with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	CV (%)	Mean (mg/dL)	CV (%)
Mean (mg/dL)	37,5	1,20	40,0	1,26
SD	1,05	0,92	1,06	2,07
CV (%)	2,79	0,77	2,65	1,65

**Sensitivity:** 1 mg/dL = 0,00180 A.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 50 samples was the following:

 Correlation coefficient (r)<sup>2</sup>: 0,98209.

Regression equation y= 1,0343x – 1,2105.

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. Glassware and distilled water must be free of ammonia and ammonium salts<sup>1</sup>.
2. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
3. Use clean disposable pipette tips for its dispensation.

**BIBLIOGRAPHY**

1. Kaplan A. Urea. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1257-1260 and 437 and 418.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Item Name UREA <b>DATA INFORMATION</b> Units mg/dL Decimals 0 <b>ANALYSIS</b> Type END W.Length 1 340 Method Urease UV <b>CORR</b> SLOPE INTER 1.000 x + 0		<b>CALIBRATION</b> TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 <b>NORMAL RANGE</b> SERUM MALE LOW HIGH FEMALE URINE	
Item Name UREA <b>ASPIRATION</b> KIND Single <input checked="" type="checkbox"/> Double VOLUME SAMPLE 3 µL REAGENT 1 240 µL REAGENT 2 60 µL Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B <b>MONITOR</b> 0 LEVEL POINT 1 SPAN 3.000		<b>DATA PROCESS</b> <b>READ</b> START END MAIN 39 40 SUB 35 36 ENDPOINT LIMIT 3 LINEAR CHECK (%) <b>FACTOR</b> Blank Correction 1.000 <b>PROZONE CHECK</b> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

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