


**Quantitative determination of cholesterol IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

The cholesterol present in the sample originates a coloured complex, according to the following reactions:

 The intensity of the color formed is proportional to the cholesterol concentration in the sample<sup>1,2</sup>.

**CLINICAL SIGNIFICANCE**

Cholesterol is a fat-like substance called a lipid that is found in all body cells. The liver makes all of the cholesterol the body needs to form cell membranes and to make certain hormones.

The determination of serum cholesterol is one of the important tools in the diagnosis a classification of lipemia.

 High blood cholesterol is one of the major risk factors for heart disease<sup>5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R</b>	PIPES pH 6.9	90 mmol/L
	Phenol	26 mmol/L
	Cholesterol esterase (CHE)	1000 U/L
	Cholesterol oxidase (CHOD)	300 U/L
	Peroxidase (POD)	650 U/L
	4 - Aminophenazone (4-AP)	0.4 mmol/L

**PREPARATION**

The reagent is ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm  $\geq 0.26$ .

**ADDITIONAL EQUIPMENT**

- Autoanalyzer Spintech 240.
- General laboratory equipment.

**SAMPLES**

 Serum or plasma<sup>1,2</sup>: Stability of the sample 7 days at 2-8°C or freezing at -20°C will keep samples stable for 3 months.

**REFERENCE VALUES**

 Risk evaluation<sup>5,6</sup>:

Less than 200 mg/dL	Normal
200-239 mg/dL	Borderline
240 mg/dL and above	High

These values are for orientation purpose; each laboratory should establish its own reference range.

**Conversion factor:** mg/dL x 0.0258= mmol/L.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTRON H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**APPLICATION SPINTECH 240**

Item Name CHOLESTEROL <b>DATA INFORMATION</b> Units mg/dL Decimals 0 <b>ANALYSIS</b> Type END W.Length 1 505 Method CHOD-PAP <b>CORR</b> SLOPE INTER 1.000 x + 0		<b>CALIBRATION</b> TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 <b>NORMAL RANGE</b> LOW HIGH SERUM MALE URINE FEMALE	
Item Name CHOLESTEROL <b>ASPIRATION</b> KIND <input checked="" type="checkbox"/> Single Double VOLUME SAMPLE 3 µL REAGENT 1 300 µL Third Mix <input checked="" type="checkbox"/> OFF ON R1 Blank Water <input checked="" type="checkbox"/> R1-B <b>MONITOR</b> 0 LEVEL POINT 1 SPAN 3.000		<b>DATA PROCESS</b> <b>ABSORBANCE LIMIT</b> <b>READ</b> START END MAIN 50 52 SUB <b>ENDPOINT LIMIT</b> 3 <b>LINEAR CHECK (%)</b> <b>FACTOR</b> Blank Correction 1.000 <b>PROZONE CHECK</b> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

**PERFORMANCE CHARACTERISTICS**
**Measuring range:** From detection limit 0,00 mg/dL to linearity limit 1000 mg/dL.

If the concentration is greater than linearity limit dilute 1/2 the sample with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	SD	Mean	SD
Mean (mg/dL)	99	201	96	197
SD	0,83	1,41	1,75	6,41
CV (%)	0,84	0,70	1,82	3,26

**Sensitivity:** 1 mg/dL = 0,0019 (A).

**Accuracy:** Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagent.

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0,99549.

 Regression equation:  $y=0,911x + 2,624$ .

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. LCF (Lipid Clearing Factor) is integrated in the reagent.
2. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
3. Use clean disposable pipette tips for its dispensation.

**BIBLIOGRAPHY**

1. Naito H.K. Cholesterol. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1194-11206 and 437.
2. Meiatini F. et al. The 4-hydroxybenzoate/4-aminophenazone Chromogenic System. Clin Chem 1978; 24 (12): 2161-2165.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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**BIBLIOGRAPHIE**

1. Naito H.K. Cholesterol. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1194-11206 and 437.
2. Meattini F. et al. The 4-hydroxybenzoate/4-aminophenazone Chromogenic System. Clin Chem 1978; 24 (12): 2161-2165.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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