**Quantitative determination of LDL cholesterol**
**IVD**




Store at 2-8°C

**PRINCIPLE OF THE METHOD**




Direct determination of serum LDLc (low-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation steps.

The assay takes place in two steps.

- 1° Elimination of lipoprotein no-LDL

- 2° Measurement of LDLc

The intensity of the color formed is proportional to the LDLc concentration in the sample.

**CLINICAL SIGNIFICANCE**

 The LDLc particle is lipoproteins that transport cholesterol to the cells. Often called "bad cholesterol" because high levels are risk factor for coronary heart disease and are associated with obesity, diabetes and nephrosis<sup>1,2,9</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b> Enzymes	PIPES Buffer	50 mmol/L
	Cholesterol esterase (CHE)	≥600 U/L
	Cholesterol oxidase (CHOD)	≥500 U/L
	Catalase	≥600 KU/L
	TOOS	2 mmol/L
<b>R 2</b> Enzymes	PIPER Buffer	50 mmol/L
	4 - Aminoantipyrine (4-AA)	4 mmol/L
	Peroxidase (POD)	≥ 4 KU/L

**PREPARATION**

 - **R 1 and R 2:** Are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

 - **R 1 and R 2:** Once opened is stable 4 weeks at 2-8°C.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.

**ADDITIONAL EQUIPMENT**

 - Autoanalyzer Spintech 240.  
 - General laboratory equipment.

**SAMPLES**

 Serum, heparinized plasma or EDTA plasma. If any sample show precipitates, centrifuge before using<sup>5</sup>.

Serum stable 6 days at 2-8°C. Do not freeze the samples.

**REFERENCE VALUES<sup>6,7,8</sup>**

Optimal	< 100 mg/dL
Near or above optimal	100-129 mg/dL
Borderline high	130-160 mg/dL
High	> 160 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

**Conversion factor:** mg/dL x 0,02586 = mmol/L

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**APPLICATION SPINTECH 240**

Item Name LDL		CALIBRATION	
<b>DATA INFORMATION</b>		TYPE Linear	
Units	mg/dL		
Decimals	1		
<b>ANALYSIS</b>		STANDARD	
Type	END	#1	* #4
W.Length 1	570	#2	#5
		#3	#6
Method Direct		<b>NORMAL RANGE</b>	
CORR SLOPE INTER 0		SERUM MALE FEMALE	LOW HIGH
Item Name LDL		DATA PROCESS	
<b>ASPIRATION</b>		<b>ABSORBANCE LIMIT</b>	
KIND	Single <input checked="" type="checkbox"/> Double	READ	LOW -3.000
VOLUME		START END	HIGH 3.000
SAMPLE	4 µL	MAIN 50 51	
REAGENT 1	225 µL	SUB 28 29	
REAGENT 2	75 µL	ENDPOINT LIMIT 3 LINEAR CHECK (%)	
Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON		<b>FACTOR</b>	
R1 Blank	Water <input checked="" type="checkbox"/> R1-B	Blank Correction 1.000	
<b>MONITOR</b>		<b>PROZONE CHECK</b>	
0 LEVEL POINT	1	START END	LIMIT (%)
SPAN	3.000	FIRST	<input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High
		SECOND	<input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High
		THIRD	

*Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.*

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 10 mg/dL to linearity limit of 976 mg/dL. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

Media (mg/dL)	Intraserie (n= 20)		Interserie (n= 20)	
		31,4	67,8	32,1
SD	0,42	1,11	0,92	2,02
CV (%)	1,35	1,64	2,87	2,97

**Sensibility:** 1 mg/dL = 0,001784 (A).

**Accuracy<sup>10,11</sup>:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

 Correlation coefficient (r)<sup>2</sup>: 0,99123.

Regression equation: y= 0,914x + 1,58283

The results of the performance characteristics depend on the analyzer used.

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Item Name LDL <b>DATA INFORMATION</b> Units mg/dL Decimals 1 <b>ANALYSIS</b> Type END W.Length 1 570 Method Direct <b>CORR</b> SLOPE INTER 1.000 x + 0		<b>CALIBRATION</b> TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 <b>NORMAL RANGE</b> LOW HIGH SERUM MALE FEMALE	
Item Name LDL <b>ASPIRATION</b> KIND Single <input checked="" type="checkbox"/> Double VOLUME SAMPLE 4 µL REAGENT 1 225 µL REAGENT 2 75 µL Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B		<b>DATA PROCESS</b> <b>ABSORBANCE LIMIT</b> <b>READ</b> START END MAIN 50 51 SUB 28 29 ENDPOINT LIMIT 3 LINEAR CHECK (%)	
<b>MONITOR</b> 0 LEVEL POINT 1 SPAN 3.000		<b>FACTOR</b> Blank Correction 1.000 <b>PROZONE CHECK</b> START END LIMIT (%) FIRST <input checked="" type="checkbox"/> Low High SECOND <input checked="" type="checkbox"/> Low High THIRD	

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