Quantitative determination of creatinine IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

In the first reaction, creatinase and sarcosine oxidase were used in the enzymatic hydrolysis of endogenous creatine to produce hydrogen peroxide, which is eliminated by catalase. In the second reaction, the catalase is inhibited by sodium azide, and creatinase and 4-aminoantipyrine (4-AA) were added, and only the creatine generated from creatinine by creatininase was hydrolyzed sequentially by creatinase and sarcosine oxidase to produce hydrogen peroxide. This newly-formed hydrogen peroxide was measured in a coupled reaction catalyzed by peroxidase, with N-ethyl-n-sulphopropyl-mtoluidine (TOPS)/4-AA as a chromogen.

CLINICAL SIGNIFICANCE

Creatinine measurements are used in the diagnosis and treatment of renal diseases, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

REAGENTS

R 1	MOPS 25 mmol/L, TOPS 0.5 mmol/L, Creatinase 10 KU/L, Sarcosine Oxidase 5 KU/L Catalase 3 KU/L, EDTA 1mmol/L, pH 7,5.
R 2	MOPS 90 mmol/L, Creatininase 30 KU/L, peroxidase 10 KU/L, pH 7,5. Sodium azide 0,5 g/L.

PREPARATION

R1 and R2 are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

R1 and R2 are stable 8 weeks after opening bottle.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

- Serum or plasma¹.
- Urine (24 h)¹: Dilute fresh urine 1/50 with distilled water.
Multiply the result by 50 (sample dilution factor).
Creatinine is stable 1 day at 2-8°C.

REFERENCE VALUES¹

Serum or plasma:

Men 0,9 - 1,3 mg/dL

Women 0,6 - 1,1 mg/dL

Urine:

Men 14- 26 mg/Kg/24 h

Women 11-20 mg/Kg/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINTECH 240

Item Name CREA			
DATA INFORMATION			
Units	mg/dL		
Decimals	2		
ANALYSIS			
Type	END		
W.Length 1	546		
Method	Trinder		
CORR			
SLOPE	INTER		
1.000 x +	0		
Item Name CREA			
ASPIRATION			
KIND	Single	✓ Double	
VOLUME			
SAMPLE	6	µL	
REAGENT 1	270	µL	
REAGENT 2	90	µL	
Third Mix	✓ OFF	ON	
R1 Blank	Water	✓ R1-B	
MONITOR			
0 LEVEL POINT	1		
SPAN	3.000		
CALIBRATION			
TYPE	Linear		
STANDARD			
#1	*	#4	
#2		#5	
#3		#6	
NORMAL RANGE			
		LOW	HIGH
SERUM			
MALE			
FEMALE			
DATA PROCESS			
READ		ABSORBANCE LIMIT	
START	END	LOW	-3.000
MAIN 50	52	HIGH	3.000
SUB 35	37		
ENDPOINT LIMIT 3			
LINEAR CHECK (%)			
FACTOR			
Blank Correction	1.000		
PROZONE CHECK			
	START	END	LIMIT (%)
FIRST			✓ Low High
SECOND			✓ Low High
THIRD			✓ Low High

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,00 mg/dL to linearity limit of 180 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	0,87	3,82	0,87	3,75
SD	0,01	0,06	0,02	0,06
CV (%)	1,63	1,44	2,31	1,72

Sensitivity: 1 mg/dL = 0,0226 (ΔA)

Accuracy: Results obtained using SPINREACT these reagents did not show systematic differences when compared with other commercial reagents or with HPLC method.

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,9730

Regression equation: y = 1,066x - 0,020.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. Calibration with an aqueous standard may cause matrix related bias, it is recommended to calibrate using a serum based calibrator.

BIBLIOGRAPHY

1. Fossati et al. Clin Chem 1983;29:1494-1496.
2. Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co.,1999.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.

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