Magnesium Xylidyl

Xylidyl Blue. Colorimetric

Quantitative determination of magnesium

IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Magnesium forms a coloured complex when reacts with Magon sulfonate in alkaline solution.

The intensity of the color formed is proportional to the magnesium concentration in the sample¹.

CLINICAL SIGNIFICANCE

Magnesium is the second more abundant intracellular cation of the human body after potassium, being essential in great number of enzymatic and metabolic processes.

Is a cofactor of all the enzymatic reactions that involve the ATP and comprises of the membrane that maintains the electrical excitability of the muscular and nervous cells.

A low magnesium level is found in malabsorption syndrome, diuretic or aminoglycoside therapy; hyperparathyroidism or diabetic acidosis.

Elevated concentration of magnesium is found in uremia, chronic renal failure, glomerulonephritis, Addisons's disease or intensive anti acid therapy^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Xylidyl Blue	0,1 mmol/L
	Thioglycolic acid	0,7 mmol/L
	DMSO	3000 mmol/L

PRECAUTIONS

R: H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles, color change and turbidity.
- Blank absorbance (A) at 546 ≥ 1,8.

ADDITIONAL EQUIPMENT

- Autoanalyser Spintech 240.
- General laboratory equipment (Note 2).

SAMPLES

- Serum, heparinized plasma¹:
Free of hemolysis and separated from cells as rapidly as possible.
Do not use oxalates or EDTA as anticoagulant.
Stability: 7 days at 2-8°C.
- Urine¹:
Should be acidified to pH 1 with HCl.
If urine is cloudy; warm the specimen to 60°C for 10 min. to dissolve precipitates.
Dilute the sample 1/10 with distilled water and multiply the result by 10.
Stability: 3 days at 2-8°C

REFERENCE VALUES¹

Serum or plasma:

1,6 – 2,5 mg/dL ≅ 0,66 – 0,03 mmol/L

Urine:

24 – 244 mg/24 h ≅ 2 – 21 mEq/L/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

Conversion factors:

mg/dL x 0,412 = mmol/L or

0,5 mmol/L = 1,0 mEq/L = 1,22 mg/dL = 12,2 mg/L¹.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINTECH 240

Item Name MG			
DATA INFORMATION			
Units	mg/dL		
Decimals	2		
ANALYSIS			
Type	END		
W.Length 1	546		
Method	Xylidyl Blue		
CORR			
SLOPE	INTER		
1.000 x +	0		
Item Name MG			
ASPIRATION			
KIND	Single	Double	
VOLUME			
SAMPLE	3	μL	
REAGENT 1	300	μL	
Third Mix	OFF	ON	
R1 Blank	Water	R1-B	
MONITOR			
O LEVEL POINT	1		
SPAN	3.000		
		CALIBRATION	
		TYPE	Linear
		STANDARD	
		#1	* #4
		#2	#5
		#3	#6
		NORMAL RANGE	
		LOW	HIGH
		SERUM MALE	
		URINE FEMALE	
		DATA PROCESS	
		READ	LOW -3.000
		START END	HIGH 3.000
		MAIN 26	27
		SUB	
		FACTOR	
		Blank Correction	1.000
		PROZONE CHECK	
		START END	LIMIT (%)
		FIRST	
		SECOND	Low High
		THIRD	Low High

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,0052 mg/dL to linearity limit of 6 mg/dL. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	1,99	3,55	1,98	3,41
SD	0,03	0,04	0,09	0,15
CV (%)	1,68	1,14	4,55	4,42

Sensitivity: 1 mg/dL = 0,5536 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

Correlation coefficient (r)²: 0,92276

Regression equation: y=1,027x + 0,102

The results of the performance characteristics depend on the analyzer used.

NOTES

1. It is recommended use disposable material to avoid magnesium contamination. If glassware is used the material should be scrupulously clean with H₂SO₄ - K₂Cr₂O₇ and then thoroughly rinsed with distilled water and dried before use.
2. Calibration with the aqueous standard may cause a systematic error in automatic procedures. It is recommended to use a serum Calibrator.
3. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Farrell E C. Magnesium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1065-1069.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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