BILIRUBIN D- DPD

Direct Bilirubin

DPD. Colorimetric

Quantitative determination of direct bilirubin IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Direct bilirubin (conjugated) couples with the diazo reagent in the presence of sulfamic acid to form azobilirubin. The intensity of color formed is proportional to the bilirubin concentration in the sample tested. The increase of absorbance at 546 nm is directly proportional to the direct bilirubin concentration.

CLINICAL SIGNIFICANCE

Bilirubin is caused by the degradation of hemoglobin and exists in two forms. Unconjugated bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted. Hyperbilirubinemia is the result of an increase of bilirubin in plasma. Possible causes: **Total bilirubin:** Increase hemolysis, genetic, neonatal jaundice, ineffective erythropoiesis and presence of drugs. **Direct bilirubin:** Hepatic cholestasis, genetic, hepatocellular damage. Clinical diagnosis should not be made based on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Sulfamic acid	100 mM
R 2	2,4-DPD Hydrochloric acid (HCl)	0,5 mM 0,3 M

PRECAUTIONS

R1: H314 - Irritation or skin corrosion. / R2: H290- Corrosive to metals.
H335 - May cause respiratory irritation. H314 - Irritation or skin corrosion.
R2: contains HCl and 2,4-DPD.
Follow the safety advice given in MSDS and product label.

PREPARATION

The reagents are provided in a ready to use format.

STORAGE AND STABILITY

The reagents are stable until the expiry date stated on the label when stored at 2-8°C, protected from light and contaminations are prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis. Protect samples from light.
Stability of the sample: 4 days at 2-8°C or 2 months at -20°C.

REFERENCE VALUES

Direct bilirubin 0 – 0,2 mg/dL (0 – 3,42µmol/L)

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210). If control values are found outside the defined range, check the instrument, reagents and calibrator for problems. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

INTERFERENCES

Interferences from hemolysis, lipemia and ascorbic acid were evaluated for this direct bilirubin method on a Spintech 240 analyzer. Two concentrations of direct bilirubin were evaluated. No interferences were observed for lipemia (Intralipid) up to 350 mg/dL and ascorbic acid up to 40 mg/L. Hemolysis causes decreased direct bilirubin values. A list of drugs and other interfering substances with bilirubin has been reported by Young et. al ^{4,5}.

APPLICATION SPINTECH 240

Item Name BILI D			
DATA INFORMATION		CALIBRATION	
Units	mg/dL	TYPE	Linear
Decimals	2		
ANALYSIS		STANDARD	
Type	END	#1 *	#4
		#2	#5
		#3	#6
W.Length 1	546	NORMAL RANGE	
W.Length 2			LOW HIGH
Method	DPD	SERUM MALE	
		URINE FEMALE	
CORR SLOPE			
1.000 x +	INTER 0		
Item Name BILI D			
ASPIRATION		DATA PROCESS	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	ABSORBANCE LIMIT
		START END	LOW -3,000
SAMPLE	VOLUME	MAIN 50 51	HIGH 3,000
REAGENT 1	240 µL	SUB 30 31	
REAGENT 2	60 µL		ENDPOINT LIMIT 3
			LINEAR CHECK (%)
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	FACTOR	
R1 Blank	<input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B	Blank Correction	1,000
MONITOR		PROZONE CHECK	
0 LEVEL POINT	1	START END LIMIT (%)	
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

The Calibration is stable until **7 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,03 mg/dL to linearity limit of 9 mg/dL. If the results obtained are greater than the linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (mg/dL)	Inter assay (n= 40)		Intra assay (n= 80)	
		0,7458	2,444	0,7458
SD	0,05868	0,0550	0,0276	0,024
CV (%)	7,9	2,2	3,7	1,0

Sensitivity: 1 mg/dL = 0,040 Abs. units

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x) on a Spintech 240 analyzer. The results obtained using 53 samples ranging from 0,06 a 9 mg/dL (1,02 to 153,9 µmol/L) were:
Correlation coefficient (r): 0,9986
Regression equation: y= 1,0056 x – 0,1046
The results of the performance characteristics depend on the analyzer used.

NOTES

- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307–328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: TK1001047	Cont.	R 1: 10 x 25 mL
		R 2: 10 x 7 mL

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BIBLIOGRAPHIE

- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Text book of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: TK1001047	Cont.	R 1: 10 x 25 mL R 2: 10 x 7 mL
----------------	-------	-----------------------------------

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