# CK-MB-LQ (Creatine kinase-MB)

Anti CK-M. Immunoinhibition. Kinetic UV. Liquid

## Quantitative determination of creatine kinase MB (CK-MB)

### IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

The procedure involves measurement of CK activity in the presence of an antibody to CK-M monomer. This antibody completely inhibits the activity of CK-MM and half of the activity of CK-MB while not affecting the B subunit activity of CK-MB and CK-BB. Then it's used the CK method to quantitatively determine CK-B activity<sup>1,2</sup>. The CK-MB activity is obtained by multiplying the CK-B activity by two.

### CLINICAL SIGNIFICANCE

CK-MB is an enzyme formed by the association of two subunits from muscle (M) and nerve cells (B). CK-MB is usually present in serum at low concentration; it is increased after an acute infarct of myocardium and later descends at normal levels. Also, is increased, rarely, in skeletal muscle damage<sup>5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

R 1	Imidazol, pH 6.7	125 mmol/L
	D-Glucose	25 mmol/L
	N-Acetyl-L-Cysteine	25 mmol/L
	Magnesium acetate	12.5 mmol/L
	NADP	2.52 mmol/L
	EDTA	2.02 mmol/L
	Hexokinase	≥8 800 U/L
Anti-human polyclonal CK-M antibody (sheep) sufficient to inhibit up to 2 000 U/L of CK-MM		
R 2	ADP	15.2 mmol/L
	AMP	25 mmol/L
	di-Adenosine-5- pentaphosphate	103 mmol/L
	Glucose-6-phosphate dehydrogenase	≥8 800 U/L
	Creatine phosphate	250 mmol/L
OPTIONAL	CK-Nac / CK-MB CONTROL	Ref: 1002260

### PREPARATION

All the reagents are ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1.60.

### ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

### SAMPLES

Serum free of hemolysis or heparin plasma<sup>1</sup>: Stability 7 days at 2-8°C, protected from light.

CK-MB activity decreases a 10% after 24 hours at 4°C or 1 hour at 25°C.

### Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1.00	1.53	2.38
30°C	0.65	1.00	1.56
37°C	0.42	0.64	1.00

### REFERENCE VALUES

The suspicion of myocardial damage is based on the three following factors:

CK-MB	25°C	30°C	37°C
	> 10 U/L	> 15 U/L	> 24 U/L
TOTAL CK	25°C	30°C	37°C
	Men, up to 80 U/L	130 U/L	195 U/L
Women, up to 70 U/L	110 U/L	170 U/L	

$\frac{\text{CK-MB Activity}}{\text{CK Total Activity}} \times 100$ : 6-25 % CK-MB Activity in the sample

These values are for orientation purpose; each laboratory should establish its own reference range.

## APPLICATION SPINTECH 240

Item Name CK-MB		CALIBRATION	
DATA INFORMATION		TYPE	Linear
Units	U/L		
Decimals	0		
ANALYSIS		STANDARD	
Type	END	#1	* #4
		#2	#5
		#3	#6
W.Length 1	340	NORMAL RANGE (37°C)	
		LOW	HIGH
Method	Kinetic UV.Liq	SERUM MALE	24 195
		FEMALE	24 170
CORR			
SLOPE	INTER		
1.000 x +	0		
Item Name CK-MB		DATA PROCESS	
ASPIRATION		ABSORBANCE LIMIT	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	LOW -3.000
		START	END
SAMPLE	VOLUME	MAIN	50 52
REAGENT 1	12 µL	SUB	35 37
REAGENT 2	240 µL		
REAGENT 2	60 µL	ENDPOINT LIMIT 3	
		LINEAR CHECK (%)	
Third Mix <input type="checkbox"/> OFF <input type="checkbox"/> ON		FACTOR	
R1 Blank <input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B		Blank Correction	1.000
MONITOR		PROZONE CHECK	
0 LEVEL POINT	1	START	END
SPAN	3.000	LIMIT (%)	
		FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

### QUALITY CONTROL

CK-Nac/CK-MB specific control sera (Ref. 1002260) are recommended to monitor the performance of assay procedures.

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of 1,9 U/L to *linearity limit* of 1000 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/1 with NaCl 9 g/L and multiply the result by 2.

### Precision:

**Sensitivity:** 1 U/L = 0,000134 (A).

Mean (U/L)	Intra-assay		Inter-assay	
	33,7	166,5	31,3	161,0
SD	1,00	3,76	1,19	3,47
CV (%)	2,96	2,26	3,81	2,15

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

Correlation coefficient (r)<sup>2</sup>: 0,999.

Regression equation: y = 0,976 x - 0,269.

The results of the performance characteristics depend on the analyzer used.

### LIMITATION OF THE PROCEDURE

- The method will also measure any CK-BB isoenzyme present in serum. The activity of the isoenzyme is usually negligible, however, if a significant amount of CK-BB activity is present the CK-MB activity will be overestimated.
- A macro form of BB (immunoglobulin complexed) has been observed which will be measured as B in the assay. If the measured CK-B activity exceeds 20% of the total CK activity, the presence of macro BB should be suspected.

### BIBLIOGRAPHY

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

### PACKAGING

Ref: TK41254  Cont. R1: 10 x 25 mL, R2: 10 x 7 mL  
 Ref: TK41255 R1: 2 x 25 mL, R2: 2 x 7 mL

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- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: TK41254	Cont.	R1: 10 x 25 mL, R2: 10 x 7 mL
Ref: TK41255		R1: 2 x 25 mL, R2: 2 x 7 mL

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### BIBLIOGRAPHIE

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: TK41254  Cont. R1: 10 x 25 mL, R2: 10 x 7 mL  
Ref: TK41255 R1: 2 x 25 mL, R2: 2 x 7 mL

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