### Quantitative determination of phosphorus

#### IVD

Store at 2-8°C

#### PRINCIPLE OF THE METHOD

Direct method for determining inorganic phosphate. Inorganic phosphate reacts in acid medium with ammonium molybdate to form a phosphomolybdate complex with yellow colour. The intensity of the color formed is proportional to the inorganic phosphorus concentration in the sample<sup>1,2</sup>.

#### CLINICAL SIGNIFICANCE

Phosphorus is an essential mineral for tissue bone formation and is required by every cell in the body for normal function. Approximately 85% of the body phosphorus is found in bone and in teeth. Low levels of phosphorus can be caused by hypervitaminosis D, primary hyperparathyroidism, renal tubular disorders, antacids or malabsorption. High levels of phosphorus can be caused by diet, bone metastases, liver disease, alcohol ingestion, diarrhea and vomiting<sup>1,5,6</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

#### REAGENTS

<b>R</b>	Ammonium molybdate	0.40 mM
Molybdic	Sulphuric acid (SO <sub>4</sub> H <sub>2</sub> )	210 mM
	Detergents	

#### PRECAUTIONS

R: H290-May be corrosive to metals. H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

#### PREPARATION

The reagent is ready to use.

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

#### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm  $\geq$  0.54.

#### ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment (Note 1).

#### SAMPLES

- Serum or plasma<sup>1,5</sup>:

Free of hemolysis. Serum or plasma should be removed from the clot as quickly as possible to avoid elevation of serum phosphorus from hydrolysis or leakage of phosphate present in erythrocytes. Stability: 7 days at 2-8°C.

- Urine<sup>1,2</sup> (24 h):

Collect the specimen into a bottle containing 10 mL of 10% v/v hydrochloric acid (HCl) to avoid phosphate precipitations. Adjust to pH 2. Dilute the sample 1/10 with distilled water. Mix. Multiply the result by 10 (dilution factor). Stability: 10 days at 2-8°C.

#### REFERENCE VALUES<sup>1</sup>

Serum or plasma:

Children	4,0 – 7,0 mg/dL	$\cong$	1,29 – 2,26 mmol/L
Adults	2,5 – 5,0 mg/dL	$\cong$	0,80 – 1,61 mmol/L

Urine:

Adults	0.4 – 1.3 g /24 h
--------	-------------------

These values are for orientation purpose; each laboratory should establish its own reference range.

**Conversion factor:** mg/dL x 0.323 = mmol/L.

#### QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

Item Name PHOS			
<u>DATA INFORMATION</u>			
Units	mg/dL		
Decimals	2		
<u>ANALYSIS</u>			
Type	END		
W.Length 1	340		
Method	UV		
<u>CORR</u>			
SLOPE	INTER		
1.000 x +	0		
Item Name PHOS			
<u>ASPIRATION</u>			
KIND	<input checked="" type="checkbox"/> Single	<input type="checkbox"/> Double	
<u>VOLUME</u>			
SAMPLE	3	$\mu$ L	
REAGENT 1	300	$\mu$ L	
Third Mix	<input checked="" type="checkbox"/> OFF	<input type="checkbox"/> ON	
R1 Blank	Water	<input checked="" type="checkbox"/> R1-B	
<u>MONITOR</u>			
0 LEVEL POINT	1		
SPAN	3.000		
		<u>CALIBRATION</u>	
		TYPE	Linear
		<u>STANDARD</u>	
		#1	*
		#2	#4
		#3	#5
		#4	#6
		<u>NORMAL RANGE</u>	
		LOW	HIGH
		SERUM	MALE
			FEMALE
		URINE (mg/24h)	
		<u>DATA PROCESS</u>	
		<u>ABSORBANCE LIMIT</u>	
		READ	LOW -3.000
		START	HIGH 3.000
		END	
		MAIN 30	31
		SUB	
		ENDPOINT LIMIT 3	
		LINEAR CHECK (%)	
		<u>FACTOR</u>	
		Blank Correction	1.000
		<u>PROZONE CHECK</u>	
		START	END
		LIMIT (%)	
		FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

#### APPLICATION SPINTECH 240

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

#### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 0,000 mg/dL to linearity limit of 35 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	4,09	7,12	4,11	7,09
SD	0,03	0,046	0,09	0,06
CV (%)	0,62	0,80	2,15	0,80

**Sensitivity:** 1 mg/dL = 0,0798 A.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,8577.

Regression equation: y = 0,724x + 0,837.

The results of the performance characteristics depend on the analyzer used.

#### NOTES

- Most of the detergents and water softening products used in the laboratories contain chelating agents and phosphates. It is recommended to rinse glassware in diluted nitric acid and water before using.
- Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- Use clean disposable pipette tips for its dispensation.

#### BIBLIOGRAPHY

- Farrell E C. Phosphorus. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1072-1074 and 418.
- Daly J A. et al. Clin Chem 1972; 18 (3): 263-265.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
- Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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Soro ou plasma:  
Crianças 4,0 – 7,0 mg/dL  $\cong$  1,29 – 2.26 mmol/L  
Adultos 2,5 – 5,0 mg/dL  $\cong$  0.80 – 1.61 mmol/L

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Item Name PHOS				
<b>DATA INFORMATION</b>			<b>CALIBRATION</b>	
Units	mg/dL		TYPE	Linear
Decimals	2		STANDARD	
<b>ANALYSIS</b>			#1	*
Type	END		#2	#5
W.Length 1	340		#3	#6
<b>METHOD</b>			<b>NORMAL RANGE</b>	
Method	UV		SERUM	MALE
<b>CORR</b>				FEMALE
SLOPE	INTER		URINE (mg/24h)	
1.000 x +	0			
Item Name PHOS				
<b>ASPIRATION</b>			<b>DATA PROCESS</b>	
KIND	Single	Double	<b>READ</b>	<b>ABSORBANCE LIMIT</b>
<b>VOLUME</b>			START	END
SAMPLE	3	$\mu$ L	MAIN 30	31
REAGENT 1	300	$\mu$ L	SUB	
<b>THIRD MIX</b>			ENDPOINT LIMIT 3	
R1 Blank	Water	R1-B	LINEAR CHECK (%)	
<b>MONITOR</b>			<b>FACTOR</b>	
0 LEVEL POINT	1		Blank Correction	1.000
SPAN	3.000		<b>PROZONE CHECK</b>	
			START	END
			LIMIT (%)	
			FIRST	
			SECOND	<input checked="" type="checkbox"/> Low High
			THIRD	<input checked="" type="checkbox"/> Low High

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