
**Quantitative determination of glucose IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

 Glucose oxidase (GOD) catalyses the oxidation of glucose to gluconic acid. The formed hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>), is detected by a chromogenic oxygen acceptor, phenol, 4 – aminophenazone (4-AP) in the presence of peroxidase (POD):

 The intensity of the color formed is proportional to the glucose concentration in the sample<sup>1,2</sup>.

**CLINICAL SIGNIFICANCE**

Glucose is a major source of energy for most cells of the body; insulin facilitates glucose entry into the cells.

 Diabetes is a disease manifested by hyperglycemia; patients with diabetes demonstrate an inability to produce insulin<sup>1,5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R</b>	TRIS pH 7.4	92 mmol/L
	Phenol	0.3 mmol/L
	Glucose oxidase (GOD)	15000 U/L
	Peroxidase (POD)	1000 U/L
	4 – Aminophenazone (4-AP)	2.6 mmol/L

**PREPARATION**

The reagent is ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

**SIGNS OF REAGENT DETERIORATION:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm  $\geq$  0.32.

**ADDITIONAL EQUIPMENT**

- Autoanalyzer Spintech 240.
- General laboratory equipment.

**SAMPLES**

 Serum or plasma, free of hemolysis<sup>1</sup>:  
 Serum should be removed from the clot as quickly as possible.  
 Stability of the sample: Glucose in serum or plasma is stable at 2-8°C for 3 days.

**REFERENCE VALUES<sup>1</sup>**

 Serum or plasma:  
 60 – 110 mg/dL  $\cong$  3.33 – 6.10 mmol/L  
 These values are for orientation purpose; each laboratory should establish its own reference range.

**Conversion factor:** mg/dL x 0.0555= mmol/L.

**QUALITY CONTROL**

Contro sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**APPLICATION SPINTECH 240**

Item Name GLU			
<b>DATA INFORMATION</b>			
Units	mg/dL		
Decimals	0		
<b>ANALYSIS</b>			
Type	END		
W.Length 1	505		
<b>METHOD</b>			
Method	GOD-PAP		
<b>CORR</b>			
SLOPE	INTER		
1.000 x +	0		
Item Name GLU			
<b>ASPIRATION</b>			
KIND	<input checked="" type="checkbox"/> Single	<input type="checkbox"/> Double	
<b>VOLUME</b>			
SAMPLE	3 $\mu$ L		
REAGENT 1	300 $\mu$ L		
<b>Third Mix</b>			
R1 Blank	<input checked="" type="checkbox"/> OFF	<input type="checkbox"/> ON	
	Water	<input checked="" type="checkbox"/> R1-B	
<b>MONITOR</b>			
0 LEVEL POINT	1		
SPAN	3.000		
		<b>DATA PROCESS</b>	
		<b>ABSORBANCE LIMIT</b>	
		READ	LOW -3.000
		START	HIGH 3.000
		MAIN 50	
		END 52	
		SUB	
		<b>ENDPOINT LIMIT 3</b>	
		<b>LINEAR CHECK (%)</b>	
<b>FACTOR</b>			
		Blank Correction	1.000
<b>PROZONE CHECK</b>			
		START	END
		LIMIT (%)	
		FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

**PERFORMANCE CHARACTERISTICS**
**Measuring range:** From detection limit 0,3709 mg/dL to linearity limit 500 mg/dL.

If the concentration is greater than linearity limit dilute 1/2 the sample with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	98,5	264,6	92,5	250
SD	0,5754	1,2733	2,76	6,44
CV (%)	0,59	0,48	2,98	2,57

**Sensitivity:** 1 mg/dL = 0,0039 (A).

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 50 samples were the following:

 Correlation coefficient (r)<sup>2</sup>: 0,99492.

Regression equation: y=1,104x - 1,249.

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

**BIBLIOGRAPHY**

1. Kaplan L.A. Glucose. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1032-1036.
2. Trinder P. Ann Clin Biochem 1969; 6 24-33.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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**APPLICATION AU SPINTECH 240**

<p>Item Name GLU</p> <p><u>DATA INFORMATION</u></p> <p>Units mg/dL</p> <p>Decimals 0</p> <p><u>ANALYSIS</u></p> <p>Type END</p> <p>W.Length 1 505</p> <p>Method GOD-PAP</p> <p><u>CORR</u></p> <p>SLOPE INTER</p> <p>1.000 x + 0</p>	<p><u>CALIBRATION</u></p> <p>TYPE Linear</p> <p>STANDARD</p> <table style="width: 100%; text-align: center;"> <tr> <td>#1</td> <td>*</td> <td>#4</td> </tr> <tr> <td>#2</td> <td></td> <td>#5</td> </tr> <tr> <td>#3</td> <td></td> <td>#6</td> </tr> </table> <p><u>NORMAL RANGE</u></p> <table style="width: 100%; text-align: center;"> <tr> <td>SERUM</td> <td>MALE</td> <td>LOW</td> <td>HIGH</td> </tr> <tr> <td>URINE</td> <td>FEMALE</td> <td></td> <td></td> </tr> </table>	#1	*	#4	#2		#5	#3		#6	SERUM	MALE	LOW	HIGH	URINE	FEMALE					
#1	*	#4																			
#2		#5																			
#3		#6																			
SERUM	MALE	LOW	HIGH																		
URINE	FEMALE																				
<p>Item Name GLU</p> <p><u>ASPIRATION</u></p> <p>KIND <input checked="" type="checkbox"/> Single Double</p> <p>3.000</p> <p><u>VOLUME</u></p> <p>SAMPLE 3 µL</p> <p>REAGENT 1 300 µL</p> <p>Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON</p> <p>R1 Blank Water <input checked="" type="checkbox"/> R1-B</p> <p><u>MONITOR</u></p> <p>0 LEVEL POINT 1</p> <p>SPAN 3.000</p>	<p><u>DATA PROCESS</u></p> <p><u>READ</u></p> <table style="width: 100%; text-align: center;"> <tr> <td>START</td> <td>END</td> <td>LOW</td> <td>HIGH</td> </tr> <tr> <td>MAIN 50</td> <td>52</td> <td></td> <td>3.000</td> </tr> </table> <p>SUB</p> <p><u>ENDPOINT LIMIT</u> 3</p> <p><u>LINEAR CHECK (%)</u></p> <p><u>FACTOR</u></p> <p>Blank Correction 1.000</p> <p><u>PROZONE CHECK</u></p> <table style="width: 100%; text-align: center;"> <tr> <td>START</td> <td>END</td> <td>LIMIT (%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td><input checked="" type="checkbox"/> Low <input type="checkbox"/> High</td> </tr> <tr> <td>THIRD</td> <td></td> <td><input checked="" type="checkbox"/> Low <input type="checkbox"/> High</td> </tr> </table>	START	END	LOW	HIGH	MAIN 50	52		3.000	START	END	LIMIT (%)	FIRST			SECOND		<input checked="" type="checkbox"/> Low <input type="checkbox"/> High	THIRD		<input checked="" type="checkbox"/> Low <input type="checkbox"/> High
START	END	LOW	HIGH																		
MAIN 50	52		3.000																		
START	END	LIMIT (%)																			
FIRST																					
SECOND		<input checked="" type="checkbox"/> Low <input type="checkbox"/> High																			
THIRD		<input checked="" type="checkbox"/> Low <input type="checkbox"/> High																			

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<table border="1"> <tr><td>Item Name</td><td>GLU</td></tr> <tr><td colspan="2"><b>DATA INFORMATION</b></td></tr> <tr><td>Units</td><td>mg/dL</td></tr> <tr><td>Decimals</td><td>0</td></tr> <tr><td colspan="2"><b>ANALYSIS</b></td></tr> <tr><td>Type</td><td>END</td></tr> <tr><td>W.Length 1</td><td>505</td></tr> <tr><td>Method</td><td>GOD-PAP</td></tr> <tr><td colspan="2"><b>CORR</b></td></tr> <tr><td>SLOPE</td><td>INTER</td></tr> <tr><td>1.000 x +</td><td>0</td></tr> </table>	Item Name	GLU	<b>DATA INFORMATION</b>		Units	mg/dL	Decimals	0	<b>ANALYSIS</b>		Type	END	W.Length 1	505	Method	GOD-PAP	<b>CORR</b>		SLOPE	INTER	1.000 x +	0	<table border="1"> <tr><td colspan="2"><b>CALIBRATION</b></td></tr> <tr><td>TYPE</td><td>Linear</td></tr> <tr><td colspan="2"><b>STANDARD</b></td></tr> <tr><td>#1</td><td>*</td><td>#4</td></tr> <tr><td>#2</td><td></td><td>#5</td></tr> <tr><td>#3</td><td></td><td>#6</td></tr> <tr><td colspan="2"><b>NORMAL RANGE</b></td></tr> <tr><td></td><td>LOW</td><td>HIGH</td></tr> <tr><td>SERUM</td><td>MALE</td><td>FEMALE</td></tr> <tr><td>URINE</td><td></td><td></td></tr> </table>	<b>CALIBRATION</b>		TYPE	Linear	<b>STANDARD</b>		#1	*	#4	#2		#5	#3		#6	<b>NORMAL RANGE</b>			LOW	HIGH	SERUM	MALE	FEMALE	URINE																								
Item Name	GLU																																																																						
<b>DATA INFORMATION</b>																																																																							
Units	mg/dL																																																																						
Decimals	0																																																																						
<b>ANALYSIS</b>																																																																							
Type	END																																																																						
W.Length 1	505																																																																						
Method	GOD-PAP																																																																						
<b>CORR</b>																																																																							
SLOPE	INTER																																																																						
1.000 x +	0																																																																						
<b>CALIBRATION</b>																																																																							
TYPE	Linear																																																																						
<b>STANDARD</b>																																																																							
#1	*	#4																																																																					
#2		#5																																																																					
#3		#6																																																																					
<b>NORMAL RANGE</b>																																																																							
	LOW	HIGH																																																																					
SERUM	MALE	FEMALE																																																																					
URINE																																																																							
<table border="1"> <tr><td>Item Name</td><td>GLU</td></tr> <tr><td colspan="2"><b>ASPIRATION</b></td></tr> <tr><td>KIND</td><td><input checked="" type="checkbox"/> Single <input type="checkbox"/> Double</td></tr> <tr><td colspan="2"><b>VOLUME</b></td></tr> <tr><td>SAMPLE</td><td>3 <math>\mu</math>L</td></tr> <tr><td>REAGENT 1</td><td>300 <math>\mu</math>L</td></tr> <tr><td>Third Mix</td><td><input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON</td></tr> <tr><td>R1 Blank</td><td>Water <input checked="" type="checkbox"/> R1-B</td></tr> <tr><td colspan="2"><b>MONITOR</b></td></tr> <tr><td>0 LEVEL POINT</td><td>1</td></tr> <tr><td>SPAN</td><td>3.000</td></tr> </table>	Item Name	GLU	<b>ASPIRATION</b>		KIND	<input checked="" type="checkbox"/> Single <input type="checkbox"/> Double	<b>VOLUME</b>		SAMPLE	3 $\mu$ L	REAGENT 1	300 $\mu$ L	Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	R1 Blank	Water <input checked="" type="checkbox"/> R1-B	<b>MONITOR</b>		0 LEVEL POINT	1	SPAN	3.000	<table border="1"> <tr><td colspan="2"><b>DATA PROCESS</b></td><td colspan="2"><b>ABSORBANCE LIMIT</b></td></tr> <tr><td colspan="2"><b>READ</b></td><td>LOW</td><td>-3.000</td></tr> <tr><td>START</td><td>END</td><td>HIGH</td><td>3.000</td></tr> <tr><td>MAIN</td><td>50 52</td><td></td><td></td></tr> <tr><td>SUB</td><td></td><td></td><td></td></tr> <tr><td colspan="2"><b>FACTOR</b></td><td colspan="2"><b>ENDPOINT LIMIT 3</b></td></tr> <tr><td>Blank Correction</td><td>1.000</td><td colspan="2">LINEAR CHECK (%)</td></tr> <tr><td colspan="2"><b>PROZONE CHECK</b></td><td colspan="2"></td></tr> <tr><td>START</td><td>END</td><td colspan="2">LIMIT (%)</td></tr> <tr><td>FIRST</td><td></td><td colspan="2"><input checked="" type="checkbox"/> Low High</td></tr> <tr><td>SECOND</td><td></td><td colspan="2"><input checked="" type="checkbox"/> Low High</td></tr> <tr><td>THIRD</td><td></td><td colspan="2"><input checked="" type="checkbox"/> Low High</td></tr> </table>	<b>DATA PROCESS</b>		<b>ABSORBANCE LIMIT</b>		<b>READ</b>		LOW	-3.000	START	END	HIGH	3.000	MAIN	50 52			SUB				<b>FACTOR</b>		<b>ENDPOINT LIMIT 3</b>		Blank Correction	1.000	LINEAR CHECK (%)		<b>PROZONE CHECK</b>				START	END	LIMIT (%)		FIRST		<input checked="" type="checkbox"/> Low High		SECOND		<input checked="" type="checkbox"/> Low High		THIRD		<input checked="" type="checkbox"/> Low High	
Item Name	GLU																																																																						
<b>ASPIRATION</b>																																																																							
KIND	<input checked="" type="checkbox"/> Single <input type="checkbox"/> Double																																																																						
<b>VOLUME</b>																																																																							
SAMPLE	3 $\mu$ L																																																																						
REAGENT 1	300 $\mu$ L																																																																						
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON																																																																						
R1 Blank	Water <input checked="" type="checkbox"/> R1-B																																																																						
<b>MONITOR</b>																																																																							
0 LEVEL POINT	1																																																																						
SPAN	3.000																																																																						
<b>DATA PROCESS</b>		<b>ABSORBANCE LIMIT</b>																																																																					
<b>READ</b>		LOW	-3.000																																																																				
START	END	HIGH	3.000																																																																				
MAIN	50 52																																																																						
SUB																																																																							
<b>FACTOR</b>		<b>ENDPOINT LIMIT 3</b>																																																																					
Blank Correction	1.000	LINEAR CHECK (%)																																																																					
<b>PROZONE CHECK</b>																																																																							
START	END	LIMIT (%)																																																																					
FIRST		<input checked="" type="checkbox"/> Low High																																																																					
SECOND		<input checked="" type="checkbox"/> Low High																																																																					
THIRD		<input checked="" type="checkbox"/> Low High																																																																					

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