


**Quantitative determination of triglycerides IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Sample triglycerides incubated with lipoproteinlipase (LPL), liberate glycerol and free fatty acids. Glycerol is converted to glycerol-3-phosphate (G3P) and adenosine-5-diphosphate (ADP) by glycerol kinase (GK) and ATP. Glycerol-3-phosphate (G3P) is then converted by glycerol phosphate oxidase (GPO) to dihydroxyacetone phosphate (DAP) and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>).

In the last reaction, hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) reacts with 4-aminophenazone (4-AP) and p-chlorophenol in presence of peroxidase (POD) to give a red colored dye:

The intensity of the color formed is proportional to the triglycerides concentration in the sample<sup>1,2,3</sup>.

**CLINICAL SIGNIFICANCE**

Triglycerides are fats that provide energy for the cell.

Like cholesterol, they are delivered to the body's cells by lipoproteins in the blood. A diet with a lot of saturated fats or carbohydrates will raise the triglyceride levels. The increases in serum triglycerides are relatively non-specific. For example, liver dysfunction resulting from hepatitis, extra hepatic biliary obstruction or cirrhosis, diabetes mellitus is associated with the increase<sup>3,6,7</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R</b>	GOOD pH 6.3	50 mmol/L
	p-Chlorophenol	2 mmol/L
	Lipoprotein lipase (LPL)	150000 U/L
	Glycerol kinase (GK)	500 U/L
	Glycerol-3-oxidase (GPO)	3500 U/L
	Peroxidase (POD)	440 U/L
	4 - Aminophenazone (4-AP)	0.1 mmol/L
	ATP	0.1 mmol/L

**PREPARATION**

The reagent is ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm  $\geq$  0.26.

**ADDITIONAL EQUIPMENT**

- Autoanalyzer Spintech 240.
- General laboratory equipment.

**SAMPLES**

Serum or plasma<sup>1</sup>.

Stability of the sample: Triglycerides are stable for 5 days at 2-8°C.

**REFERENCE VALUES**

Men 40 – 160 mg/dL  
 Women 35 – 135 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

**Conversion factor:** mg/dL x 0.0113= mmol/L.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**APPLICATION SPINTECH 240**

Item Name TG <b>DATA INFORMATION</b> Units mg/dL Decimals 0 <b>ANALYSIS</b> Type END W.Length 1 505 Method GPO-POD <b>CORR</b> SLOPE INTER 1.000 x + 0		<b>CALIBRATION</b> TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 <b>NORMAL RANGE</b> SERUM MALE LOW HIGH FEMALE	
Item Name TG <b>ASPIRATION</b> KIND <input checked="" type="checkbox"/> Single Double VOLUME SAMPLE 3 µL REAGENT 1 300 µL Third Mix <input checked="" type="checkbox"/> OFF ON R1 Blank Water <input checked="" type="checkbox"/> R1-B <b>MONITOR</b> 0 LEVEL POINT 1 SPAN 3.000		<b>DATA PROCESS</b> <b>ABSORBANCE LIMIT</b> <b>READ</b> LOW -3.000 START END HIGH 3.000 MAIN 50 52 SUB ENDPOINT LIMIT 3 LINEAR CHECK (%) <b>FACTOR</b> Blank Correction 1.000 <b>PROZONE CHECK</b> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than **40 days**.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From *detection limit* 0,000 mg/dL to *linearity limit* 1200 mg/dL. If the concentration is greater than linearity limit, dilute 1/2 the sample with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	SD	CV (%)	CV (%)
Mean (mg/dL)	109	224	111	224
SD	0,64	1,01	3,74	7,90
CV (%)	0,58	0,45	3,38	3,52

**Sensitivity:** 1 mg/dL = 0,0013 (A).

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99810.

Regression equation: y= 0,9178x - 0,5426

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. LCF (Lipid Clearing Factor) is integrated in the reagent.
2. Calibration with the aqueous Standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
3. Use clean disposable pipette tips for its dispensation.

**BIBLIOGRAPHY**

1. Buccolo G et al. Quantitative determination of serum triglycerides by use of enzymes. Clin Chem 1973; 19 (5): 476-482.
2. Fossati P et al. Clin. Chem 1982; 28(10): 2077-2080.
3. Kaplan A et al. Tryglycerides. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 437 and Lipids 1194-1206.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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1. Buccolo G et al. Quantitative determination of serum triglycerides by use of enzymes. Clin Chem 1973; 19 (5): 476-482.
2. Fossati P et al. Clin. Chem 1982; 28(10): 2077-2080.
3. Kaplan A et al. Tryglicerides. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 437 and Lipids 1194-1206.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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