Total Bilirubin

DPD. Colorimetric

Quantitative determination of total bilirubin IVD

Store at 2-8°C

INTENDED USE

Bilirubin T-DPD is an in vitro diagnostic medical device intended to be used by healthcare professionals for the quantitative determination of total bilirubin in serum or plasma.

PRINCIPLE OF THE METHOD

Bilirubin (both conjugated and unconjugated) couples with the diazo reagent in the presence of a surfactant to form azobilirubin. The intensity of color formed is proportional to the bilirubin concentration in the sample tested. The increase of absorbance at 546 nm is directly proportional to the total bilirubin concentration.

CLINICAL SIGNIFICANCE

Bilirubin is caused by the degradation of hemoglobin and exists in two forms. Unconjugated bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted. Hyperbilirubinemia is the result of an increase of bilirubin in plasma.

Possible causes:

Total bilirubin: Increase hemolysis, genetic alteration, neonatal anemia, erythropoiesis alterations and presence of drugs.

Direct Bilirubin: cholestasis liver, liver abnormalities and genetic.

Clinical diagnosis should not be made based on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Surfactants	<1%
	Hydrochloric acid (HCl)	160 mM
R 2	2,4-DPD	≥2 mM
	Hydrochloric acid (HCl)	120 mM
	Surfactant	<1%

PRECAUTIONS

R1/R2: H314-Causes severe skin burns and eye damage. Follow the safety advice given in MSDS and product label.

PREPARATION

The reagents are provided in a ready to use format.

STORAGE AND STABILITY

The reagents are stable until the expiry date stated on the label when stored at 2-8°C, protected from light and contaminations are prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

The use of fresh serum and plasma free of hemolysis is recommended. Protect samples from light.

The samples can be stored in the refrigerator (2-8°C) for 4 days maximum. For longer storage periods they must be kept frozen at -20°C for 2 months. Frozen samples should be totally thawed and brought to room temperature before testing. Avoid repeated freezing and thawing of the samples.

REFERENCE VALUES

Total bilirubin 0,2-1,2 mg/dL (3,4 – 20,5 µmol/L)

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210). If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINTECH 240

Item Name BILI T			
DATA INFORMATION			
Units	mg/dL		
Decimals	2		
ANALYSIS			
Type	END		
W.Length 1	546		
W.Length 2			
Method	DPD		
CORR			
SLOPE	INTER		
1.000 x +	0		
Item Name BILI T			
ASPIRATION			
KIND	Single	<input checked="" type="checkbox"/> Double	
VOLUME			
SAMPLE	12	µL	
REAGENT 1	240	µL	
REAGENT 2	60	µL	
Third Mix	<input checked="" type="checkbox"/> OFF	ON	
R1 Blank	<input checked="" type="checkbox"/> Water	R1-B	
MONITOR			
0 LEVEL POINT	1		
SPAN	3.000		
		DATA PROCESS	
		ABSORBANCE LIMIT	
		READ	LOW -3,000
		START	END
		MAIN 50	51
		SUB 30	31
		ENDPOINT LIMIT 3	
		LINEAR CHECK (%)	
		FACTOR	
		Blank Correction	1,000
		PROZONE CHECK	
		START	END
		FIRST	LIMIT (%)
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

The Calibration is stable until **7 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

Measuring range: From quantification limit of 0,1 mg/dL to linearity limit of 30 mg/dL.

If the results obtained were greater than the linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Inter assay (n= 20)		Intra assay (n= 20)	
	Mean (mg/dL)	SD	Mean (mg/dL)	SD
Mean (mg/dL)	1,24	4.4	1,22	4.33
SD	0,01	0,03	0,02	0,08
CV (%)	1,0	0,7	2,0	1,9

Sensitivity: 1 mg/dL = 0,026 Abs. units.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x) on a Spintech 240 analyzer. The results obtained using 60 samples were:

Correlation coefficient (r): 0,996

Regression equation: $y = 0,9836x + 0,1644$

The results of the performance characteristics depend on the analyzer used.

NOTES

- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2: 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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1. Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
2. Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: TK1001046	Cont.	R 1: 10 x 25 mL
		R 2: 10 x 7 mL

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BIBLIOGRAPHIE

- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Text book of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: TK1001046	Cont.	R 1: 10 x 25 mL R 2: 10 x 7 mL
----------------	-------	-----------------------------------

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BIBLIOGRAFIA

- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: TK1001046

Cont.

R 1: 10 x 25 mL
R 2: 10 x 7 mL