α₁-Acid Glycoprotein

Turbidimetry

Quantitative determination of α₁-Acid Glycoprotein (α₁-Ac GLY) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Anti- α₁-Ac Gly antibodies when mixed with samples containing α₁-Ac Gly, form insoluble complexes. These complexes cause an absorbance change, dependent upon the α₁-Ac Gly concentration of the patient sample, that can be quantified by comparison from a calibrator of known α₁-Ac Gly concentration.

CLINICAL SIGNIFICANCE

α₁-Ac Glycoprotein (also known as orosomucoid) is a glycoprotein synthesized by hepatic parenchymal cells, but granulocytes and monocytes may also contribute significantly to plasma levels in sepsis. It has long been known to bind a large number of basic and lipophilic compounds (progesterone and related hormones).

It is an acute phase response protein that shows a 3 to 4-fold increase in most conditions associated with inflammation or tissue necrosis, and may be one of the most reliable indicators of clinical activity of ulcerative colitis. Levels also are increased by glucocorticoid effect.

Synthesis and plasma levels are decreased by estrogens.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH, 8,3. Preservative.
R 2 Antibody	Goat serum, anti-human α ₁ -Ac Glycoprotein, pH 7.5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the α₁-Ac Glycoprotein calibrator by the corresponding factor stated in table below to obtain the α₁-Ac Glycoprotein concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (μL)	--	10	25	50	75	100
NaCl 9 g/L (μL)	100	90	75	50	25	-
Factor	0	0,1	0,25	0,5	0,75	1,0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

REFERENCE VALUES

Between 50 - 120 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

SPINTECH 240 / 300 APPLICATION

Item Name GLYC		<u>CALIBRATION</u>	
<u>DATA INFORMATION</u>		TYPE Spline	
Units	mg/dL		
Decimals	0		
<u>ANALYSIS</u>		STANDARD	
Type	END	#1 0.10 x Cal.Val.	#4 0.75 x Cal.Val.
		#2 0.25 x Cal.Val.	#5 1.00 x Cal.Val.
W.Length 1	340	#3 0.50 x Cal.Val.	#6
Method Turbidimetry		<u>NORMAL RANGE</u>	
<u>CORR SLOPE</u>		SERUM	LOW HIGH
1.000 x +	INTER 0	MALE	
		FEMALE	
Item Name GLYC		<u>DATA PROCESS</u>	
<u>ASPIRATION</u>		<u>ABSORBANCE LIMIT</u>	
KIND	Single <input checked="" type="checkbox"/> Double	<u>READ</u>	LOW -3.000
		START END	HIGH 3.000
	VOLUME**	MAIN 50 51	
SAMPLE	2.5 μL	SUB 30 31	
REAGENT 1	300 μL	ENDPOINT LIMIT 3	
REAGENT 2	75 μL	LINEAR CHECK (%)	
Third Mix <input checked="" type="checkbox"/> OFF ON		<u>FACTOR</u>	
R1 Blank <input checked="" type="checkbox"/> Water R1-B		Blank Correction	1.000
<u>MONITOR</u>		<u>PROZONE CHECK</u>	
0 LEVEL POINT	1	START END LIMIT (%)	
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days.

PERFORMANCE CHARACTERISTICS

- Measurement range:** Up to 250 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection Limit:** Values less than 12,9 mg/dL give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 1000 mg/dL.
- Sensitivity:** Δ 5.0 mA / mg/dL.
- Precision:** The reagent has been tested for 20 days, using two levels of serum in a EP5-based study.

EP5	CV (%)	
	56,4 mg/dL	112,07 mg/dL
Total	3,3%	3,1%
Within Run	1,1%	1,6%
Between Run	3%	2,1%
Between Day	0,7%	1,6%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the method Immage from Beckman. 51 samples ranging from 50 to 120 mg/dL of α₁-Ac Gly were assayed. The correlation coefficient (r) was 0.95 and the regression equation

$$y = 0,9304x + 6,5367.$$

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Bienvenue J et al. Clin Chem Clin Biochem 1981; 27: 721-726.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACCPres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACCPres, 1997.

PRESENTATION

Ref: TK1102154 Cont. R1. Diluent: 2 x 24 mL
R2. Antibody: 2 x 6 mL

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BIBLIOGRAPHIE

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Bienvenue J et al. Clin Chem Clin Biochem 1981; 27: 721-726.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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Item Name GLYC DATA INFORMATION Units mg/dL Decimals 0 ANALYSIS Type END W.Length 1 340 Method Turbidimetry CORR SLOPE INTER 1.000 x + 0		CALIBRATION TYPE Spline STANDARD #1 0.10 x Cal.Val. #4 0.75 x Cal.Val. #2 0.25 x Cal.Val. #5 1.00 x Cal.Val. #3 0.50 x Cal.Val. #6 NORMAL RANGE SERUM MALE LOW HIGH FEMALE	
Item Name GLYC ASPIRATION KIND Single <input checked="" type="checkbox"/> Double VOLUME** SAMPLE 2.5 µL REAGENT 1 300 µL REAGENT 2 75 µL Third Mix <input checked="" type="checkbox"/> OFF ON R1 Blank <input checked="" type="checkbox"/> Water R1-B MONITOR 0 LEVEL POINT 1 SPAN 3.000		DATA PROCESS ABSORBANCE LIMIT READ START END LOW -3.000 HIGH 3.000 MAIN 50 51 SUB 30 31 ENDPOINT LIMIT 3 LINEAR CHECK (%) FACTOR Blank Correction 1.000 PROZONE CHECK START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

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