# Transferrin

## Turbidimetry

### Quantitative determination of Transferrin (TRF)

#### IVD

Store 2 - 8°C.

#### PRINCIPLE OF THE METHOD

Anti-transferrin antibodies when mixed with samples containing TRF, form insoluble complexes. These complexes cause an absorbance change, dependent upon the TRF concentration of the patient sample, that can be quantified by comparison from a calibrator of known TRF concentration.

#### CLINICAL SIGNIFICANCE

Transferrin is a plasma protein that contains a single polypeptide chain with approximately 6% carbohydrate. It is synthesized in the liver and transfers iron through the serum.

Evaluation of plasma TRF levels is useful for the differential diagnosis of anemia and for monitoring its treatment. In the hypochromic anemia (iron deficiency), the TRF level is increased. On the other hand, if the anemia is due to a failure to incorporate iron into erythrocytes, the TRF level is normal or low but the protein is highly saturated with iron. In iron overload, the TRF concentration is normal but saturation exceeds 55% and may be as great as 90%.

TRF concentration may, in fact, be used for assessing nutritional status. In the congenital defect atransferrinemia, very low level of TRF is accompanied by iron overload and a severe hypochromic anemia that is resistant to iron therapy. High levels of TRF occur in pregnancy and during estrogen administration.

#### REAGENTS

<b>R 1</b>	Tris buffer 20 mmol/L, PEG 8000, pH 8.3.
Diluent	Preservative.
<b>R 2</b>	Goat serum, anti-human transferrin, pH 7.5.
Antibody	Preservative.
<b>Optional</b>	Ref: 1102003 PROT CAL.

#### CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

#### PREPARATION

**Reagents:** Ready to use.

**Calibration Curve:** Prepare the following PROT CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the transferrin calibrator by the corresponding factor stated in table below to obtain the transferrin concentration of each dilution:

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	6,25	12,5	25	50	100
NaCl 9 g/L (µL)	100	93,75	87,5	75	50	-
Factor	0	0,0625	0,125	0,25	0,5	1,0

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

#### ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

#### SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

#### REFERENCE VALUES

Between 200 - 360 mg/dL. Each laboratory should establish its own reference range.

#### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### SPINTECH 240 / 300 APPLICATION

Item Name TRF		CALIBRATION	
<u>DATA INFORMATION</u>		TYPE Spline	
Units	mg/dL		
Decimals	1		
<u>ANALYSIS</u>		STANDARD	
Type	END	#1 0.0625 x Cal.Val.	#4 0.5 x Cal.Val.
		#2 0.125 x Cal.Val.	#5 1.00 x Cal.Val.
		#3 0.25 x Cal.Val.	#6
W.Length 1	340	<u>NORMAL RANGE</u>	
		LOW HIGH	
Method	Turbidimetry	SERUM	MALE FEMALE
SLOPE	INTER		
1.000 x +	0		
Item Name TRF		<u>DATA PROCESS</u>	
<u>ASPIRATION</u>		<u>ABSORBANCE LIMIT</u>	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	LOW -3.000
		START END	HIGH 3.000
SAMPLE	VOLUME**	MAIN 44 45	
REAGENT 1	240 µL	SUB 30 31	
REAGENT 2	60 µL		ENDPOINT LIMIT 3
			LINEAR CHECK (%)
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	<u>FACTOR</u>	
R1 Blank	Water <input checked="" type="checkbox"/> R1-B <input type="checkbox"/>	Blank Correction	1.000
<u>MONITOR</u>		<u>PROZONE CHECK</u>	
0 LEVEL POINT	1	START END LIMIT (%)	
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

\*\* Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days.

#### PERFORMANCE CHARACTERISTICS

- Measurement range:** Up to the calibrator value (aprox. 800 mg/dL) under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection Limit:** 0 mg/dL
- Prozone effect:** No prozone effect was detected upon 2000 mg/dL
- Sensitivity:** Δ 3,0 mA / mg/dL (94 mg/dL).
- Precision:** The reagent has been tested for 20 days, using three levels of serum in an EP5-based study.

EP5	CV (%)		
	77,02 mg/dL	206,99 mg/dL	377 mg/dL
Total	5,4%	2,5%	5,4%
Within Run	1%	0,8%	1,2%
Between Run	1,7%	1,3%	2,1%
Between Day	5%	2%	4,9%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the Immage from Beckman. 100 samples ranging from 50 to 700 mg/dL of TRF were assayed. The correlation coefficient (r) was 0,95 and the regression equation  $y = 1,046x + 3,843$ .

The results of the performance characteristics depend on the used analyzer.

#### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

#### BIBLIOGRAPHY

- McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutz HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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### BIBLIOGRAFIA

- McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutzer HJJ. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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**APPLICATION AU SPINTECH 240 / 300**

Item Name TRF <b>DATA INFORMATION</b> Units mg/dL Decimals 1 <b>ANALYSIS</b> Type END W.Length 1 340 Method Turbidimetry SLOPE 1.000 x + INTER 0		<b>CALIBRATION</b> TYPE Spline STANDARD #1 0.0625 x Cal.Val. #4 0.5 x Cal.Val. #2 0.125 x Cal.Val. #5 1.00 x Cal.Val. #3 0.25 x Cal.Val. #6 <b>NORMAL RANGE</b> SERUM MALE LOW HIGH FEMALE	
Item Name TRF <b>ASPIRATION</b> KIND Single <input checked="" type="checkbox"/> Double VOLUME** SAMPLE 3 µL REAGENT 1 240 µL REAGENT 2 60 µL Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B <b>MONITOR</b> 0 LEVEL POINT 1 SPAN 3.000		<b>DATA PROCESS</b> <b>READ</b> START END MAIN 44 45 SUB 30 31 ENDPOINT LIMIT 3 LINEAR CHECK (%) <b>FACTOR</b> Blank Correction 1.000 <b>PROZONE CHECK</b> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

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**BIBLIOGRAPHIE**

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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EP5	CV (%)		
		77,02 mg/dL	206.99 mg/dL
Total	5,4%	2,5%	5,4%
Within Run	1%	0,8%	1,2%
Between Run	1,7%	1,3%	2,1%
Between Day	5%	2%	4,9%

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**BIBLIOGRAFIA**

- McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutzer H.JH. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS.Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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Ref: TK1102134	Cont.	R1.Diluyente:	2 x 24 mL
		R2.Anticorpo:	2 x 6 mL