Quantitative determination of complement C3 (C3) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-human C3 antibodies when mixed with samples containing C3, form insoluble complexes. These complexes cause an absorbance change, dependent upon the C3 concentration of the patient sample, that can be quantified by comparison from a calibrator of known C3 concentration.

CLINICAL SIGNIFICANCE¹

C3 is the functional link between classical and alternative pathways of activation and it is the most concentrate component of the complement system in human plasma. Hepatic cells synthesize C3, although bacterial endotoxins induce synthesis by monocytes and fibroblasts. Concentration C3 increases as a consequence of an acute-phase response (trauma, surgery or inflammatory process), biliary obstruction and focal glomerulosclerosis. Decreasing C3 levels are consequence of a genetic deficiency that may increase the risk of infections particularly with encapsulated bacteria, or acquired deficiency that causes vascular disorders and severe infections.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
Antibody (R2)	Goat serum, anti-human C3, pH 7.5. Preservative.
Optional	Cod: 1102003 PROT CAL.

CALIBRATION

The assay is calibrated to the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the C3 calibrator by the corresponding factor stated in table below to obtain the C3 concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not use. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged. Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES⁵

Neonates: Between 70 - 196 mg/dL.

Adults: Between 90 - 180 mg/dL.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

1. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

SPINTECH 240 / 300 APPLICATION

Item Name C3		CALIBRATION	
DATA INFORMATION		TYPE Spline	
Units	mg/dL		
Decimals	0		
ANALYSIS		STANDARD	
Type	ENDPOINT	#1 0.10 x Cal. Val.	#4 0.75 x Cal. Val.
W.Length 1	340	#2 0.25 x Cal. Val.	#5 1.00 x Cal. Val.
		#3 0.50 x Cal. Val.	#6
Method Turbidimetry		NORMAL RANGE	
		LOW	HIGH
SLOPE 1.000 x +		SERUM	MALE FEMALE
	INTER 0	URINE	
Item Name C3		DATA PROCESS	
ASPIRATION		ABSORBANCE LIMIT	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	LOW -3.000
		START END	HIGH 3.000
	VOLUME**	MAIN 38 39	
SAMPLE	3 µL	SUB 30 31	
REAGENT 1	240 µL	ENDPOINT LIMIT 3	
REAGENT 2	60 µL	LINEAR CHECK (%)	
Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON		FACTOR	
R1 Blank	Water <input checked="" type="checkbox"/> R1-B	Blank Correction 1.000	
MONITOR		PROZONE CHECK	
0 LEVEL POINT	1	START END LIMIT (%)	
SPAN	3.000	FIRST	<input checked="" type="checkbox"/> Low High
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days.

PERFORMANCE CHARACTERISTICS

1. **Measurement range:** Up to 400 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and re-tested again. The linearity limit and measurement range depends on the sample to reagent / ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. **Detection Limit:** Values less than 1 mg/dL give non-reproducible results.
3. **Prozone effect:** No prozone effect was detected upon 1500 mg/dL.
4. **Sensitivity:** Δ 8,86 mA. mg/dL (23,8 mg/dL), Δ 4,3 mA. mg/dL (190 mg/dL).
5. **Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	42,98 mg/dL	118,96 mg/dL	229,5 mg/dL
Total	6,6%	2,3%	3,1%
Within Run	0,9%	0,8%	0,8%
Between Run	3,7%	2,2%	1,8%
Between Day	5,4%	0%	2,4%

6. **Accuracy:** Results obtained using this reagent (y) were compared to those obtained using an immunoturbidimetric method from Bayer. 48 samples ranging from 50 to 200 mg/dL of C3 were assayed. The correlation coefficient (r) was 0,96 and the regression equation $y = 1,1x - 0,6$.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Carrol MC. Annual Review of Immunology 1998; 16: 545-568.
3. Lambris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1: 16-45.
4. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
5. Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
6. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
7. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAPHIE

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Carrol MC. Annual Review of Immunology 1998; 16: 545-568.
- Lambris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1: 16-45.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Carroll MC. Annual Review of Immunology 1998; 16: 545-568.
- Lambris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1: 16-45.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997

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Ref.: TK1102094 Cont. R1. Diluente: 2 x 24 mL
R2. Antibody: 2 x 6 mL