Quantitative determination of complement C4 (C4) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-human C4 antibodies when mixed with samples containing C4, form insoluble complexes. These complexes cause an absorbance change, dependent upon the C4 concentration of the patient sample, that can be quantified by comparison from a calibrator of known C4 concentration.

CLINICAL SIGNIFICANCE¹

C4 is the second component reacting in the classical pathway cascade. Most synthesis occurs in the hepatic parenchymal cells, although some may be synthesized by monocytes or other tissues.

C4 levels in plasma rise modestly after trauma or inflammation and tissue necrosis (acute phase process).

Inherited primary deficiency of C4 is associated with a high prevalence of autoimmune or collagen vascular disease, particularly Systemic Lupus Erythematosus (SLE). Also, levels of C4 are more commonly depressed because of consumption as a consequence of formed immune-complexes.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
Antibody (R2)	Goat serum, anti-human C4, pH 7.5. Preservative.
Optional	Cod: 1102003 PROT CAL.

CALIBRATION

The assay is calibrated to the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the C4 calibrator by the corresponding factor stated in table below to obtain the C4 concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not use.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant.

Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES⁵

Neonates: Between 13 - 38 mg/dL.

Adults: Between 10 - 40 mg/dL.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

1. The linearity depends on the calibrator concentration.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

SPINTECH 240 / 300 APPLICATION

Item Name C4		CALIBRATION	
<u>DATA INFORMATION</u>		TYPE	Spline
Units	mg/dL		
Decimals	0		
<u>ANALYSIS</u>		STANDARD	
Type	ENDPOINT	#1 0.10 x Cal. Val	#4 0.75 x Cal. Val
		#2 0.25 x Cal. Val	#5 1.00 x Cal. Val
W.Length 1	340	#3 0.50 x Cal. Val	#6
		<u>NORMAL RANGE</u>	LOW HIGH
Method	Turbidimetry	SERUM MALE	
		URINE FEMALE	
SLOPE	INTER		
1.000 x +	0		
Item Name C4		DATA PROCESS	
<u>ASPIRATION</u>		<u>ABSORBANCE LIMIT</u>	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	LOW -3.000
		START END	HIGH 3.000
SAMPLE	VOLUME**	MAIN 38 39	
REAGENT 1	240 µL	SUB 30 31	
REAGENT 2	60 µL		ENDPOINT LIMIT 3
			LINEAR CHECK (%)
		<u>FACTOR</u>	1.000
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON		
R1 Blank	Water <input checked="" type="checkbox"/> R1-B <input type="checkbox"/>		
<u>MONITOR</u>		<u>PROZONE CHECK</u>	
0 LEVEL POINT	1	START END LIMIT (%)	
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

*** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.*

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days.

PERFORMANCE CHARACTERISTICS

1. **Measurement range:** Up to 75 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and re-tested again. The linearity limit and measurement range depends on the sample to reagent / ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. **Detection Limit:** Values less than 1 mg/dL give non-reproducible results.
3. **Prozone effect:** No prozone effect was detected upon 500 mg/dL.
4. **Sensitivity:** Δ 23,6 mA. mg/dL (5 mg/dL), Δ 12,9 mA. mg/dL (37mg/dL).
5. **Precision:** The reagent has been tested for 20 days, using three levels of serum in an EP5-based study.

EP5	CV (%)		
	8,57 mg/dL	22,46 mg/dL	42,98 mg/dL
Total	3,9%	2,4%	1,9%
Within Run	1,6%	1%	1%
Between Run	2,2%	1,6%	1,1%
Between Day	2,8%	1,4%	1,2%

6. **Accuracy:** Results obtained using this reagent (y) were compared to those obtained using an immunoturbidimetric method from Bayer. 46 samples ranging from 9 to 60 mg/dL of C4 were assayed. The correlation coefficient (r) was 0,97 and the regression equation $y = 1,16x - 1,9$. The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Yang Y et al. Curr Dir Autoimmun 2004; 7: 98-132.
3. Borque L et al. Clin Biochem 1983; 16: 330-333.
4. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
5. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
6. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
7. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Yang Y et al. Curr Dir Autoimmun 2004; 7: 98-132.
- Borque L et al. Clin Biochem 1983; 16: 330-333.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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