Quantitative determination of Haptoglobin (HAPTO) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-haptoglobin antibodies when mixed with samples containing haptoglobin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the haptoglobin concentration of the patient sample, that can be quantified by comparison from a calibrator of known haptoglobin concentration.

CLINICAL SIGNIFICANCE

The haptoglobin is an α_2 -glycoprotein synthesized in the liver that binds hemoglobin irreversibly. The hapto-hemoglobin complexes, as well as free haptoglobin itself, play significant roles in the iron storage and prevents of possible renal damage as a consequence of hemoglobin excretion. As an acute-phase protein, haptoglobin is increased in the presence of acute inflammatory process, tissue necrosis or malignancy.

Haptoglobin deficiency in plasma is a consequence of hemolysis "in vivo", presence of estrogens in pregnancy and oral contraceptive therapy, as well as most forms of acute or chronic hepatocellular disease, including viral hepatitis.

Haptoglobin test is mainly used for the determination and monitoring the hemolytic disorders. Under normal circumstances, approximately 1% of circulating red blood cells are destroyed from circulation. An increase of only 2% destruction per day will completely deplete plasma haptoglobin in the absence of a stimulus to production, such as acute inflammation or corticosteroids therapy.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
Antiserum (R2)	Goat serum, anti-human haptoglobin pH 7.5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL Calibrator to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the haptoglobin calibrator by the corresponding factor stated in table below to obtain the haptoglobin concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (μ L)	--	10	25	50	75	100
NaCl 9 g/L (μ L)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES?

Between 30 - 200 mg/dL. Each laboratory should establish its own reference range

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

1. Linearity depends on the calibrator concentration.
2. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

SPINTECH 240 / 300 APPLICATION

Item Name HAPTO		CALIBRATION	
DATA INFORMATION		TYPE	Spline
Units	mg/dL		
Decimals	1		
ANALYSIS		STANDARD	
Type	END	#1 0.10 x Cal.Val.	#4 0.75 x Cal.Val.
W.Length 1	340	#2 0.25 x Cal.Val.	#5 1.00 x Cal.Val.
		#3 0.50 x Cal.Val.	#6
Method		NORMAL RANGE	
	Turbidimetry	LOW	HIGH
SLOPE		SERUM MALE	
1.000 x +	INTER	FEMALE	
	0	URINE	
Item Name HAPTO		DATA PROCESS	
ASPIRATION		ABSORBANCE LIMIT	
KIND	Single <input checked="" type="checkbox"/> Double	READ	LOW -3.000
		START END	HIGH 3.000
	VOLUME**	MAIN 40 41	
SAMPLE	3 μ L	SUB 28 29	
REAGENT 1	240 μ L		ENDPOINT LIMIT 3
REAGENT 2	60 μ L		LINEAR CHECK (%)
Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON		FACTOR	
R1 Blank	Water <input checked="" type="checkbox"/> R1-B	Blank Correction	1.000
MONITOR		PROZONE CHECK	
0 LEVEL POINT	1	START END LIMIT (%)	
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days.

PERFORMANCE CHARACTERISTICS

1. **Measurement range:** Up to 300 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. **Detection Limit:** Values less than 1,3 mg/dL give non-reproducible results.
3. **Prozone effect:** No prozone effect was detected upon 1200 mg/dL
4. **Sensitivity:** Δ 4,96 mA / mg/dL (100 mg/dL).
5. **Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	39,25 mg/dL	97,35 mg/dL	191,5 mg/dL
Total	8%	3,2%	2,3%
Within Run	1,5%	0,9%	1,2%
Between Run	6,7%	2,3%	1,2%
Between Day	4%	2%	1,5%

6. **Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the method from Beckman (System Array 360 CE). 35 samples ranging from 10 to 400 mg/dL of Haptoglobin were assayed. The correlation coefficient (r) was 0,95 and the regression equation $y = 0,88x + 4,8$. The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Kluthe R et al. Nature 1965; 205: 93-94.
5. Tarukosky PH et al. Scand J Clin Lab Invest 1966; 18: 80-86.
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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