Immunoglobulin M

Turbidimetry

Quantitative determination of Immunoglobulin M (IgM)

IVD

Store at 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-human IgM antibodies when mixed with samples containing IgM, form insoluble complexes. These complexes cause an absorbance change, dependent upon the IgM concentration of the patient sample, that can be quantified by comparison from a calibrator of known IgM concentration.

CLINICAL SIGNIFICANCE

IgM is the only immunoglobulin that a neonate normally synthesizes, and in adults, it represents the 5-10% of the total immunoglobulins. Its structure is a pentamer of five IgG molecules and its high molecular weight (900.000 daltons) prevents its passage into extravascular spaces.

IgM concentration is decreased in diseases related with hereditary or acquired deficiencies of the immunoglobulin production. Polyclonal increases in serum immunoglobulins are the normal response to infections. The IgM generally increases as a primary response to virus infections and blood stream infections such as malaria and primary biliary cirrhosis. In multiple myeloma, if the paraprotein proves to be IgM, the diagnosis is probably Waldenström macroglobulinemia.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
R 2 Antibody	Goat serum, anti-human IgM, pH 7.5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay is calibrated to the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL Calibrator to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the IgM calibrator by the corresponding factor stated in table below to obtain the IgM concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0,1	0,25	0,5	0,75	1,0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Between 40 - 230 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (cod.:1102004). Each laboratory should establish its own quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

SPINTECH 240 / 300 APPLICATION

Item Name IgM		<u>CALIBRATION</u>	
<u>DATA INFORMATION</u>		TYPE	Spline
Units	mg/dL		
Decimals	1		
<u>ANALYSIS</u>		STANDARD	
Type	END	#1 0.10 x Cal. Val.	#4 0.75 x Cal. Val.
		#2 0.25 x Cal. Val.	#5 1.00 x Cal. Val.
		#3 0.50 x Cal. Val.	#6
W.Length 1	340	<u>NORMAL RANGE</u>	
Method	Turbidimetry	LOW	HIGH
<u>CORR</u>		SERUM	MALE
SLOPE	INTER		FEMALE
1.000 x +	0		
Item Name IgM		<u>DATA PROCESS</u>	
<u>ASPIRATION</u>		<u>ABSORBANCE LIMIT</u>	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	LOW -3.000
		START END	HIGH 3.000
	VOLUME**	MAIN 40 41	
SAMPLE	2 µL	SUB 28 29	
REAGENT 1	240 µL		
REAGENT 2	60 µL		ENDPOINT LIMIT 3
			LINEAR CHECK (%)
		<u>FACTOR</u>	Blank Correction 1.000
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	<u>PROZONE CHECK</u>	
R1 Blank	Water <input checked="" type="checkbox"/> R1-B <input type="checkbox"/>		
<u>MONITOR</u>			
0 LEVEL POINT	1	START END LIMIT (%)	
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

**Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days.

PERFORMANCE CHARACTERISTICS

- 1. Measurement range:** Up to 300 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and re-tested again. The linearity limit and measurement range depends on the sample to reagent / ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- 2. Detection Limit:** Values less than 1 mg/dL give non-reproducible results.
- 3. Prozone effect:** No prozone effect was detected upon 2000 mg/dL.
- 4. Sensitivity:** Δ 2,4 mA. mg/dL at 30 mg/dL.
- 5. Precision:** The reagent has been tested for 20 days, using two levels of serum in a EP5-based study.

EP5	CV (%)	
	68,67 mg/dL	143,3 mg/dL
Total	5,7%	2,8%
Within Run	1,1%	0,7%
Between Run	3,8%	2,3%
Between Day	4,2%	1,3%

- 6. Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the Elecsys system from Roche. 100 samples ranging from 50 to 210 mg/dL of IgM were assayed. The correlation coefficient (r) was 0.958 and the regression equation $y = 0,974x + 1,296$.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
5. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAPHIE

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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