Quantitative determination of α₁-antitrypsin (α₁-ATRYP)

IVD
Store at 2-8°C

PRINCIPLE OF THE METHOD

Anti- α₁-antitrypsin antibodies when mixed with samples containing α₁-antitrypsin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the α₁-antitrypsin concentration of the patient sample, that can be quantified by comparison from a calibrator of known α₁-antitrypsin concentration.

CLINICAL SIGNIFICANCE

α₁-antitrypsin is a glycoprotein synthesized in the hepatic parenchyma cells that circulates in the bloodstream. It is the second highest proteinase inhibitor in plasma after α₂-macroglobulin. α₁-antitrypsin is a strong reactor with elastase, skin collagenase, chymotrypsin, plasmin, and thrombin, and also shows inhibitory activity against fungal and leukocytic proteases.

α₁-antitrypsin deficiency is an inherited disorder, and occurs when both parents pass on an abnormal gene (PIZ) to their child. This deficiency is associated with a very risk for development of pulmonary emphysema and diseases of the liver - neonatal cholestasis, hepatitis, cirrhosis and hepatocellular carcinoma.

α₁-antitrypsin increases in inflammatory or necrosis process. Serum levels begin to rise after approximately 24 hours and peak at 3 or 4 days if the insult is acute and short-lived.

REAGENTS

R 1	
Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
R 2	
Antibody	Goat serum, anti-human α ₁ -antitrypsin, pH 7.5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the α₁-antitrypsin calibrator by the corresponding factor stated in table below to obtain the α₁-antitrypsin concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (μL)	--	10	25	50	75	100
NaCl 9 g/L (μL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

REFERENCE VALUES

Newborn: Between 124 - 348 mg/dL.

Adults: 90 - 200 mg/dL.

Each laboratory should establish its own reference range.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data

SPINTECH 240 / 300 APPLICATION

<p>Item Name ATRY</p> <p>DATA INFORMATION</p> <p>Units mg/dL Decimals 0</p> <p>ANALYSIS</p> <p>Type END</p> <p>W.Length 1 340</p> <p>Method Turbidimetry</p> <p>CORR</p> <p>SLOPE INTER 1.000 x + 0</p>	<p>CALIBRATION</p> <p>TYPE Spline</p> <p>STANDARD</p> <table style="width: 100%; border: none;"> <tr> <td>#1 0.10 x Cal.Val.</td> <td>#4 0.75 x Cal.Val.</td> </tr> <tr> <td>#2 0.25 x Cal.Val.</td> <td>#5 1.00 x Cal.Val.</td> </tr> <tr> <td>#3 0.50 x Cal.Val.</td> <td>#6</td> </tr> </table> <p>NORMAL RANGE</p> <table style="width: 100%; border: none;"> <tr> <td>SERUM</td> <td>MALE</td> <td>LOW</td> <td>HIGH</td> </tr> <tr> <td></td> <td>FEMALE</td> <td>90</td> <td>200 U/L</td> </tr> <tr> <td></td> <td></td> <td>90</td> <td>200 U/L</td> </tr> </table>	#1 0.10 x Cal.Val.	#4 0.75 x Cal.Val.	#2 0.25 x Cal.Val.	#5 1.00 x Cal.Val.	#3 0.50 x Cal.Val.	#6	SERUM	MALE	LOW	HIGH		FEMALE	90	200 U/L			90	200 U/L													
#1 0.10 x Cal.Val.	#4 0.75 x Cal.Val.																															
#2 0.25 x Cal.Val.	#5 1.00 x Cal.Val.																															
#3 0.50 x Cal.Val.	#6																															
SERUM	MALE	LOW	HIGH																													
	FEMALE	90	200 U/L																													
		90	200 U/L																													
<p>Item Name ATRY</p> <p>ASPIRATION</p> <p>KIND Single <input checked="" type="checkbox"/> Double</p> <p>VOLUME**</p> <p>SAMPLE 2 μL REAGENT 1 240 μL REAGENT 2 60 μL</p> <p>Third Mix <input checked="" type="checkbox"/> OFF ON R1 Blank <input checked="" type="checkbox"/> Water R1-B</p> <p>MONITOR</p> <p>0 LEVEL POINT 1 SPAN 3.000</p>	<p>DATA PROCESS</p> <p>ABSORBANCE LIMIT</p> <table style="width: 100%; border: none;"> <tr> <td>READ</td> <td>START</td> <td>END</td> <td>LOW</td> <td>-0.1</td> </tr> <tr> <td>MAIN</td> <td>40</td> <td>41</td> <td>HIGH</td> <td>3.000</td> </tr> <tr> <td>SUB</td> <td>29</td> <td>30</td> <td></td> <td></td> </tr> </table> <p>ENDPOINT LIMIT 3 LINEAR CHECK (%)</p> <p>FACTOR</p> <p>Blank Correction 1.000</p> <p>PROZONE CHECK</p> <table style="width: 100%; border: none;"> <tr> <td></td> <td>START</td> <td>END</td> <td>LIMIT (%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td><input checked="" type="checkbox"/> Low High</td> </tr> <tr> <td>THIRD</td> <td></td> <td></td> <td><input checked="" type="checkbox"/> Low High</td> </tr> </table>	READ	START	END	LOW	-0.1	MAIN	40	41	HIGH	3.000	SUB	29	30				START	END	LIMIT (%)	FIRST				SECOND			<input checked="" type="checkbox"/> Low High	THIRD			<input checked="" type="checkbox"/> Low High
READ	START	END	LOW	-0.1																												
MAIN	40	41	HIGH	3.000																												
SUB	29	30																														
	START	END	LIMIT (%)																													
FIRST																																
SECOND			<input checked="" type="checkbox"/> Low High																													
THIRD			<input checked="" type="checkbox"/> Low High																													

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days.

PERFORMANCE CHARACTERISTICS

1. **Measurement range:** Up to 500 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. **Detection Limit:** Values less than 16 mg/dL give non-reproducible results.
3. **Prozone effect:** No prozone effect was detected upon 1000 mg/dL
4. **Sensitivity:** Δ 3,4 mA / mg/dL.
5. **Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	34,39 mg/dL	92,7 mg/dL	181,8 mg/dL
Total	4,2%	2,6%	2,8%
Within Run	0,8%	1,1%	1,6%
Between Run	3,8%	2,4%	2,3%
Between Day	1,6%	0%	0%

6. **Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the Beckman Array 360 CE. 35 samples ranging from 70 to 250 mg/dL of α₁antitrypsin were assayed. The correlation coefficient (r) was 0.92 and the regression equation y = 0,84x +26,5.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Sharp HL. Hospital Practice; May 1971: 83-96
5. Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

PACKAGING

Ref: TK1102054 Cont. R1.Diluent: 2 x 24 mL
R2.Antibody: 2 x 6 mL

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Sharp HL. Hospital Practice; May 1971: 83-96
- Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
- Young DS.Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clin. laboratory tests, 3tn ed. AACC Pres, 1997.

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BIBLIOGRAPHIE

1. Clinical Guide to Laboratory Tests, Edited by NWTietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. PesceAJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Sharp HL. Hospital Practice; May 1971: 83-96
5. Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

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Ref: TK1102054	Cont.	R1.Diluant:	2 x 24 mL
		R2.Anticorp:	2 x 6 mL

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