# Ceruloplasmin

## Turbidimetry

### Quantitative determination of Ceruloplasmin IVD

Store at 2-8°C.

#### INTENDED USE

The Ceruloplasmin reagent is a quantitative turbidimetric test for the measurement of Ceruloplasmin in human serum or plasma.

#### PRINCIPLE OF THE METHOD

Anti-human Ceruloplasmin antibodies when mixed with samples containing Ceruloplasmin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the Ceruloplasmin concentration of the patient sample, that can be quantified by comparison from a calibrator of known Ceruloplasmin concentration.

#### CLINICAL SIGNIFICANCE

Caeruloplasmin is an  $\alpha_2$ -globulin that contains approximately 95% of total serum copper. Each molecule of caeruloplasmin contain six to eight copper atoms. The high content of copper ions gives caeruloplasmin a blue color. Caeruloplasmin can also bind, and possible transport, other cations such as magnesium. The molecule of caeruloplasmin has a single polypeptide chain and carbohydrate, and results a molecular mass of 132 kD. Caeruloplasmin is synthesized primarily by the hepatic cells and small quantities by macrophages and lymphocytes.

Caeruloplasmin is most often quantified as a screening test for Wilson's disease. However, it is important to realize that several other factors, including diet, hormone levels, and other genetic disorders, can influence plasma levels.

Synthesis of caeruloplasmin is increased modestly in the acute-phase response, peaking at 4 to 20 days after a single, acute insult. Synthesis is also stimulated by estrogens, and during pregnancy.

Low plasma caeruloplasmin levels are due to lack of incorporation of  $Cu^{2+}$  into the molecule during synthesis. The causes are, the dietary insufficiency (including malabsorption), inability to release  $Cu^{2+}$  from gastrointestinal epithelium into circulation, or inability to insert  $Cu^{2+}$  into developing caeruloplasmin molecule. Levels may also be low in blood loss or gastrointestinal or renal proteinlosing syndromes.

#### REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
Antibody (R2)	Goat serum, anti-human Ceruloplasmin, pH 7.5. Preservative.
Optional	Cod: 1102003 PROT CAL.

#### CALIBRATION

The assay is calibrated to the Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements). It must be used the PROT CAL to calibrate the reagent.

#### PREPARATION

**Reagents:** Ready to use.

**Calibration Curve:** Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the Ceruloplasmin calibrator by the corresponding factor stated in table below to obtain the Ceruloplasmin concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator ( $\mu$ L)	--	10	25	50	75	100
NaCl 9 g/L ( $\mu$ L)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not use. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

#### ADDITIONAL EQUIPMENT

- Laboratory equipment.
- Spintech 240 autoanalyzer

#### SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

#### REFERENCE VALUES

Between 15 – 60 mg/dL. Each laboratory should establish its own reference range.

#### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

Spinreact has instructions sheets available for several automatic analyzers.

### SPINTECH 240 APPLICATION

Item Name CERUL			
<u>DATA INFORMATION</u>		<u>CALIBRATION</u>	
Units	mg/dL	TYPE	Spline
Decimals	0		
<u>ANALYSIS</u>		STANDARD	
Type	END	#1	0.10 x Cal. Val.
		#2	0.25 x Cal. Val.
		#3	0.50 x Cal. Val.
		#4	0.75 x Cal. Val.
		#5	1.00 x Cal. Val.
		#6	
W.Length 1	340	<u>NORMAL RANGE</u>	
			LOW HIGH
Method	Turbidimetry	SERUM	MALE FEMALE
		URINE	
<u>CORR</u>			
SLOPE	INTER		
1.000 x +	0		
Item Name CERUL			
<u>ASPIRATION</u>		<u>DATA PROCESS</u>	
KIND	Single <input checked="" type="checkbox"/> Double	<u>READ</u>	<u>ABSORBANCE LIMIT</u>
		START	END
		MAIN	40 41
		SUB	28 29
VOLUME**			
SAMPLE	3 $\mu$ L		
REAGENT 1	240 $\mu$ L		
REAGENT 2	60		
			ENDPOINT LIMIT 3
			LINEAR CHECK (%)
		<u>FACTOR</u>	
		Blank Correction	1.000
		<u>PROZONE CHECK</u>	
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON		
R1 Blank	Water <input checked="" type="checkbox"/> R1-B		
<u>MONITOR</u>			
0 LEVEL POINT	1		
SPAN	3.000		
		START	END
		FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

\*\* Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days

#### PERFORMANCE CHARACTERISTICS

**Measurement range:** Up to 91 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

**Detection Limit:** Values less than 1,12 mg/dL give non-reproducible results.

**Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study (NCCLS).

EP5	CV (%)		
	28.96mg/dL	55.47 mg/dL	76.54 mg/dL
Total	4%	2.3%	2%
Within Run	2.2%	1.5%	1%
Between Run	3.1%	1.1%	1.5%
Between Day	1.1%	1.3%	0.8%

**Accuracy:** Results obtained using this reagent (y) were compared to those obtained with a Bayer immunoturbidimetric method. 45 samples ranging from 20 to 80 mg/dL of Ceruloplasmin were assayed. The correlation coefficient (r) was 0.96 and the regression equation  $y = 0.896x + 10.57$ .

The results of the performance characteristics depend on the used analyzer.

#### INTERFERENCES

Hemoglobin (16 g/L), bilirubin (40 mg/dL), lipemia (< 2.5 g/L), and rheumatoid factor (800 IU/mL) do not interfere. Other substances may interfere.

#### BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
3. Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
4. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
5. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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### BIBLIOGRAPHIE

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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### BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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Ref.: TK1102064

Cont.

R1. Solvente: 2 x 24 mL  
R2. Anticorpo: 2 x 6 mL