Lp(a)-turbilatex

Latex turbidimetry

Quantitative determination of Lipoprotein (a) (Lp(a)) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The Lp(a)-turbilatex is a quantitative turbidimetric test for the measurement of Lp(a) in human serum or plasma.

Latex particles coated with antibodies anti-Lp(a) are agglutinated when mixed with samples containing Lp(a). The agglutination causes an absorbance change, dependent upon the Lp(a) contents of sample that can be quantified by comparison from a calibrator of known Lp(a) concentration.

CLINICAL SIGNIFICANCE

Lp(a) is a low density lipoprotein-like particle containing apolipoprotein B-100 disulphide-linked to one large glycoprotein called apolipoprotein (a). Many investigators have confirmed that a high Lp(a) concentration represents an indicator of risk for cardiovascular disease, especially when serum LDL-cholesterol or Apo B are elevated. The quantification of Lp(a) in serum or plasma is important for identification of individuals at risk for developing atherosclerosis.

REAGENTS

Diluent (R1)	Glycine buffer 50 mmol/L. Preservative.
Latex (R2)	Latex particles coated with rabbit polyclonal anti-human Lp(a). Preservative.
Optional	Ref.: 1107022 Lp(a) Calibrator. Ref.: 1107024 Lp(a) Control.

CALIBRATION

The sensitivity of the assay and the target value of the calibrator have been standardized against an Internal Reference Material. It is not recommended the use of other commercially available Lp(a) calibrators.

PREPARATION

Calibration Curve Prepare the following Lp(a) calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the Lp(a) calibrator by the corresponding factor stated in table below to obtain the Lp(a) concentration of each dilution.

Calibrator dilution	1	2	3	4	5
Lp(a) Calibrator (µL)	--	25	50	75	100
NaCl 9 g/L (µL)	100	75	50	25	-
Factor	0	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen latex and diluent could change the functionality of the test.

Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1).

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Normal values up to 30 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Lp(a) Control Ref.: 1107024.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

INTERFERENCES

Hemoglobin (1000 mg/dL), bilirubin (60 mg/dL) and lipemia (10 g/L) do not interfere. Other substances may interfere⁹.

SPINTECH 240 APPLICATION

Item Name Lp (a)		CALIBRATION	
DATA INFORMATION		TYPE	
Units	mg/dL	Spline	
Decimals	1		
ANALYSIS		STANDARD	
Type	END	#1 *	#4 *
W.Length 1	570	#2 *	#5 *
		#3 *	#6 *
Method		NORMAL RANGE (37°C)	
Method	Turbilatex	SERUM MALE	LOW HIGH
		FEMALE	0 30
		URINE (mg/24h)	0 30
CORR		ASPIRATION	
SLOPE	INTER	DATA PROCESS	
1.000 x +	0	ABSORBANCE LIMIT	
ASPIRATION		READ	
KIND	Single <input checked="" type="checkbox"/> Double	START	END
		MAIN 47	48
		SUB 35	36
VOLUME**		ENDPOINT LIMIT 3	
SAMPLE	4,5 µL	LINEAR CHECK (%)	
REAGENT 1	240 µL		
REAGENT 2	60 µL		
Third Mix		FACTOR	
	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	Blank Correction 1.000	
R1 Blank	<input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B	PROZONE CHECK	
MONITOR		START END LIMIT (%)	
0 LEVEL POINT	1	FIRST	
SPAN	3.000	SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days

PERFORMANCE CHARACTERISTICS

- Linearity:** Up to 100 mg/dL, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit and measurement range depends on the sample to reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Quantification limit:** Values less than 3,17 mg/dL give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 650 mg/dL.
- Precision:**

	Level 1	Level 2	Level 3
Mean (mg/dL)	18,81	34,12	47,60
SD	0,59	0,48	0,33
CV	3,15	1,42	0,69
- Accuracy:** The Spinreact method was compared with another manufacturer. The study was performed with 35 serum samples. Both tests were performed on a Spintech240. Both tests were calibrated with their respective calibrators. The correlation coefficient (r) was 0,9815 and the regression equation $y = 1,5734x + 1,5873$.

The results of the performance characteristics depend on the analyzer used.

NOTES

- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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Item Name Lp (a) DATA INFORMATION Units mg/dL Decimals 1 ANALYSIS Type END W.Length 1 570 Method Turbidilatex CORR SLOPE INTER 1.000 x + 0		CALIBRATION TYPE Spline STANDARD #1 * #4 * #2 * #5 #3 * #6 NORMAL RANGE LOW HIGH SERUM MALE 0 30 FEMALE 0 30 URINE (mg/24h)	
Item Name Lp (a) ASPIRATION KIND Single <input checked="" type="checkbox"/> Double VOLUME** SAMPLE 4.5 µL REAGENT 1 240 µL REAGENT 2 60 µL Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank <input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B MONITOR 0 LEVEL POINT 1 SPAN 3.000		DATA PROCESS ABSORBANCE LIMIT READ START END LOW HIGH MAIN 47 48 -3.000 3.000 SUB 35 36 ENDPOINT LIMIT 3 LINEAR CHECK (%) FACTOR Blank Correction 1.000 PROZONE CHECK START END LIMIT (%) FIRST <input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High SECOND <input checked="" type="checkbox"/> Low <input checked="" type="checkbox"/> High THIRD	

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<p>Item Name Lp (a)</p> <p>DATA INFORMATION</p> <p>Units mg/dL</p> <p>Decimals 1</p> <p>ANALYSIS</p> <p>Type END</p> <p>W.Length 1 570</p> <p>Method Turbidilatex</p> <p>CORR</p> <p>SLOPE INTER</p> <p>1.000 x + 0</p>		<p>CALIBRATION</p> <p>TYPE Spline</p> <p>STANDARD</p> <p>#1 * #4 *</p> <p>#2 * #5</p> <p>#3 * #6</p> <p>NORMAL RANGE (37°C)</p> <p>SERUM MALE 0 30</p> <p>FEMALE 0 30</p> <p>URINE (mg/24h)</p>	
<p>Item Name Lp (a)</p> <p>ASPIRATION</p> <p>KIND Single <input checked="" type="checkbox"/> Double</p> <p>VOLUME**</p> <p>SAMPLE 4,5 µL</p> <p>REAGENT 1 240 µL</p> <p>REAGENT 2 60 µL</p> <p>Third Mix <input checked="" type="checkbox"/> OFF ON</p> <p>R1 Blank <input checked="" type="checkbox"/> Water R1-B</p> <p>MONITOR</p> <p>0 LEVEL POINT 1</p> <p>SPAN 3.000</p>		<p>DATA PROCESS</p> <p>ABSORBANCE LIMIT</p> <p>READ</p> <p>START END LOW HIGH</p> <p>MAIN 47 48 0 3.000</p> <p>SUB 35 36</p> <p>ENDPOINT LIMIT 3</p> <p>LINEAR CHECK (%)</p> <p>FACTOR</p> <p>Blank Correction 1.000</p> <p>PROZONE CHECK</p> <p>START END LIMIT (%)</p> <p>FIRST</p> <p>SECOND <input checked="" type="checkbox"/> Low High</p> <p>THIRD <input checked="" type="checkbox"/> Low High</p>	

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