**Quantitative determination of immunoglobulin E (IgE)**
**IVD**

Store 2 - 8°C.

**PRINCIPLE OF THE METHOD**

The IgE-Turbilatex is a quantitative turbidimetric test for the measurement of IgE in human serum or plasma.

Latex particles coated with mouse IgG anti-human IgE are agglutinated when mixed with samples containing IgE. The agglutination causes an absorbance change, dependent upon the IgE contents of the patient sample that can be quantified by comparison from a calibrator of known IgE concentration.

**CLINICAL SIGNIFICANCE**

IgE is an immunoglobulin with a molecular weight of approximately 190,000 normally present in trace amounts. Continual production of IgE antibodies in response to common naturally occurring allergens, however, often results in elevated serum levels and in the development of such clinically important Type I allergic reactions as asthma, hay fever, dermatitis and food allergies. Elevated IgE levels are also seen in parasitic diseases, IgE myeloma, and in hepatitis. The measurement of IgE in human serum is thus considered to be useful in the diagnosis, treatment assessment of disease progression, or postoperative prognosis for such conditions.

**REAGENTS**

<b>Diluent (R1)</b>	Glycine buffer, pH 8.3. Preservative.
<b>Latex (R2)</b>	Latex particles coated with mouse IgG anti-human IgE, pH 7.3. Preservative.
<b>Optional</b>	Ref: 1107052 IgE CAL (Calibrator) Ref: 1107053 IgE Control.

**CALIBRATION**

The sensitivity of the assay and the target value of the calibrator have been standardized against the International Reference Preparation of IgE 75/502 (2<sup>nd</sup> IRP 1981) from WHO. It is not recommended the use of other commercially available IgE calibrators.

**PREPARATION**

**Reagents:** Ready to use.

**Calibration curve:** Prepare the following IgE calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the IgE calibrator by the corresponding factor stated in table below to obtain the IgE concentration of each dilution.

Calibrator dilution	1	2	3	4	5
IgE Calibrator (µL)	--	12.5	25	50	100
NaCl 9 g/L (µL)	100	87.5	75	50	--
Factor	0	0.125	0.25	0.5	1.0

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Frozen Latex and Diluent could change the functionality of the test.

**Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

**ADDITIONAL EQUIPMENT**

- Spintech 240 autoanalyzer
- Laboratory equipment.

**SAMPLES**

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

**REFERENCE VALUES**

Up to 350 IU/mL. Each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control Sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact Control Serum IgE is available (Ref.: 1107053). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**NOTES**

Clinical diagnosis should be not made on findings of a single test result but should integrate both clinical and laboratory data.

**SPINTECH 240 APPLICATION**

Item Name IgE <b>DATA INFORMATION</b> Units UI/mL Decimals 0 <b>ANALYSIS</b> Type END W.Length 1 570 Method Turbilatex <b>CORR</b> SLOPE INTER 1.000 x + 0		<b>CALIBRATION</b> TYPE Spline STANDARD #1 0.125 x Cal.Val. #4 1.00 x Cal.Val. #2 0.250 x Cal.Val #5 #3 0.500 x Cal.Val. #6 <b>NORMAL RANGE</b> SERUM MALE LOW HIGH FEMALE	
Item Name IgE <b>ASPIRATION</b> KIND Single <input checked="" type="checkbox"/> Double VOLUME** SAMPLE 4 µL REAGENT 1 175 µL REAGENT 2 95 µL Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B		<b>DATA PROCESS</b> <b>ABSORBANCE LIMIT</b> <b>READ</b> START END LOW -3.000 MAIN 51 52 HIGH 3.000 SUB 34 35 <b>ENDPOINT LIMIT 2</b> <b>LINEAR CHECK (%)</b> <b>FACTOR</b> Blank Correction 1.000 <b>PROZONE CHECK</b> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	
<b>MONITOR</b> 0 LEVEL POINT 1 SPAN 3.000			

**\*\* Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio. Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days**

**PERFORMANCE CHARACTERISTICS**

- 1. Measurement range:** Up to 1000 IU/mL, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit and measurement range depend on the sample to reagent / ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- 2. Limit detection:** Values less than 3,04 IU/mL give non-reproducible results.
- 3. Prozone effect:** No prozone effect was detected upon 60.000 IU/mL.
- 4. Sensitivity:** Δ 0,5 mA. IU/mL.
- 5. Precision:**

	Within run		
	Mean (IU/mL)	SD (IU/mL)	CV (%)
n= 21	64,0	1,07	1,68
	121,0	1,43	1,18
	295,4	1,20	0,41
	491,3	4,01	0,82
	739,8	2,63	0,36

	Between days		
	Mean (IU/mL)	SD (IU/mL)	CV (%)
n= 21	61,0	2,9	4,76
	406	7,8	1,92

- 6. Accuracy:** Results obtained using this reagents (y) were compared to those obtained using a nephelometric method (x). 92 samples were assayed by both methods. The correlation coefficient (r) was 0,9748 and the regression line equation  $y = 1,0321x + 8,7369$ .

The results of the performance characteristics depend on the analyzer used.

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