Quantitative determination of C-Reactive Protein (CRP)
IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

CRP-Turbilatex is a quantitative turbidimetric test for the measurement of C-reactive protein (CRP) in human serum or plasma.

Latex particles coated with specific anti-human CRP are agglutinated when mixed with samples containing CRP. The agglutination causes an absorbance change, dependent upon the CRP contents of the patient sample that can be quantified by comparison from a calibrator of known CRP concentration.

CLINICAL SIGNIFICANCE

CRP is an acute-phase protein present in normal serum, which increases significantly after most forms of tissue injuries, bacterial and virus infections, inflammation and malignant neoplasia. During tissue necrosis and inflammation resulting from microbial infections, the CRP concentration can rise up to 300 mg/L in 12-24 hours.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, pH 8.2. Preservative.
Latex (R2)	Latex particles coated with goat IgG anti-human CRP, pH 7.3. Preservative.
CRP-CAL	Calibrator. C-Reactive protein concentration is stated on the vial label.
Optional	Ref.:1102114 Control serum ASO/CRP/RF Level L Ref.:110211 Control serum ASO/CRP/RF Level H.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Use CRP Calibrator Reference 1107002.

The sensitivity of the assay and the target value of the calibrator have been standardized against the Reference Material ERM-DA 474/IFCC.

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready to use.

CRP Calibrator: Reconstitute (→) with 1.0 mL of distilled water. Mix gently and incubate 10 minutes at room temperature before use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1).

CRP Calibrator: Stable for 1 month at 2-8°C or 3 months at -20°C.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Normal values up to 6 mg/L.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used SPINREACT Controls ASO/CRP/RF Level L (Ref.:1102114) and Level H (Ref.: 1102115).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

INTERFERENCES

Bilirubin (20 mg/dL), lipemia (10 g/L) and rheumatoid factors (300 IU/mL) do not interfere. Hemoglobin (≥ 5 g/L), interferes. Other substances may interfere.

SPINTECH 240 APPLICATION

Item Name	CRP		
<u>DATA INFORMATION</u>		<u>CALIBRATION</u>	
Units	mg/L	TYPE	Linear
Decimals	0		
<u>ANALYSIS</u>		STANDARD	
Type	END	#1 *	#4
		#2	#5
		#3	#6
W.Length 1	546	<u>NORMAL RANGE</u>	
		LOW	HIGH
Method	Turbilatex	SERUM	MALE
			FEMALE
<u>CORR</u>		URINE	
SLOPE	INTER		
1.000 x +	0		
Item Name	CRP		
<u>ASPIRATION</u>		<u>DATA PROCESS</u>	<u>ABSORBANCE LIMIT</u>
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	<u>READ</u>	LOW -3.000
		START END	HIGH 3.000
	VOLUME**	MAIN 41 42	
SAMPLE	3 μ L	SUB 33 34	
REAGENT 1	240 μ L		ENDPOINT LIMIT 3
REAGENT 2	60 μ L		LINEAR CHECK (%)
		<u>FACTOR</u>	
Third Mix	OF <input type="checkbox"/> ON <input checked="" type="checkbox"/>	Blank Correction	1.000
R1 Blank	Water <input type="checkbox"/> R1-B <input checked="" type="checkbox"/>		
<u>MONITOR</u>		<u>PROZONE CHECK</u>	
0 LEVEL POINT	1	START END LIMIT (%)	
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days

PERFORMANCE CHARACTERISTICS

1. Linearity limit: Up to 150 mg/L, under the described assay conditions.

Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

2. Detection limit: Values less than 1 mg/L give non-reproducible results.

3. Prozone effect: No prozone effect was detected upon 800 mg/L.

4. Sensitivity: Δ 4,2 mA.mg/L.

5. Precision: The reagent has been tested for 20 days, using three different CRP concentrations in a EP5-based study.

EP5	CV (%)		
	9,2 mg/L	16,8 mg/L	57,97 mg/L
Total	7,3%	6,9%	5,9%
Within Run	2,8%	3,1%	2,9%
Between Run	6,1%	4,7%	3,9%
Between Day	3,0%	4,0%	3,4%

6. Accuracy: Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 50 samples of different concentrations of CRP were assayed. The correlation coefficient (r^2) was 0.99 and the regression equation $y = 1,101x + 2,518$.

The results of the performance characteristics depend on the analyzer used.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

- Lars-Olof Hanson et al. Current Opinion in Infect Diseases 1997; 10: 196-201.
- Chetana Vaishnavi. Immunology and Infectious Diseases 1996; 6: 139 - 144.
- Yoshitsugy Hokama et al. Journal of Clinical Lab. Status 1987; 1: 15 - 27.
- Kari Pulki et al. Sacand J Clin Lab Invest 1986; 46: 606 - 607.
- Werner Müller et al. Journal of Immunological Methods 1985; 80: 77 - 90.
- Shogo Otsuji et al. Clin Chem 1982; 28/10: 2121 - 2124.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AAC Press, 1995.

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1. Lars-Olof Hanson et al. Current Opinion in Infectious diseases 1997; 10: 196-201.
2. Chetana Vaishnavi. Immunology and Infectious Diseases 1996; 6: 139 - 144.
3. Yoshitsugy Hokama et al. Journal of Clinical Laboratory Status 1987; 1: 15 - 27.
4. Kari Pulki et al. Sacand J Clin Lab Invest 1986; 46: 606 - 607.
5. Werner Müller et al. Journal of Immunological Methods 1985; 80: 77 - 90.
6. Shogo Otsuji et al. Clin Chem 1982; 28/10: 2121 - 2124.
7. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.

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APPLICATION AU SPINTECH 240

Item Name CRP DATA INFORMATION Units mg/L Decimals 0 ANALYSIS Type END W.Length 1 546 Method Turbilatex CORR SLOPE INTER 1.000 x + 0		CALIBRATION TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 NORMAL RANGE LOW HIGH SERUM MALE FEMALE URINE	
Item Name CRP ASPIRATION KIND Single <input checked="" type="checkbox"/> Double VOLUME** SAMPLE 3 µL REAGENT 1 240 µL REAGENT 2 60 µL Third Mix OF <input checked="" type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B MONITOR 0 LEVEL POINT 1 SPAN 3.000		DATA PROCESS ABSORBANCE LIMIT READ START END MAIN 41 42 SUB 33 34 ENDPOINT LIMIT 3 LINEAR CHECK (%) FACTOR Blank Correction 1.000 PROZONE CHECK START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

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BIBLIOGRAPHIE

- Lars-Olof Hanson et al. Current Opinion in Infectious diseases 1997; 10: 196-201.
- Chetana Vaishnavi. Immunology and Infectious Diseases 1996;6:139- 144.
- Yoshitsugu Hokama et al. Journal of Clinical Laboratory Status 1987;1: 5-27.
- Kari Pulki et al. Sacand J Clin Lab Invest 1986; 46: 606 - 607.
- Werner Müller et al. Journal of Immunological Methods 1985; 80: 77 - 90.
- Shogo Otsuji et al. Clin Chem 1982; 28/10: 2121 - 2124.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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Item Name CRP DATA INFORMATION Units mg/L Decimals 0 ANALYSIS Type END W.Length 1 546 Method Turbilatex CORR SLOPE INTER 1.000 x + 0		CALIBRATION TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 NORMAL RANGE LOW HIGH SERUM MALE FEMALE URINE	
Item Name CRP ASPIRATION KIND Single <input checked="" type="checkbox"/> Double VOLUME** SAMPLE 3 µL REAGENT 1 240 µL REAGENT 2 60 µL Third Mix OF <input checked="" type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B MONITOR 0 LEVEL POINT 1 SPAN 3.000		DATA PROCESS ABSORBANCE LIMIT READ START END LOW -3.000 HIGH 3.000 MAIN 41 42 SUB 33 34 ENDPOINT LIMIT 3 LINEAR CHECK (%) FACTOR Blank Correction 1.000 PROZONE CHECK START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

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BIBLIOGRAFIA

1. Lars-Olof Hanson et al. Current Opinion in Infectious diseases 1997;10:196-201.
2. Chetana Vaishnavi. Immunology and Infectious Diseases 1996; 6: 139 – 144.
3. Yoshitsugy Hokama et al. Journal of Clinical Laboratory Status 1987; 1: 15 – 27.
4. Kari Pulki et al. Sacand J Clin Lab Invest 1986; 46: 606 – 607.
5. Werner Müller et al. Journal of Immunological Methods 1985; 80: 77 – 90.
6. Shogo Otsuji et al. Clin Chem 1982; 28/10: 2121 – 2124.
7. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCPress, 1995.

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Ref: TK1107001

Cont.

 R1. Diluíente: 2 x 24 mL
 R2. Latex: 1 x 12 mL
 CRP-CAL: 1 x 1 mL