Quantitative determination of anti-streptolysin O (ASO) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The ASO-Turbilatex is a quantitative turbidimetric test for the measurement of ASO in human serum or plasma.

Latex particles coated with recombinant streptolysin O (SLO) are agglutinated when mixed with samples containing ASO. The agglutination causes an absorbance change, dependent upon the ASO contents of the patient sample that can be quantified by comparison from a calibrator of known ASO concentration.

CLINICAL SIGNIFICANCE

SLO is a toxic immunogenic exoenzyme produced by β -hemolytic Streptococci of groups A, C and G. Measuring the ASO antibodies are useful for the diagnostic of rheumatoid fever, acute glomerulonephritis and streptococcal infections. Rheumatic fever is an inflammatory disease affecting connective tissue from several parts of human body as skin, heart, joints etc... and acute glomerulonephritis is a renal infection that affects mainly to renal glomerulus.

REAGENTS

| | |
|---------------------|--|
| Diluent (R1) | Tris buffer 20 mmol/L, pH 8.2. Preservative. |
| Latex (R2) | Latex particles coated with recombinant Streptolysin-O, pH 8.2. Preservative. |
| ASO-CAL | Calibrator. Human serum. ASO concentration is stated on the vial label. |
| Optional | Ref.:1102114 Control ASO/CRP/RF Level L Ref.:1102115 Control ASO/CRP/RF Level H |

PRECAUTIONS

R1: H317 - May cause an allergic skin reaction.

Contains 2-Methylisothiazol-3(2H)-one (Proclin 950).

R2: H360 – May damage fertility or the unborn child. Contains boric acid (H3BO3)

Follow the precautionary statements given in MSDS and label of the product.

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Use ASO Calibrator Reference 1107012.

The sensitivity of the assay and the target value of the calibrator have been standardized against the ASO International Calibrator from NIBSC ASO.

The calibration in the SPINLAB 180 is stable for 3 weeks.

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready to use.

ASO Calibrator: Reconstitute (→) with 1.0 mL of distilled water. Mix gently and incubate at room temperature for 10 minutes before use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Frozen Latex and Diluent could change the functionality of the test.

Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1).

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Normal values up to 200 IU/mL (adults) and 100 IU/mL (children < 5 years old)⁶.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Controls ASO/CRP/RF Level L (Ref.: 1102114) and Level H (Ref.: 1102115).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

INTERFERENCES

Bilirubin (20 mg/dL), hemoglobin (10 g/L), lipemia (10 g/L) and rheumatoid factors (600 IU/mL), do not interfere. Other substances may interfere⁷.

SPINTECH 240 APPLICATION

| | | | |
|-------------------------|--|-------------------------|--|
| Item Name ASO | | CALIBRATION | |
| <u>DATA INFORMATION</u> | | TYPE Linear | |
| Units | IU/mL | | |
| Decimals | 0 | | |
| <u>ANALYSIS</u> | | STANDARD | |
| Type | END | #1 * | #4 |
| | | #2 | #5 |
| | | #3 | #6 |
| W.Length 1 | 546 | <u>NORMAL RANGE</u> | |
| | | LOW | HIGH |
| Method | Turbilatex | SERUM MALE | |
| | | FEMALE | |
| | | URINE | |
| CORR SLOPE | INTER | | |
| 1.000 x + | 0 | | |
| Item Name ASO | | DATA PROCESS | |
| <u>ASPIRATION</u> | | <u>ABSORBANCE LIMIT</u> | |
| KIND | Single <input type="checkbox"/> Double <input checked="" type="checkbox"/> | READ | LOW -3.000 |
| | | START END | HIGH 3.000 |
| SAMPLE | VOLUME** | MAIN 41 42 | |
| REAGENT 1 | 3 μ L | SUB 33 34 | |
| REAGENT 2 | 240 μ L | | ENDPOINT LIMIT 3 |
| REAGENT 2 | 60 μ L | | LINEAR CHECK (%) |
| Third Mix | <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON | <u>FACTOR</u> | |
| R1 Blank | Water <input checked="" type="checkbox"/> R1-B | Blank Correction | 1.000 |
| <u>MONITOR</u> | | <u>PROZONE CHECK</u> | |
| 0 LEVEL POINT | 1 | START END LIMIT (%) | |
| SPAN | 3.000 | FIRST | |
| | | SECOND | <input checked="" type="checkbox"/> Low High |
| | | THIRD | <input checked="" type="checkbox"/> Low High |

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days

PERFORMANCE CHARACTERISTICS

1. Linearity limit: Up to 800 IU/mL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/3 in NaCl 9 g/L and retested again. The linearity limit depends on the sample-reagent ratio, as well the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

2. Detection limit: Values less than 20 IU/mL give non-reproducible results.

3. Prozone effect: No prozone effect was detected up to 1000 IU/mL.

4. Sensitivity: Δ 0,73 mA. IU/mL.

5. Precision: The reagent has been tested for 20 days, using three different ASO concentrations in a EP5-based study.

| EP5 | CV (%) | | |
|-------------|---------------|---------------|---------------|
| | +/- 100 IU/mL | +/- 200 IU/mL | +/- 400 IU/mL |
| Total | 6,4% | 5,7% | 5,1% |
| Within Run | 2,4% | 1,7% | 1,4% |
| Between Run | 3,6% | 4,2% | 4,9% |
| Between Day | 4,7% | 3,5% | 0,7% |

6. Accuracy: Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 60 samples of different concentrations of ASO were assayed. The correlation coefficient (r) was 0.99 and the regression equation $y = 0,915x - 4,844$.

The results of the performance characteristics depend on the analyzer used.

NOTES

Clinical diagnosis should be not made on findings of a single test result, but should integrate both clinical and laboratory data.

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