Quantitative determination of Ferritin IVD

Store at 2 - 8°C.

RECOMMENDED USE

Turbidimetric immunoassay for the quantitative determination of ferritin in human serum or plasma.

PRINCIPLE OF THE METHOD

Ferritin-turbilatex is a quantitative turbidimetric test for the measurement of ferritin in human serum or plasma.

Latex particles coated with specific anti-human ferritin are agglutinated when mixed with samples containing ferritin. The agglutination causes an absorbance change, dependent upon the ferritin contents of the sample that can be quantified by comparison from a calibrator of known ferritin concentration.

CLINICAL SIGNIFICANCE

Serum ferritin concentration usually reflects body iron stores and is considered one of the most reliable indicators of iron status of patients. Whereas low serum concentrations of ferritin are always indicative of an iron deficiency, elevated concentrations can occur for variety of reasons. Thus, although elevated concentrations often indicate an excessive iron intake, they are also caused by liver disease, chronic inflammation and malignancies. Pregnant women, blood donors, hemodialysis patients, adolescents and children are groups particularly at risk.

REAGENTS

Diluent (R1)	Tris Buffer 20 mmol/L, pH 8,2. Preservative
Latex (R2)	Latex particles coated with rabbit IgG anti-human ferritin pH 8,2. Preservative.
FERR-CAL	Calibrator. Ferritin concentration is stated on the vial.
Opcional	Ref: 1107044 Ferritin Control.

PRECAUTIONS

R1 and R2: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary statements given in MSDS and label of the product.

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

PREPARATION

Ferritin Calibrator: Reconstitute (→) with 3.0 mL of distilled water. Mix gently and incubate at room temperature for about 10 minutes before testing.

CALIBRATION

Use Ferritin Calibrator Reference 1107042.

The sensitivity of the assay and the target value of the calibrator have been standardized against the 3rd International Standard of Ferritin (94/572, 2008 WHO). Recalibrate when control results are out of specified values; when using a different lot of reagent and when the instrument is adjusted.

Calibration curve: Prepare the following dilutions of the FERR Calibrator using NaCl 9 g/L. To obtain the concentration of each dilution, multiply using the dilution factor shown in the next table:

Calibrator dilution	1	2	3	4	5	6
Calibrator FERR (µL)	--	25	50	100	200	400
NaCl 9 g/L (µL)	400	375	350	300	200	--
Dilution Factor	0	1/16	1/8	1/4	1/2	1,0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1).

ADDITIONAL EQUIPMENT

- Laboratory equipment
- Spintech 240 analyzer.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Men: 30 - 220 µg/L.

Women: 20 - 110 µg/L.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Ferritin Control (Ref.: 1107044).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

SPINTECH 240 APPLICATION

Item Name FERR DATA INFORMATION Units ug/L Decimals 0 ANALYSIS Type END W.Length 1 546 Method CORR SLOPE INTER 1.000 x + 0		CALIBRATION TYPE Spline STANDARD #1 0,0625 x Cal. Val #4 0,50 x Cal. Val #2 0,125 x Cal. Val #5 1,00 x Cal. Val #3 0,250 x Cal. Val. #6 NORMAL RANGE LOW HIGH SERUM MALE FEMALE URINE	
Item Name FERR ASPIRATION KIND Single <input checked="" type="checkbox"/> Double VOLUME** SAMPLE 18 µL REAGENT 1 160 µL REAGENT 2 40 µL Third Mix <input checked="" type="checkbox"/> OF ON R1 Blank <input checked="" type="checkbox"/> Water R1-B		DATA PROCESS ABSORBANCE LIMIT READ LOW -3.000 START END HIGH 3.000 MAIN 51 52 SUB 33 34 ENDPOINT LIMIT 3 LINEAR CHECK (%) FACTOR Blank Correction 1.000 PROZONE CHECK 0 LEVEL POINT 1 START END LIMIT (%) SPAN 3.000 FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 30 days

INTERFERENCES

Bilirubin (40 mg/dL), hemoglobin (5 g/L), y and rheumatoid factor (750 UI/mL), do not interfere. Lipids (≥ 2,5 g/L) do interfere. Other substances may interfere ⁵.

NOTE

Clinical diagnosis should not be based on findings of a single test result, but should integrate both clinical and laboratory data.

PERFORMANCE CHARACTERISTICS

Measuring range: Up to 600 µg/L. Samples with higher values should be diluted 1/5 in NaCl 9 g/L and retested. The upper linearity limit increases as the sample volume and the sensitivity decrease.

Detection limit: 5,04 µg/L.

Quantification limit: Values under 6,6 µg/L may give non-reproducible results.

Prozone effect: No prozone effect was detected at least up to 9000 µg/L.

Precision: According to the EP5-A2 standards (CLSI), the reagent has been tested for 20 days, measuring each level per duplicate twice a day (n=80):

Mean (µg/L)	Intra-assay (n= 80)			Total (n= 80)		
	33,4	114,5	289,8	33,4	114,5	289,8
SD	1,7	1,4	2,4	2,1	3,4	7,5
CV (%)	5,1	1,2	0,8	6,3	2,9	2,6

Method comparison: The reagent was compared to another commercially available Ferritin reagent by testing 144 samples (male and female), with concentrations between 6,97 and 730 µg/L. The coefficient of correlation (r) was 0,988, and the equation $y = 0,96x + 1,15$

Performance characteristics depend on the analyzer used.

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