Quantitative determination of glycated hemoglobin (HbA_{1c}) in human blood IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

This method utilizes the interaction of antigen and antibody to directly determine the HbA_{1c} in whole blood. Total hemoglobin and HbA_{1c} have the same unspecific absorption rate to latex particles. When mouse anti-human HbA_{1c} monoclonal antibody is added (R2), latex-HbA_{1c}-mouse anti human HbA_{1c} antibody complex is formed. Agglutination is formed when goat anti-mouse IgG polyclonal antibody interacts with the monoclonal antibody. The amount of agglutination is proportional to the amount of HbA_{1c} absorbed on to the surface of latex particles. The amount of agglutination is measured as absorbance. The HbA_{1c} value is obtained from a calibration curve.

CLINICAL SIGNIFICANCE

Throughout the circulatory life of the red cell, Hemoglobin A_{1c} is formed continuously by the addition of glucose to the N-terminal of the hemoglobin beta chain. This process, which is non-enzymatic, reflects the average exposure of hemoglobin to glucose over an extended period. In a classical study, Trivelli et al¹ showed Hemoglobin A_{1c} in diabetic subjects to be elevated 2-3 fold over the levels found in normal individuals. Several investigators have recommended that Hemoglobin A_{1c} serve as an indicator of metabolic control of the diabetic, since Hemoglobin A_{1c} levels approach normal values for diabetics in metabolic control.^{2,3,4} Hemoglobin A_{1c} has been defined operationally as the "fast fraction" hemoglobins (HbA_{1a}, A_{1b}, A_{1c}) that elute first during column chromatography with cation-exchange resins. The non-glycosylated hemoglobin, which consists of the bulk of the hemoglobin has been designated HbA₀. The present procedure utilizes an antigen and antibody reaction to directly determine the concentration of the HbA_{1c}.

REAGENTS

R1	Latex 0.13%, Buffer, stabilizer.
R2	Mouse anti-human HbA _{1c} monoclonal antibody 0.05mg/ml, goat anti-mouse IgG polyclonal antibody 0.08mg/dl, Buffer, stabilizers.
R3 (Hemolysis reagent)	Water and stabilizers
Optional	Ref: 43105 HbA _{1c} CAL. HbA _{1c} (4) Calibrators. Ref: 43106 HbA _{1c} CONTROL. HbA _{1c} (2) Controls.

PRECAUTIONS

All human specimens should be regarded as potentially biohazardous. Therefore, universal precautions should be used in specimen handling (gloves, lab garments, avoid aerosol production, etc.).

PREPARATION

R1, R2 and R3 are ready to use. Mix gently before use.

CALIBRATION

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date. Do not use reagents over the expiration date.

R1 and R2 are stable for at least one month after opening stored at 2-8°C.

Hemoglobin A_{1c} in whole blood collected with EDTA is stable for one week at 2-8°C.⁵

Reagent deterioration: Alterations in the physical appearance of the reagents or values of control materials outside of the manufacturer's acceptable range may be an indication of reagent instability.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Special preparation of the patient is unnecessary. Fasting specimens are not required. No special additives or preservatives other than anticoagulants are required. Collect venous blood with EDTA using aseptic technique.

To determine HbA_{1c}, a hemolysate must be prepared for each sample:

1. Dispense 1 mL Hemolysis Reagent into tubes labeled: Calibrator, Control, Patients, etc. Note: Plastic or glass tubes of appropriate size are acceptable.
2. Place 20 µL of well mixed whole blood into the appropriately labeled lyse reagent tube. Mix.
3. Allow to stand for 5 minutes or until complete lysis is evident. Hemolysates may be stored up to 10 days at 2-8°C.

REFERENCE VALUES

Recommended Values: less than 6% for a non-diabetic, less than 7% for glycemic control of a person with diabetes.

Each laboratory should establish its own expected values. In using Hemoglobin A_{1c} to monitor diabetic patients, results should be interpreted individually. That is, the patient should be monitored against him or herself. There is a 3-4 week time lag before Hemoglobin A_{1c} reflects changes in blood glucose level.

QUALITY CONTROL

HbA_{1c} Control (ref: 43106) is recommended to monitor the performance of manual and automated assay procedures. **Controls require hemolysis pretreatment after being reconstituted.**

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

SPINTECH 240 APPLICATION

Item Name HbA _{1c} -d		CALIBRATION*	
<u>DATA INFORMATION</u>		TYPE Spline	
Units	%	STANDARD	
Decimals	2	#1 Cal.Val.1	#4 Cal.Val.4
<u>ANALYSIS</u>		#2 Cal.Val.2	#5
Type	END	#3 Cal.Val.3	#6
W.Length 1	600	<u>NORMAL RANGE</u>	
W.Length 2		SERUM MALE LOW HIGH	
Method	Turbilatex	FEMALE	
<u>CORR</u>		URINE	
SLOPE	INTER		
1.000 x +	0		
Item Name HbA _{1c} -d		<u>DATA PROCESS</u>	
<u>ASPIRATION</u>		<u>ABSORBANCE LIMIT</u>	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	READ	LOW -3.000
VOLUME**		START END	HIGH 3.000
SAMPLE	6 µL	MAIN 51 52	
REAGENT 1	225 µL	SUB	
REAGENT 2	75 µL	ENDPOINT LIMIT 3	
		LINEAR CHECK (%)	
<u>MONITOR</u>		<u>FACTOR</u>	
Third Mix	<input checked="" type="checkbox"/> OF <input type="checkbox"/> ON	Blank Correction	1.000
R1 Blank	<input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B	<u>PROZONE CHECK</u>	
<u>0 LEVEL POINT</u>		START END LIMIT (%)	
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 2% to linearity limit of 16%.

Precision:

	Intra-assay (n=20)	
Mean (%)	5,95	12,15
SD	0,19	0,18
CV (%)	3,20	1,47

	Inter-assay (n=20)	
	5,97	12,21
	0,14	0,15
	2,31	1,24

Sensitivity: 1% = 0,056 (A)

Accuracy: Results obtained using SPINREACT reagents (y) did not show differences when compared with other commercial reagent (x). The results obtained using 40 samples were the following:

Correlation coefficient (r)² 0,995 and Regression equation: y= 0,989x -0,047

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES AND LIMITATIONS

1. Bilirubin to 50mg/dL, ascorbic acid to 50mg/dL, triglycerides to 2000mg/dL, carbamylated Hb to 7.5mmol/L and acetylated Hb to 5.0mmol/L do not interfere in this assay.
2. This assay should not be used as the unique test for the diagnosis of diabetes mellitus. Other tests must be also considered to establish a correct diagnosis.
3. Patient specimens should always be assayed using a calibration curve.
4. It has been reported that results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, alcoholism, ingest large doses of aspirin.^{6,7,8,9}
5. It has been reported that elevated levels of HbF may lead to underestimation of HA_{1c} and, that uremia does not interfere with HbA_{1c} determination by immunoassay.¹⁰ It has been reported that labile intermediates (Schiff base) are not detected and therefore, do not interfere with HbA_{1c} determination by immunoassay.⁵
6. It has been determined that Hemoglobin variants HbA₂, HbC and HbS do not interfere with this method.
7. Other very rare variants of hemoglobin (e.g. HbE) have not been assessed.

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