		<div style="border: 1px solid black; padding: 2px; display: inline-block;">RF-LATEX</div>
RF-Latex		
Slide agglutination		

Qualitative determination of Rheumatoid Factors (RF) IVD

Store at 2-8°C.

PRINCIPLE OF THE METHOD

The RF-latex is a slide agglutination test for the qualitative and semi-quantitative detection of RF in human serum. Latex particles coated with human gammaglobulin are agglutinated when mixed with samples containing RF.

CLINICAL SIGNIFICANCE

Rheumatoid factors are a group of antibodies directed to determinants in the Fc portion of the immunoglobulin G molecule. Although rheumatoid factors are found in a number of rheumatoid disorders, such as systemic lupus erythematosus (SLE) and Sjögren's syndrome, as well as in nonrheumatic conditions, its central role in clinic lies its utility as an aid in the diagnosis of rheumatoid arthritis (RA). A study of the "American College of Rheumatology" shows that the 80,4% of RA patients were RF positive.

PRECAUTIONS

Control +/-: H317: May cause an allergic skin reaction. Follow the precautionary advice indicated on the SDS and product label. Components from human origin havebeen tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The RF-latex sensitivity is calibrated against the RF International Standard from NIBSC 64/002.

STORAGE AND STABILITY

All the kit components are ready to use, and will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze: frozen reagents could change the functionality of the test. Mix reagents gently before use.

Reagents deterioration: Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Mechanical rotator with adjustable speed at 80-100 r.p.m.
- Vortex mixer.
- Pipettes 50 µL.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at –20°C. Samples with presence of fibrin should be centrifuged before testing. Do not use highly haemolized or lipemic samples.

PROCEDURE

Qualitative method

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Place 50 µL of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
- Mixthe RF-latex reagent rigorously or on a vortex mixerbefore using and add one drop (50 µL) next to the sample to be tested.
- Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
- Place the slide on a mechanical rotator at 80-100 r.p.m. for 2 minutes. False positive results could appear if the test is read later than two minutes.

Semi-quantitative method

- Make serial two fold dilutions of the sample in 9 g/L saline solution.
- Proceed for each dilution as in the qualitative method.

READING AND INTERPRETATION

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator.

The presence of agglutination indicates a RF concentration equal or greater than 8 IU/mL (Note 1). The titer, in the semi-quantitative method, is defined as the highest dilution showing a positive result.

CALCULATIONS

The approximate RF concentrationin the patient sample is calculated as follows:

8 x RF Titer = IU/mL

QUALITY CONTROL

Positive and Negative controls are recommended to monitor the performance of test procedure, as well as a comparativepattern for a better results interpretation. All result different from the negative control result, will be considered as a positive.

REFERENCE VALUES

Up to 8 IU/mL. Each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

- Analytical sensitivity:8 (6-16) IU/mL, under the described assay conditions.
- Prozone effect: No prozone effect was detected up to 1500IU/mL.
- Diagnostic sensitivity: 100%.
- Diagnostic specificity: 100%.

The diagnostic sensitivity and specificity have been obtained using 139 samples compared with the same method of a competitor.

INTERFERENCES

Bilirubin (20 mg/dL), hemoglobin (10 g/L), and lipids (10 g/L), do not interfere. Other substances may interfere⁶.

LIMITATIONS OF PROCEDURE

- Diagnosis should not be solely based on the results of latex method but also should be complemented witha Waaler Rose test along with the clinical examination.

NOTES

- Results obtained with a latex method do not compare with those obtained with Waaler Rose test. Differences in the results between methods do not reflect differences in the ability to detect rheumatoid factors.

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