RPR-CARBON

RPR-carbon
Slide agglutination

RPR-CARBON

RPR-carbon
Agglutination en porta**Qualitative determination of plasma reagins****IVD**

Store at 2 -8°C.

PRINCIPLE OF THE METHOD

The RPR-carbon is a non-treponemal slide agglutination test for the qualitative and semi-quantitative detection of plasma reagins in human serum. Carbon particles coated with a lipid complex are agglutinated when mixed with samples containing reagins of patient affected by syphilis.

CLINICAL SIGNIFICANCE

Reagins are a group of antibodies against some components of the damage tissues from patients infected by *Treponema pallidum*, the agent which causes the syphilis. This microorganism produces some damage to the liver and heart, releasing some tissue fragments. Immunological patient system reacts producing reagins, antibodies against these fragments. The assay is useful to follow the antibiotic therapy answer.

PRECAUTIONS

Control +: H319-Causes serious eye irritation.

Follow the precautionary statements given in MSDS and label of the product.

CALIBRATION

The sensitivity is calibrated against the International Reference WHO (1st International Standard for Syphilitic Human Serum, ref. 05/132).

PREPARATION

RPR-carbon: Homogenize the reagent before use. Open the RPR-carbon vial, place the micropipette to the dispensing vial and draw by suction the required volume of RPR-carbon. Once the test is completed, return the reagent to the original vial and rinse the micropipette and vial with distilled water.

STORAGE AND STABILITY

All the kit components will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze: frozen reagents could change the functionality of the test.

Mix reagents gently before use.

Reagents deterioration: Presence of particles and turbidity.**ADDITIONAL EQUIPMENT-**

- Mechanical rotator with adjustable speed at 80-100 r.p.m.

- Vortex mixer.

- Pipettes 20 and 50 µL.

SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolized or lipemic samples.

PROCEDURE**Qualitative method**

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50 µL of the sample and one drop of Positive control into separate circles on the slide test.
3. Once the reagent has been homogenized as indicated in the preparation, deposit 20 µL of this reagent next to each of the previous drops.
4. Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
5. Place the slide on a mechanical rotator at 80-100 r.p.m. for 8 min (Note 1). False positive results could appear if the test is read later than 8 minutes.

Semi-quantitative method

1. Make serial two fold dilutions of the sample in 9 g/L saline solution.
2. Proceed for each dilution as in the qualitative method.

READING AND INTERPRETATION

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide test from the rotator. Rotate the slide twice by hand before reading.

Interpretation

Agglutination	Reading	Report
Medium or large clumps	R	Reactive
Small clumps	W	Weakly reactive
No clumping or very slight "roughness"	N	Non Reactive

The titer, in the semi-quantitative method, is defined as the highest dilution showing a positive result.

QUALITY CONTROL

Positive and Negative controls are recommended to monitor the performance of procedure, as well as a comparative pattern for a better result interpretation. All result different from the negative control result, will be considered as a positive.

PERFORMANCE CHARACTERISTICS

1. **Analytical sensitivity:** Accurate titer determination of the Reference Material, under the described assay conditions (see calibration).
2. **Prozone effect:** No prozone effect was detected up to titers $\geq 1/128$.
3. **Diagnostic sensitivity:** 100%.
4. **Diagnostic specificity:** 100%.

INTERFERENCES

Bilirubin (20 mg/dL), hemoglobin (10 g/L) and lipids (10 g/L), do not interfere. Rheumatoid factors (300 IU/mL), interfere. Other substances may interfere⁵.

NOTES

1. During the 8 minutes of reaction time do not expose the slide to a source of heat or intense light in order to reduce evaporation. Such evaporation could cause a false agglutination and therefore false positive results.

LIMITATIONS OF THE PROCEDURE

- RPR carbon test is non-specific for syphilis. All Reactive samples should be retested with treponemal methods such as TPHA and FTA-Abs to confirm the results.
- A Non Reactive result by itself does not exclude a diagnosis of syphilis. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.
- False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.

BIBLIOGRAPHY

1. George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
2. Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
3. Sandra Larsen et al. A manual of Test for Syphilis American Public Health Association 1990: 1-192.
4. Joseph Earle Moore et al. Gastrointestinal Haemorrhage 1952; 150(5): 467-473.
5. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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