Lp(a)-CONTROL

Lp(a) Control / Lp(a) Control
Turbidimetry / Turbidimetría**Quantitative determination of Lipoprotein (a) Lp(a) IVD**
Store at 2 - 8°C.**PRODUCT CHARACTERISTICS**

The control is a lyophilised serum from human origin containing Lp(a) used to evaluate the precision and accuracy of Lp(a) determinations in human serum by turbidimetry methods.

REAGENTS

Lp(a) Control	Human Lp(a) serum. With a low concentration level of Lp(a). The Lp(a) concentration is stated on the vial label.
---------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against an Internal Reference Material.

PREPARATION

Reconstitute with 1.0 mL of distilled water. Mix gently, avoiding foam forming. Bring to room temperature for about 15 minutes before use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 2 weeks at 2-8°C

PROCEDURE

To be used in Lp(a) turbidimetric assay.

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 1 x 1 mL

Lp(a)-CONTROL

Lp(a) Control / Lp(a) Control
Turbidimetry / Turbidimetría**Quantitative determination of Lipoprotein (a) Lp(a) IVD**
Store at 2 - 8°C.**PRODUCT CHARACTERISTICS**

The control is a lyophilised serum from human origin containing Lp(a) used to evaluate the precision and accuracy of Lp(a) determinations in human serum by turbidimetry methods.

REAGENTS

Lp(a) Control	Human Lp(a) serum. With a low concentration level of Lp(a). The Lp(a) concentration is stated on the vial label.
---------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against an Internal Reference Material.

PREPARATION

Reconstitute with 1.0 mL of distilled water. Mix gently, avoiding foam forming. Bring to room temperature for about 15 minutes before use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date. After reconstitution, stable for 2 weeks at 2-8°C

PROCEDURE

To be used in turbidimetric Lp(a) assay.

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