Quantitative determination of fecal occult blood IVD

Store at 2-8 °C

PRINCIPLE OF THE METHOD

FOB-Turbilatex is a quantitative turbidimetric test for the measurement of human hemoglobin in human feces.

Latex particles coated with specific anti-human hemoglobin antibody are agglutinated when mixed with samples containing hemoglobin or hemoglobin/haptoglobin complexes. The agglutination causes an absorbance change, dependent upon the hemoglobin contents of the sample that can be quantified by comparison from a calibrator of known hemoglobin concentration.

CLINICAL SIGNIFICANCE

The presence of hemoglobin in stool is associated with lower intestinal tract pathologies like colon polyps, Crohn's disease, ulcerative colitis and colorectal cancer. An early diagnosis is the key for a successful therapy. FOB-Turbilatex allows detecting amounts of hemoglobin that are invisible to the naked eye. The test can be performed with an automatic analyzer using directly a stool sample collection tube. These features make this the ideal test for screening programs.

REAGENTS

Diluent (R1)	Buffered solution. Preservative
Latex (R2)	Latex particles coated with rabbit polyclonal antibody anti-human hemoglobin. Preservative.
FOB-CAL (opcional)	Calibrator. Hemoglobin concentration is stated on each vial label (6 x 1 mL)

PRECAUTIONS

R1 and R2: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary statements given in MSDS and label of the product.

Contains products of animal origin. Handle cautiously as potentially infectious.

CALIBRATION

Use FOB Calibrator Ref. 1107082. The value of the calibrator is traceable to Reference Material WHO International Standard Haemoglobin cyanide NIBSC (98/708). Recalibrate when control results are out of specified tolerances, when using a different lot of reagents and when the instrument is adjusted.

PREPARATION

The reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date. Do not freeze. Reagent deterioration: presence of particles (R1, R2) and turbidity (R1).

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E autoanalyzer.
- Laboratory equipment.
- FOB Calibrator (6 x 1 mL)

SAMPLES

Use only Spinreact's sample collection tube, following the instructions included in the package. Store stool samples in collection tube during maximum 7 days at 2-8°C.

REFERENCE VALUES

Each laboratory should determine its own normal reference value. The sensitivity and specificity of the results with the cutoff of 100ng/mL was 76% and 95%, respectively⁴.

QUALITY CONTROL

Control solutions are recommended to monitor the performance of assay procedures (Control 1 and Control 2). If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

The clinical diagnosis should not be based on a single test result; the complete medical history of the patient should be taken into consideration.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS

Test	FOB	R1	180
No.	**	R2	60
Full Name	FOB	Sample volume	8.5
Standard No.	6	R1 Blank	
Reac. Type	Fixed Time	Mixed Rgt Blank	
Pri. Wavelength	510/505 nm	Linearity Range	0 ng/mL 1100 ng/mL
Sec. Wavelength		Linearity Limit	*
Direction	Increase	Substrate Limit	*
Reac. Time	1 17	Factor	
Incuba. Time		Prozone check	*
Unit	ng/mL	q1	q2
Precision	0.01	q3	q4
		PC	Abs

CALIBRATION (Cal + Rgt Blk)

Rule	Spline	Calib. Level	CONCENTRATION
Sensitivity	1	Cal 0	Stated on the vial label
Replicates	2	Cal 1	
Interval (days)	0	Cal 2	
		Cal 3	
		Cal 4	
		Cal 5	

6 CALIBRATION POINTS

With this parameter the blank must be performed in order to get good results in CALIB screen from main menu. The calibration is stable until 30 days. After this period a new calibration must be performed

PERFORMANCE CHARACTERISTICS

1. Linearity limit: Up to the maximum value of the calibrator (aprox. 1000 ng/mL), under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in sample diluent (obtained from sample collection tube). The linearity limit depends on the sample/reagent ratio and the analyzer used.

2. Quantification limit: Values less than 10 ng/mL give non-reproducible results.

3. Prozone effect: No prozone effect was detected up to 10000 ng/mL.

4. Precision: The reagent has been tested during 19, using two different hemoglobin concentrations in a EP5 (CLSI)-based study.

	CV (%)	
	75 ng/mL	232 ng/mL
Total	5,2	5,0
Within Run	2,4	0,9
Between Run	3,4	2,2
Between Days	3,0	4,5

5. Accuracy: 108 fecal samples were tested, including 14 positives. All results were reported correctly, at the cutoff of 100 ng/mL. Results were compared against a immunocromatographic method from another manufacturer.

The results of the performance characteristics depend on the analyzer used.

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