Quantitative determination of β₂-microglobulin (β₂-m) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The β₂-m Turbilatex is a quantitative turbidimetric test for the measurement of β₂-microglobulin (β₂-m) in human serum, plasma or urine. Latex particles coated with anti-human β₂-m are agglutinated when mixed with samples containing β₂-m. The agglutination causes an absorbance change, dependent upon the β₂-m contents of the patient sample that can be quantified by comparison from a calibrator of known concentration.

CLINICAL SIGNIFICANCE

β₂-m is a protein located on the surface of human lymphocytes and other nucleated cells. Free β₂-m is filtered by the glomerulus and subsequently reabsorbed in the proximal tubular cells. Increased urinary excretion of β₂-m is a sensitive indicator of renal insufficiency. Also, the β₂-m level in serum is a useful marker of other diseases including carcinomas, lymphoid tumors, rheumatoid arthritis and AIDS.

REAGENTS

β₂-m Diluent (R1)	Tris buffer 20 mmol/L, pH 8.2. Preservative.
β₂-m Latex (R2)	Particles coated with goat IgG anti-human β ₂ -m, pH 7.5. Preservative.
β₂-m CAL	Calibrator. β ₂ -m concentration is stated on the vial label.
Optional:	Ref: 1107034 β ₂ -m Control.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

Use β₂-Microglobulin Calibrator Reference 1107032. The sensitivity of the assay and the target value of the calibrator have been standardized against the 1st International ERM-DA470k/IFCC. The calibration in the SPINLAB 180 is stable for 1 month. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION
β₂-m Calibrator:

Serum method: Reconstitute (→) with 1.0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.

Urine method: Dilute reconstituted calibrator 1/6 with NaCl 9 g/L (50 μL calibrator + 250 μL NaCl 9 g/L).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1).

β₂-m Calibrator: Stable for 1 month at 2-8°C or 3 months at -20°C.

ADDITIONAL EQUIPMENT

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C with a 540 nm filter.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Fresh urine. Adjust samples to pH 7-8 by the addition of K₂HPO₄. Stable 2 days at 2-8°C or 2 months at -20°C.

The samples with particles or fibrin should be centrifuged before testing. Do not use hemolyzed or lipemic samples.

PROCEDURE

- Bring the Reagents and the photometer (cuvette holder) to 37°C.
- Assay conditions:
 - Wavelength: 540 nm (530-550).
 - Temperature: 37°C
 - Cuvette light path: 1 cm.
- Adjust the instrument to zero with distilled water.

4. Pipette into a cuvette:

Diluent R1	800 μL
Latex R2	200 μL
Calibrator or sample	10 (serum), 50 (urine) μL

- Mix and read the absorbance immediately (A₁) and after 3 minutes (A₂) of the sample addition.

Spinreact has instruction sheets available for several automatic analyzers. Instructions for many of them are available on request.

CALCULATIONS
Serum:

$$\frac{(A_2 - A_1)_{\text{sample}}}{(A_2 - A_1)_{\text{Calibrator}}} \times \text{Calibrator concentration} = \text{mg/L } \beta_2\text{-m}$$

Urine:

$$\frac{(A_2 - A_1)_{\text{sample}}}{(A_2 - A_1)_{\text{Calibrator}}} \times \frac{\text{Calibrator concentration}}{6} = \text{mg/L } \beta_2\text{-m}$$

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT β₂-m Control (Ref: 1107034). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Serum: from 1.0 to 3.0 mg/L.

Urine: from 0.1 to 0.3 mg/L.

Each laboratory should establish its own reference range

PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 18 mg/L (serum) and 3 mg/L (urine), under the described assay conditions. Samples higher results should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity depends on the sample-reagent ratio, as well as the analyzer used. It will be higher by decreasing sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 0.09 mg/L (serum) and 0.01 mg/L (urine) give non-reproducible results.
- Prozone effect:** No prozone effect was detected up to 60 mg/L (serum) and 20 mg/L (urine).
- Sensitivity:** Δ 0.048 A. mg/L (serum) and Δ 0.228 A. mg/L (urine).
- Precision:** The reagent has been tested for 20 days, using three different β₂m concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 1 mg/L	+/- 3.2 mg/L	+/- 8.5 mg/L
Total	4.0%	3.4%	1.7%
Within Run	2.8%	2.0%	1.2%
Between Run	1.7%	1.5%	1.2%
Between Day	2.2%	2.4%	0.0%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 36 samples of different concentrations of β₂-m were assayed. The correlation coefficient (r) was 0.97 and the regression equation y = 1.709x - 2.627. The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Serum method: bilirubin (20 mg/L), hemoglobin (10 g/L) and lipids (10 g/L), do not interfere. Rheumatoid factors (150 IU/mL), interfere.

Urine method: urea (urine) (50 g/L), uric ac. (20 g/L) and glucose (100 g/L), do not interfere.

Other substances may interfere⁷.

BIBLIOGRAPHY

- Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
- Malaguarnera M et al. Digestive Diseases and Sciences 1997; 42: 762-766.
- Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93.
- Wibell L et al. Nephron 1973; 10: 320-331.
- Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
- Davey P G et al. Clin Chem 1982; 28/6: 1330-1333.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Pres, 1995.

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