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- Pulkki K. et al. Scand J Clin Lab Invest 1986; 46: 605 - 607.
- Müller W. et al. Journal of Immunological Methods 1985; 80: 77 - 90.
- Strang F. et al., Mediators Inflamm., 2014.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AAC Press, 1995.

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Ref. 1107103	Cont.	R1 Diluyente: 1 x 40 mL
		R2 Látex: 1 x 10 mL
		U-CRP CAL: 1 x 2 mL
		CRP CAL high level: 1 x 2 mL

## Quantitative determination of C-Reactive Protein (CRP)

### IVD

Store 2 - 8°C.

### PRINCIPLE OF THE METHOD

CRP-Turbilatex is a quantitative turbidimetric test for the measurement of C-reactive protein (CRP) in human serum or plasma.

Latex particles coated with specific anti-human CRP are agglutinated when mixed with samples containing CRP. The agglutination causes an absorbance change, dependent upon the CRP contents of the patient sample that can be quantified by comparison from a calibrator of known CRP concentration. The reagent in the kit, by using a different calibrator (also included in the kit) and a different application, can be used for a high-sensitive-CRP test (CRP-ultra).

### CLINICAL SIGNIFICANCE

CRP is an acute-phase protein present in normal serum, which increases significantly (up to 300mg/L in 12-24h) after acute conditions such as bacterial/fungal infections and pelvic inflammatory disease. The test is also used for monitoring some inflammatory chronic conditions such as certain types of arthritis and autoimmune diseases. Apart from this main use of the test, the hs-CRP test is used to quantify low levels of CRP (0,5-10mg/L) which, combined with other tests like lipid profile, can be used as a marker for cardiovascular disease.

### REAGENTS

<b>Diluent (R1)</b>	Tris buffer 20 mmol/L, pH 8,2. Preservative.
<b>Latex (R2)</b>	Latex particles coated with goat IgG anti-human CRP, pH 8,2. Preservative.
<b>CRP-CAL (x2)</b>	CRP-Ultra (U-CRP) Calibrator; CRP high level Calibrator C-Reactive protein concentration is stated on the vial label.
<b>Optional</b>	Ref.: 43036 CRP-Ultra (U-CRP) control. Ref.: 1102114 Control serum ASO/CRP/RF Level L. Ref.: 1102115 Control serum ASO/CRP/RF Level H.

### PRECAUTIONS

R1, R2 and CRP-CAL: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950).

Follow the precautionary statements given in MSDS and label of the product.

Components from human origin have been tested and found negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

### CALIBRATION

Use U-CRP Calibrator and/or CRP Calibrator

The sensitivity of the assay and the target value of the calibrator have been standardized against the Reference Material ERM-DA 474/IFCC.

The calibration on the Spintech240 is stable for at least 1 month.

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

### PREPARATION

**U-CRP Calibrator:** Liquid ready to use.

**CRP Calibrator:** Liquid ready to use.

**Calibration Curve:** Prepare the following U-CRP/CRP calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the calibrator by the corresponding factor stated in table below to obtain the concentration of each dilution:

Calibration dilution	0	1/16	1/8	1/4	1/2	1
Calibrator (µL)	--	25	50	100	200	400
NaCl 9 g/L (µL)	400	375	350	300	200	--
Factor	0	0,062	0,125	0,25	0,5	1,0

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

**Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

### ADDITIONAL EQUIPMENT

- Thermostatic bath at 37°C.

- Spectrophotometer or photometer at 37°C with 546 nm and 660nm filters.

### SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

### PROCEDURE

1. Bring the reagents and the photometer (cuvette holder) to 37°C.
2. Assay conditions: Wavelength: 546nm (CRP-ultra); 660 nm (CRP); cuvette light path: 1 cm.
3. Adjust the instrument to zero with distilled water.
4. Pipette into a cuvette:

Diluent R1	800 µL
Latex R2	200 µL
Calibrator or sample	CRP-ultra: 33 µL; CRP: 10 µL

5. Mix and read the absorbance immediately ( $A_1$ ) and 4 minutes after ( $A_2$ ) of the sample addition.

**SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

### CALCULATIONS

Calculate the absorbance difference ( $A_2-A_1$ ) of each point of the calibration curve and plot the values obtained against the CRP concentration of each calibrator dilution. CRP in the sample is calculated by interpolation in the calibration curve.

### QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES

Normal values up to 3 mg/L (CRP-ultra) and 6 mg/L (CRP).

Each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Linearity limit:** CRP-ultra: up to 10 mg/L; CRP: up to the maximum level of the calibrator, approx. 150 mg/L, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again.

**Detection limit:** CRP-ultra: <0,05 mg/L; CRP: <0,2 mg/L give non-reproducible results.

**Prozone effect:** No prozone effect was detected up to 150 mg/L (CRP-ultra); 450 mg/L (CRP).

**Precision:** The reagent has been tested for 20 days, using three different CRP concentrations (mg/L) in an EP5-based study (%CV):

CRP (mg/L)	2,9 (hsCRP)	8,9 (CRP)	37,8 (CRP)
Total	2,4	4,2	4,1
Within Run	0,8	1,3	0,9
Between Run	1,4	3,4	3,6
Between Day	1,8	2,0	1,7

**Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. Sixty-one (CRP-ultra) samples were tested. CRP-ultra:  $r=0,99$ ;  $y=1,15x+0,1$ . Forty-eight (CRP) samples were tested. CRP:  $r=0,999$ ;  $y=0,95x+0,45$ . The results of the performance characteristics depend on the analyzer used.

### INTERFERENCES

Bilirubin (9mg/dL CRP-ultra; 25mg/dL CRP), lipids (5 g/L), rheumatoid factors (300 IU/mL) and hemoglobin (10g/L) do not interfere. Other substances may interfere.

### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

### BIBLIOGRAPHY

1. Pulkki K. et al. Scand J Clin Lab Invest 1986; 46: 605 – 607.
2. Müller W. et al. Journal of Immunological Methods 1985; 80: 77 – 90.
3. Strang F. et al., Mediators Inflamm., 2014.
4. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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