# Stool collection tube

## Stool collection tube for FOB-TURBI assay IVD

For professional in vitro diagnostic use only.  
Store at 2 - 30 °C

### INTENDED PURPOSE

FOB-TURBI STOOL COLLECTION TUBE is a device intended to be used by professional users with the FOB-TURBI medical device to collect, contain and preserve fecal samples to be analyzed with turbidimetric assay analysers.

### TEST PRINCIPLE

Stool collection tube for fecal occult blood testing to be used with Spinreact's FOB-turbidimetric assay. The buffered solution contained in the tube stabilizes the hemoglobin present in the sample during up to 7 days; time enough from sampling to testing. The tube design standardizes the amount of stool collected (17mg) for an easier interpretation of the results.

### REAGENTS

<b>Sample Buffer</b>	1,7 mL of buffered solution, stabilizers, preservatives.	
<b>Related products</b>	FOB-TURBI	Ref.: MI1107185 (200 tests)
		Ref.: MI1107186 (1000 tests)
	FOB Control H & Control L	Ref.: 1107083
	FOB-CAL	Ref.: 1107082

### PRECAUTIONS

Contains sodium azide which can react with lead and copper plumbing, to form potentially explosive azide. When disposing of such reagents, flush with large volumes of water to prevent azide build up.

### STORAGE AND IN-USE STABILITY

The buffered solution contained in the tube is stable until the expiry date printed on the tube, if stored at 2 - 30 °C. The hemoglobin obtained from the sample is stable 7 days if the tube is stored at 2 - 8 °C.

### PROCEDURE

1.

Place toilet paper in the toilet bowl.

2.


Defecate directly onto the toilet paper.

3.

Unscrew the blue cap of the tube.

4.

Dip the stick into 4 different zones before placing it back into the tube.

5.

Screw the blue dip-stick back into the tube.

6.

Identify the tube.

7.

Keep the tube in the refrigerator. Do not freeze. Return to the laboratory as soon as possible.

### NOTES

- Before testing, shake the tube vigorously in a vortex mixer until the fecal sample is completely detached from the sampler. Then, let the tube stand for 1 hour in an upright position with the white cap on top.
- Place the tube in the analyser after releasing the white cap.
- Verify that the tube does not collide with the sampling probe, if necessary, transfer the sample to the instrument sample cup.
- If the tube seal is totally broken discard the sample.

**Any serious incident related to the product shall be reported to Spinreact and to the Competent Authority of the Member State in which the user and/or patient is established.**

### LIMITATIONS

- Avoid sample collection during the menstrual period, and in case of hemorrhoids or blood in urine.
- Blood, serum, or other body fluids containing traces of blood processed with the same analyser may interfere with the assay.
- Liquid samples must be added volumetrically into the collection tub with a pipette (17µL).

### PACKAGING

Ref.: 1107185

Cont.

1 x 100 tubes

### SYMBOLLOGY

"CONTENT". This symbol is not homologated but is placed on the packaging to identify reagents and the volume included in the IVD device kit.

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