FOB-CONTROL

FOB Control
Turbidimetry / Turbidimetria**Quantitative determination of FOB**

IVD

Store at 2 – 8 °C.

PRODUCT CHARACTERISTICS

The control is a liquid solution containing hemoglobin used to evaluate the precision and accuracy of fecal occult blood determinations in human samples by turbidimetric method.

REAGENTS

FOB- Control	Human hemoglobin. Buffered solution. Preservative. Hemoglobin concentration is stated on the vial label.
--------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Calibrator has been standardized against the WHO International Standard Haemiglobincyanide NIBSC 98/708.

PREPARATION

Ready for use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8 °C and contaminations are prevented during their use. Do not use reagents over the expiration date.

PROCEDURE

To be used in Spinreact's FOB-turbidimetric assay.

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Store at 2 – 8 °C.

PRODUCT CHARACTERISTICS

The control is a liquid solution containing hemoglobin used to evaluate the precision and accuracy of fecal occult blood determinations in human samples by turbidimetric method.

REAGENTS

FOB- Control	Human hemoglobin. Buffered solution. Preservative. Hemoglobin concentration is stated on the vial label.
--------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Calibrator has been standardized against the WHO International Standard Haemiglobincyanide NIBSC 98/708.

PREPARATION

Ready for use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

PROCEDURE

To be used in Spinreact's FOB-turbidimetric assay.

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