# Stool collection tube

**ENGLISH**
**Stool collection tube for Calprotectin turbidimetric assay**
**IVD**

For professional in vitro diagnostic use only.  
Store at 2-30 °C

**INTENDED PURPOSE**

CALPRO-TURBI STOOL COLLECTION TUBE is a device intended to be used by professional users with the CALPRO-TURBI medical device to collect, contain and preserve fecal samples to be analyzed with turbidimetric assays analysers.

**TEST PRINCIPLE**

Stool collection tube for calprotectin testing to be used with Spinreact's Calprotectin turbidimetric assay. The buffered solution contained in the tube stabilizes the calprotectin present in the sample during up to 28 days: time enough from sampling to testing. The tube design standardizes the amount of stool collected (10 mg) for an easier interpretation of the results.

**REAGENTS**

<b>Sample Buffer</b>	1.7 mL of buffered solution, stabilizers, preservatives.	
<b>Related products</b>	Calpro-turbi reagent	Ref. MI1107090
	Controls (Low/High)	Ref: 1107093
	Calibrator	Ref: 1107092

**PRECAUTIONS**

H317 - May cause an allergic skin reaction.  
Contains 2-Methylisothiazol-3(2H)-one (Proclin 950).  
Follow the precautionary statements given in MSDS and label of the product.


**STORAGE AND IN-USE STABILITY**

The buffered solution contained in the tube is stable until the expiry date printed on the tube, if stored at 2 - 30 °C. The calprotectin obtained with the sample is stable 28 days if the tube is stored at 2 - 30 °C.

**PROCEDURE**
**1.**

**Place toilet paper in the toilet bowl.**
**2.**

**Defecate directly onto the toilet paper.**
**3.**

**Unscrew the sampler.**
**4.**

**Dip the sampler into 4 different zones before placing it back into the tube.**
**5.**

**Screw the sampler back into the tube.**
**6.**

**Identify the tube.**
**7.**

**Keep the tube at 2 – 30 °C. Do not freeze. Return to the laboratory as soon as possible.**
**NOTES FOR THE LABORATORY**

- Before testing, shake vigorously the tube and leave gently shaking at least 1 hour to help the feces sample dissolve in the buffered solution. If not dissolved, shake vigorously again until the fecal sample is totally detached from the sampler.
- Centrifuge the tube for 10 min at 1000 x g and place the tube in the autoanalyzer after releasing the white cap.
- Verify that the tube does not collide with the sampling probe, if necessary, transfer the sample to the instrument sample cup.
- If the tube seal is totally broken discard the sample.

**Any serious incident related to the product shall be reported to Spinreact and to the Competent Authority of the Member State in which the user and/or patient is established.**

**LIMITATIONS**

- Liquid samples must be added volumetrically into the collection tube with a pipette (10 µL).

**PACKAGING**

Ref.: 1107195

Cont.

1 x 100 tubes

**SYMBOLGY**

"CONTENT". This symbol is not homologated but is placed on the packaging to identify reagents and the volume included in the IVD device kit.

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