HbA_{1c}

Glycated Hemoglobin

Latex turbidimetry

Quantitative determination of glycated hemoglobin (HbA_{1c}) in human blood

IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

This method utilizes the interaction of antigen and antibody to directly determine the HbA_{1c} in whole blood. Total hemoglobin and HbA_{1c} have the same unspecific absorption rate to latex particles. When mouse antihuman HbA_{1c} monoclonal antibody is added (R2), latex-HbA_{1c}-mouse anti human HbA_{1c} antibody complex is formed. Agglutination is formed when goat anti-mouse IgG polyclonal antibody interacts with the monoclonal antibody. The amount of agglutination is proportional to the amount of HbA_{1c} absorbed on to the surface of latex particles. The amount of agglutination is measured as absorbance. The HbA_{1c} value is obtained from a calibration curve.

CLINICAL SIGNIFICANCE

Throughout the circulatory life of the red cell, Hemoglobin A_{1c} is formed continuously by the adduction of glucose to the N-terminal of the hemoglobin beta chain. This process, which is non-enzymatic, reflects the average exposure of hemoglobin to glucose over an extended period. In a classical study, Trivelli et al¹ showed Hemoglobin A_{1c} in diabetic subjects to be elevated 2-3 fold over the levels found in normal individuals. Several investigators have recommended that Hemoglobin A_{1c} serve as an indicator of metabolic control of the diabetic, since Hemoglobin A_{1c} levels approach normal values for diabetics in metabolic control.^{2,3,4}

Hemoglobin A_{1c} has been defined operationally as the "fast fraction" hemoglobins (HbA, A₁B, A₁C) that elute first during column chromatography with cation-exchange resins. The non-glycosylated hemoglobin, which consists of the bulk of the hemoglobin has been designated HbA₀. The present procedure utilizes an antigen and antibody reaction to directly determine the concentration of the HbA_{1c}.

REAGENTS

R1	Latex 0,13%, Buffer, stabilizer.
R2	Mouse anti-human HbA _{1c} monoclonal antibody 0,05mg/ml, goat anti-mouse IgG polyclonal antibody 0,08mg/dl, Buffer, stabilizers.
R3 (Hemolysis reagent)	Water and stabilizers

PRECAUTIONS

All human specimens should be regarded as potentially biohazardous. Therefore, universal precautions should be used in specimen handling (gloves, lab garments, avoid aerosol production, etc.).

PREPARATION

R1, R2 and R3 are ready to use. Mix gently before use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

R1 and R2 are stable for at least one month after opening stored at 2-8°C.

Reagent deterioration: Alterations in the physical appearance of the reagents or values of control materials outside of the manufacturer's acceptable range may be an indication of reagent instability.

SAMPLES

Special preparation of the patient is unnecessary. Fasting specimens are not required. No special additives or preservatives other than anticoagulants are required. Collect venous blood with EDTA using aseptic technique. HbA_{1c} in whole blood collected with EDTA is stable for one week at 2-8°C⁵.

To determine HbA_{1c}, a hemolysate must be prepared for each sample:

- Dispense 0,5 mL Hemolysis Reagent into tubes labeled: Calibrator, Control, Patients, etc. Note: Plastic or glass tubes of appropriate size are acceptable.
- Place 10 µL of well mixed whole blood into the appropriately labeled lyse reagent tube. Mix.
- Allow to stand for 5 minutes or until complete lysis is evident. Hemolysates may be stored up to 10 days at 2-8°C.

QUALITY CONTROL

HbA_{1c} Control (ref: 43106) is recommended to monitor the performance of manual and automated assay procedures. **Controls require hemolysis pretreatment after being reconstituted.**

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹¹

Recommended Values: less than 6% for a non-diabetic, less than 7% for glycemic control of a person with diabetes.

Each laboratory should establish its own expected values. In using Hemoglobin A_{1c} to monitor diabetic patients, results should be interpreted individually. That is, the patient should be monitored against him or herself. There is a 3-4 weeks time lag before Hemoglobin A_{1c} reflects changes in blood glucose level.

NOTES

- In order to avoid contamination, it is recommended to use disposable material.
- Use clean disposable pipette for its dispensation.
- SPINREACT has instruction sheets for several automatic analyzers.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	180 µL
Test	HbA _{1c}	Vol. R2	60 µL
Full Name	Glycated Hemoglobin	Vol. R3	
Standard nº	5	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	5 µL	Decimal	0.01
Vol. Sample Increas.		Unit	%
Vol. Sample Dec		Slope	1
		Inter.	0
REACTION PARAMETERS			
Reac. Type	End Point	Direction	Increase
Pri. Wave.	660	Reagent Blank	0-0
Sec. Wave.		React. Time	76-77
CALIBRATION DILUTION (4 LEVEL CAL + 1 WATER BLANK)			
CALIBRATION TYPE: SPLINE			

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	HbA _{1c} -d	No.	**
Full name	HbA _{1c} -d	Print name	HbA _{1c} -d
Reac. Type	End Point	Direction	Increase
Pri. Wave.	660	Sec. Wave.	
Unit	%	Decimal	0.01
Reagent Blank	0-0	React. Time	81 - 82
Vol. Sample	5 µL	R1	180 µL
Increased		R2	60 µL
Decreased		R3	
Sample blank		R4	
CALIBRATION			
Spline	Water	0	
Calibrator	CAL LEV. 1		
Calibrator	CAL LEV. 2		
Calibrator	CAL LEV. 3		
Calibrator	CAL LEV. 4		

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 2% to linearity limit of 16%.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (%)	5,95	12,15	5,97	12,21
SD	0,19	0,18	0,14	0,15
CV (%)	3,20	1,47	2,31	1,24

Sensitivity: 1% = 0,056 (A)

Accuracy: Results obtained using SPINREACT reagents (y) did not show differences when compared with other commercial reagent (x). The results obtained using 40 samples were the following:

Correlation coefficient (r)² 0,995

Regression equation: y = 0,989x - 0,047

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

- Bilirubin to 50 mg/dL, ascorbic acid to 50 mg/dL, triglycerides to 2000 mg/dL, carbamylated Hb to 7,5 mmol/L and acetylated Hb to 5,0 mmol/L do not interfere in this assay.
- It has been reported that results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, alcoholism, ingest large doses of aspirin.^{6,7,8,9}
- It has been reported that elevated levels of HbF may lead to underestimation of HA_{1c} and, that uremia does not interfere with HbA_{1c} determination by immunoassay.¹⁰ It has been reported that labile intermediates (Schiff base) are not detected and therefore, do not interfere with HbA_{1c} determination by immunoassay.⁵
- It has been determined that Hemoglobin variants HbA₂, HbC and HbS do not interfere with this method.
- Other very rare variants of hemoglobin (e.g. HbE) have not been assessed.

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APPLICATION AU SPIN640

TEST INFORMATION		REAGENT VOLUME	
N°	**	Vol. R1	180 µL
Test	HbA1C	Vol. R2	60 µL
Full Name	Glycated Hemoglobin	Vol. R3	
Standard n°	5	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	5	Decimal	0.01 Slope 1
Vol. Sample Increas.		Unit	% Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	End Point	Direction	Increase
Pri. Wave.	660	Reagent Blank	0-0
Sec. Wave.		React. Time	76-77
CALIBRATION DILUTION (4 LEVEL CAL + 1 WATER BLANK)			
CALIBRATION TYPE: SPLINE			

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