# Total Bilirubin

DPD. Colorimetric

## Quantitative determination of total bilirubin IVD

Store at 2-8°C

### INTENDED USE

Bilirubin T-DPD is an in vitro diagnostic medical device intended to be used by healthcare professionals for the quantitative determination of total bilirubin in serum or plasma.

### PRINCIPLE OF THE METHOD

Bilirubin (both conjugated and unconjugated) couples with the diazo reagent in the presence of a surfactant to form azobilirubin. The intensity of color formed is proportional to the bilirubin concentration in the sample tested. The increase of absorbance at 546 nm is directly proportional to the total bilirubin concentration.

### CLINICAL SIGNIFICANCE

Bilirubin is caused by the degradation of hemoglobin and exists in two forms. Unconjugated bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted. Hyperbilirubinemia is the result of an increase of bilirubin in plasma. Possible causes:

**Total bilirubin:** Increase hemolysis, genetic alteration, neonatal anemia, erythropoiesis alterations and presence of drugs.

**Direct Bilirubin:** cholestasis liver, liver abnormalities and genetic.

Clinical diagnosis should not be made based on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

|          |                         |             |
|----------|-------------------------|-------------|
| R 1      | Surfactants             | <1%         |
|          | Hydrochloric acid (HCl) | 160 mM      |
| R 2      | 2,4-DPD                 | ≥2 mM       |
|          | Hydrochloric acid (HCl) | 120 mM      |
|          | Surfactant              | <1%         |
| Optional | <b>BILIRUBIN CAL</b>    | Ref:1002250 |

### PRECAUTIONS

R1/ R2: H290- Corrosive to metals. H314 - Irritation or skin corrosion.  
R2: contains 2,4-DPD.

Follow the safety advice given in MSDS and product label.

### PREPARATION

The reagents are provided in a ready to use format.

### STORAGE AND STABILITY

The reagents are stable until the expiry date stated on the label when stored at 2-8°C, protected from light and contaminations are prevented during their use. Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.

### ADDITIONAL EQUIPMENT

- Spectrophotometer or analyzer capable of measuring absorbance at 546 nm.
- Cuvettes 1.0 cm light path.
- General laboratory equipment.

### SAMPLES

The use of fresh serum and plasma free of hemolysis is recommended. Protect samples from light.

The samples can be stored in the refrigerator (2-8°C) for 4 days maximum. For longer storage periods they must be kept frozen at -20°C for 2 months. Frozen samples should be totally thawed and brought to room temperature before testing. Avoid repeated freezing and thawing of the samples.

### PROCEDURE

- Assay conditions:  
Wavelength: ..... 546 nm (530-580)  
Cuvette: ..... 1 cm light path  
Temperature: ..... 37°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

|                 | Calibrator blank | Sample blank |
|-----------------|------------------|--------------|
| R 1 (μL)        | 800              | 800          |
| Calibrator (μL) | 40               | -            |
| Sample (μL)     | -                | 40           |

- Mix and incubate for **5 minutes** at 37°C.
- Read the absorbance (A1) of the sample and calibrator.
- Add:

|          | Calibrator | Sample |
|----------|------------|--------|
| R 2 (μL) | 200        | 200    |

- Mix and incubate for **5 minutes** at 37°C.
- Read the absorbance (A2) of the sample and calibrator against the blank.
- Calculate the increase of the absorbance:  $\Delta A = A2 - A1$ .

### CALCULATIONS:

- **With calibrator:**

$$\frac{(\Delta A)_{\text{Sample}}}{(\Delta A)_{\text{Calibrator}}} \times \text{Calibrator conc.} = \text{mg/dL of bilirubin in the sample}$$

- **With Factor:**  $(\Delta A)_{\text{Sample}} \times \text{Factor}^* = \text{mg/dL bilirubin in the sample}$

$$\text{*Factor: } \frac{\text{Calibrator concentration}}{(\Delta A)_{\text{Calibrator}}}$$

**Conversion factor:** mg/dL x 17,1 = μmol/L.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210). If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES

Total bilirubin 0,2-1,2 mg/dL (3,4 – 20,5 μmol/L)

These values are for orientation purpose; each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *quantification limit* of 0,1 mg/dL to *linearity limit* of 30 mg/dL. If the results obtained were greater than the linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

### Precision:

|              | Inter assay (n= 20) |      | Intra assay (n= 20) |      |
|--------------|---------------------|------|---------------------|------|
|              | Mean (mg/dL)        | SD   | CV (%)              |      |
| Mean (mg/dL) | 1,24                | 4.4  | 1,22                | 4.33 |
| SD           | 0,01                | 0,03 | 0,02                | 0,08 |
| CV (%)       | 1.0                 | 0.7  | 2.0                 | 1,9  |

**Sensitivity:** 1 mg/dL = 0,026 Abs. units.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x) on a Spintech 240 analyzer. The results obtained using 60 were:

Correlation coefficient: (r) 0,996

Regression equation:  $y = 0,9836x + 0,1644$

The results of the performance characteristics depend on the analyzer used.

### INTERFERENCES

No interferences were observed for lipemia (Intralipid) up to 2000 mg/dL, hemoglobin up to 1000 mg/dL and ascorbic acid up to 40 mg/L. A list of drugs and other interfering substances with bilirubin has been reported by Young et. al <sup>4,5</sup>.

### NOTES

- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

### BIBLIOGRAPHY

- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2: 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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### BIBLIOGRAPHIE

- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Text book of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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