Fructose-1,6-diphosphate aldolase

GDH - UV

Quantitative determination of Aldolase IVD


Store at 2-8°C

INTENDED USE

For the quantitative *in vitro* determination of Aldolase in serum or plasma.

PRINCIPLE OF THE METHOD

Aldolase converts fructose-1,6-diphosphate (F-1,6-DP) to glyceraldehyde-3-phosphate (GAP) and dihydroxyacetone phosphate (DAP). The addition of triosephosphate isomerase (TIM), glycerolphosphate dehydrogenase (GDH) and NADH converts the dihydroxyacetone phosphate to glycerol-1-phosphate. The rate of the aldolase reaction is measured by the decrease in absorbance at 340 nm as a consequence of the conversion of NADH to NAD⁺.

CLINICAL SIGNIFICANCE

Levels of Aldolase are elevated in cases of muscle or liver damage. Increased serum aldolase is also found with myocardial infarction and is used in the monitoring of patients with muscular dystrophy.

REAGENTS

R 1	Collidine buffer, pH 7.4	51 mmol/L
	Mono-iodoacetate	0.27 mmol/L
	F-1,6-DP	2.7 mmol/L
R 2	NADH	0.23 mmol/L
R3	GDH	≥ 326 mU/mL
	TIM	≥ 4.35 U/mL
	LDH	≥ 616 mU/mL
	Ammonium Sulphate	> 35 %

PREPARATION AND STABILITY

R1. Buffer/Substrate

Reconstitute one vial of Buffer/Substrate R1 with **20 mL** of distilled water. Stable for two weeks at +2 to +8°C.

R2. NADH

Reconstitute one vial of NADH R2 with **1 mL** of distilled water. Stable for 4 weeks at +2 to +8°C.

R3. GDH/TIM/LDH

Contents ready for use. Stable up to the expiry date specified when stored at +2 to +8°C.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 37° C (± 0.1°C)
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Serum, heparinized plasma or EDTA plasma.

PROCEDURE

1. Assay conditions:
 Wavelength: 340nm
 Cuvette: 1 cm light path
 Constant temperature 37°C
2. Adjust the instrument to zero with distilled water.
3. Pipette into a test tube.

	Sample blank	Calibrator	Sample
Sample (mL)	0.2	0.2	0.2
R1 (mL)	-	2.5	2.5
NaCl 9 g/L (mL)	-	-	-
R2 (mL)	-	0.05	0.05
R3 (mL)	-	0.01	0.01

4. Mix sample, R1, R2 and R3, and incubate at 37°C for 5 minutes.
5. Read absorbance A₁ against sample blank. Allow to stand at 37°C for exactly 20 minutes after first reading and then measure absorbance A₂ against blank.
6. If A₁ < 0.95, dilute 1+1 with 0.9% NaCl and reassay. Multiply the result by 2.

CALCULATIONS

To calculate the Aldolase activity use the following formula:

$$\frac{(A_1 - A_{sb}) - (A_2 - A_{cb})}{(A_1 - A_{sb}) - (A_2 - A_{cb})} \times \text{Conc. of calibrator}$$

$$(A_1 - A_{cb}) - (A_2 - A_{cb})$$

Where: A_{sb} = A sample blank, A_{cb} = A calibrator blank.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: Aldolase 2 levels control ref.1002222. Two levels of controls should be assayed at least once a day. Values obtained should fall within the specified range,

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Serum: up to 7.6 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Linearity: The method is linear up to a concentration of 106U/L. If the sample concentration exceeds this value, dilute the sample 1/5 with NaCl 9 g/L and multiply the result by 5.

Sensitivity: The minimum detectable concentration of Aldolase with an acceptable level of precision was determined as 1.73 U/L.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	CV (%)	CV (%)
Mean (U/L)	6.36	19.1	6.36	19.1
SD	0.295	0.854	0.402	1.35
CV (%)	4.64	4.47	6.32	7.09

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results were the following:

Correlation coefficient (r): 0.9917.

Regression equation: y = 0.9815x + 0.183.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Haemolysis interferences with the test.

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