**Quantitative determination of Fibrinogen IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Fibrinogen in presence of an excess of thrombin concentration, changes into Fibrin.

The time for clot formation in dilute plasma is inversely proportional to the fibrinogen concentration in the sample.

The thrombin clotting time fibrinogen assay is based on the method originally described by Clauss.<sup>1</sup> In the presence of high concentrations of thrombin, the time required for clot formation in dilute plasma is inversely proportional to the fibrinogen concentration.

**CLINICAL SIGNIFICANCE**

Fibrinogen (Factor I), protein synthesized by the liver, is the substance used in the blood to form a clot. Its determination is used to evaluate abnormal blood clotting.

Elevated Fibrinogen levels are observed in acute inflammations and in pregnancy; low values are observed in trombotic therapy, in hepatic disease, in the congenital non fibrinogen, in DIC (Disseminated Intravascular Coagulation) and in pancreatitis (low values)<sup>1</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b>	Bovine thrombin	≈ 100 NIH u /mL
<b>R 2</b>	Imidazole Buffer , Sodium azide	
<b>R3</b>	Caolin Solution	
<b>Optional</b>	COAGULATION CAL	REF: 1709101
	CONTROL NORMAL	REF: 1709104
	CONTROL PATHOLOGIC	REF: 1709106

**PRECAUTIONS**

R2: H290- May be corrosive to metals. H302-Harmful if swallowed. H360-May damage fertility or the unborn child. H412-Harmful to aquatic life with long lasting effects

Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

R1: Dissolve (→) the contents with 2.0 mL of distilled water. Cap vial and mix gently to dissolve contents. Stability: 7 days at 2-8°C or 1 month at –20°C, if immediately frozen and stored in the original container. Do not re-freeze.

R2: Mix before use. R3: Ready for use reagent.

**TRACEABILITY**

The Coagulation Calibrator is traceable for Fibrinogen to the WHO standard, 2nd International Standard for Fibrinogen, plasma (98/612).

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Quality control values outside established ranges.
- Product colour variations.

**ADDITIONAL EQUIPMENT**

- Coagulometer or stopwatch and bath at 37°C ± 0.5°C.
- General laboratory equipment<sup>(Note 1)</sup>.

**SAMPLES**

Plasma from venous puncture diluted 1/10 in trisodium citrate solution 3.8% (105 mmol/L).

Mixing immediately the blood with anticoagulant. Avoid foaming the specimen.

Centrifuge the sample at 3000 x g for 10 min and transfer the plasma to siliconized glass or plastic containers.

Turbid, icteric, lipemic or hemolyzed samples may generate erroneous results.

The sample is stable for 4 hours at room temperature (15-25°C) or 28 days if immediately frozen at below 20°C.

**NOTES**

1. All labware must be clean and free of trace amounts of detergents.
2. Always follow instrument manufacturer's instructions; the results must be validated by the test laboratory.

**PROCEDURE**

The reagent can be used by manual procedure, mechanical, photo-optical or other means of end clot detection<sup>(Note 2)</sup>.

1. Dilute the citrated plasma and Control 1/10 with Imidazole buffer:  
50 µL plasma + 450 µL Imidazole buffer.  
The diluted sample must be processed in 1 hour.
2. Prepare the following dilutions of the Calibrator in Imidazole buffer.

Calibrator Dilution	1/40	1/30	1/20	1/10	1/5
IMIDAZOLE BUFFER (mL)	3.9	2.9	1.9	0.9	0.4
COAGULATION CAL (mL)	0.1	0.1	0.1	0.1	0.1
Factor	10/40*	10/30*	10/20*	10/10*	10/5*
Concentration (mg/dL)	0.25* x c	0.33* x c	0.5* x c	1* x c	2* x c

(c = Calibrator value)

3. Add 20 µL of R3 to 0.2 mL of each dilution, and allow to reach 37°C for 4-6 minutes.
4. Add 0.1 mL of R1 and time clot formation. Do not prewarm thrombin R1.

**Calculations**

1. Calculate the mean of duplicate clotting times immediately after reaction. Use all five of the calibrator points to construct a log-log curve that plots fibrinogen concentration (mg/dL) vs. clotting time (s).
2. Draw the straight line of best fit. Examine the curve and, if necessary, omit non-linear points. The final curve must consist of at least three consecutive points. Constructing the curve with only the most linear points will produce the best recovery on control and patient samples.
3. The following curve is **only orientative**. It will change with lot and concentration of the calibrator, a well as, with the instrument used.

Time (s)	Concentration (mg/dL)	Concentration (g/L)
18.1	608	6.08
26.4	304	3.04
49.6	152	1.52
84.7	76	0.76
153	38	0.38

4. Find the clotting time of quality control and patient samples on the curve and read the corresponding fibrinogen value.
5. If clotting times for the 1/10 dilution fall outside the linear curve, prepare 1/5 or 1/20 dilutions as needed. If the sample is diluted 1/5, divide the result from the standard curve by 2; if the sample was diluted 1/20, multiply the curve result by 2 to get the final result.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures:

CONTROL NORMAL REF: 1709104

CONTROL PATHOLOGIC REF: 1709106

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES**

200 - 400 mg/dL<sup>1</sup> (2.00 – 4.00 g/L)

These values are for orientation purpose; each laboratory should establish its own reference range.

**PERFORMANCE CHARACTERISTICS**

**Precision:**

	Inter-assay (n=30)		
	144	294	488
Mean (U/L)			
CV (%)	5.9	3.4	2.9

**Accuracy:** Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents.

**INTERFERENCES AND LIMITATIONS**

Has been observed interferences in samples with fibrinogen degradation.

Acute inflammatory reactions can elevate circulating fibrinogen. Hemolysis can cause clotting factor activation and end point detection interference. High paraprotein levels, and drugs that activate the fibrinolytic system can interfere with fibrinogen assays. A list of drugs and other interfering substances with the determination has been reported<sup>2,3</sup>.

**BIBLIOGRAPHY**

1. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.

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