## Quantitative measurement of D-dimer in human plasma or serum

### IVD

Store 2 - 8°C.

### PRINCIPLE OF THE METHOD

D-dimer in the sample reacts with the anti-human D-dimer mouse monoclonal antibody-coated latex, resulting in agglutination and elevation of turbidity. The resulting turbidity changes are then measured using a spectrophotometer, allowing quantitative measurement of the D-dimer concentration in the sample.

### CLINICAL SIGNIFICANCE

D-dimer is a type of fibrin degradation products (FDP) composed of stable fibrin degraded by plasmin. Stable fibrin is crosslinked by the action of coagulation factor XIII in the blood coagulation and fibrinolysis system. Various types of D-dimer molecules exist in the blood, including YY/DXD, YD/DY, DD/E and DD complex.

Increasing in the level of D-dimer in the blood proves thrombus production, as well as the efficacy of fibrinolysis. Increasing in the level of D-dimer is also known to be associated with various diseases, including malignant tumors, obstetric diseases, vascular lesions, and DIC (disseminated intravascular coagulation syndrome).

### REAGENTS

<b>Diluent (R1)</b>	Tris (hydroxymethyl aminomethane) buffer (pH 8.5)
<b>Latex (R2)</b>	Anti-human D-dimer mouse monoclonal antibody-coated latex (2.8 mg/mL)
<b>Optional</b>	Ref. 1709111 D-Dimer CAL (Calibrator)
<b>Optional</b>	Ref.1709114 D-Dimer Control

### PREPARATION

**Reagents:** Ready to use.

### PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

The reagents contain Proclin 300 as a preservative. Therefore, if the reagents come in contact with skin or clothes, rinse immediately with plenty of water, and consult the doctor if skin irritation develops.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Frozen Latex and Diluent could change the functionality of the test. Protect the reagents from direct sunlight.

**Reagent deterioration:** Presence of particles and turbidity.

### ADDITIONAL EQUIPMENT

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C with a 570 nm filter (560 – 580 nm).

### SAMPLES

Serum or citrated plasma can be used. When using serum, whole blood should be collected with the specific tube for FDP containing thrombin and aprotinin.

Stable 1 day at 2-8°C and 1 month at -80°C.

Do not repeat the freeze-thaw cycle for samples.

### PROCEDURE

1. Bring the reagents and the photometer (cuvette holder) to 37°C.

2. Assay conditions:

Wavelength : 570 nm (560-580 nm)  
 Temperature : 37 °C  
 Cuvette lighth path : 1cm

3. Adjust the instrument to zero with distilled water.

4. Pipette into a cuvette:

R1. Diluent (µL)	400
Calibrator or Sample (µL)	48

5. Incubate for 5 min at 37°C.

6. Pipette:

R2. Latex (µL)	400
----------------	-----

7. Mix and read the absorbance after 30 sec (A1) and 1,4 minutes (A2) of the R2 addition.

**Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

### CALCULATIONS

Calculate the absorbance difference ( $A_2 - A_1$ ) of each point of the calibration curve and plot the values obtained against the D-Dimer concentration of each calibrator dilution. D-Dimer concentration in the sample is calculated by interpolation of its ( $A_2 - A_1$ ) in the calibration curve.

### QUALITY CONTROL

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES

Normal Range

1.0 µg/mL or lower

There may be non-specific reaction or interfering reactions. When a plasma sample is inappropriately collected, false values which are higher than the actual values may be obtained. If measurement results appear unreliable, repeat the measurement (if necessary after dilution) or try another analytical measurement.

### PERFORMANCE CHARACTERISTICS

- Measurement range:** Up to 60 µg/mL, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit and measurement range depends on the sample to reagent / ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Limit detection:** Values less than 0,5 µg/mL give non-reproducible results.
- Sensitivity:**  $\Delta$  0.01-0.05/min per 10 µg/mL of D-Dimer.
- Accuracy:** Results obtained using this reagents (y) were compared to those obtained using a Latex turbidimetric method (x). 120 samples were assayed by both methods. The correlation coefficient (r) was 0.99 and the regression line equation  $y = 0.97x + 0.81$ .

The results of the performance characteristics depend on the analyzer used.

### INTERFERENCES

Hemoglobin (500 mg/dL), bilirubin-C (21 mg/dL), bilirubin-F (17 mg/dL), and rheumatoid factors (500 IU/mL), do not interfere. Other substances may interfere.

### NOTES

Clinical diagnosis should be not made on findings of a single test result, but should integrate both clinical and laboratory data.

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