HbA_{1c} Multicalibrator / Multicalibrador

Latex turbidimetry / Turbidimetria Látex

Quantitative determination of glycosylated hemoglobin (HbA_{1c})

IVD

Store 2 - 8°C.

INTENDED USE

This product is for the purpose of calibrating results in the quantitative determination of human hemoglobin A1c (HbA_{1c}) in blood by automated immunoassay

PRODUCT CHARACTERISTICS

The lyophilized hemoglobin A1c calibrators are a hemolysate prepared from packed human erythrocytes. Stabilizers are added to maintain hemoglobin in the reduced state for the accurate calibration of the hemoglobin A1c procedure.

REAGENTS

HbA _{1c} Calibrator 1	Human blood. Contains preservatives. The HbA _{1c} concentration is indicated on the vial label.
HbA _{1c} Calibrator 2	
HbA _{1c} Calibrator 3	
HbA _{1c} Calibrator 4	

PREPARATION

Reconstitute vials with 0.5mL deionized water. Gently mix for 10 minutes, or until all material has dissolved.

The calibrators should be treated in the same manner as patient specimens regarding the hemolysate procedure.

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

Avoid contact with skin and mucous membranes.

CALIBRATION

The setpoint value of the calibrators was obtained by calibrating with materials referenced to NGSP values using the Spinreact HbA_{1c} Reagent Set.

MATERIAL REQUIRED AND NOT PROVIDED

Hemoglobin A1c Reagent kit, from Spinreact Ref: 43100 or 43101
Hemoglobin A1c Control kit, from Spinreact Ref:43106

STORAGE AND STABILITY

1. Store calibrators at 2-8°C. Stable until expiration date if sealed tightly. PROTECT FROM LIGHT AND HEAT.
2. The reconstituted calibrators should be stored refrigerated (2-8°C) and sealed tightly. The calibrator retains its assigned value for at least 30 days at 2-8°C.

PROCEDURE

To be used in turbidimetric assays.

The lyophilized hemoglobin A1c calibrator kit will produce a calibration curve that will be stable for at least 7 days on most analyzers.

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PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious. Avoid contact with skin and mucous membranes.

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STORAGE AND STABILITY

- Store calibrators at 2-8°C. Stable until expiration date if sealed tightly. PROTECT FROM LIGHT AND HEAT.
- The reconstituted calibrators should be stored refrigerated (2-8°C) and sealed tightly. The calibrator retains its assigned value for at least 30 days at 2-8°C.

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