**Quantitative determination of glycosylated hemoglobin (HbA<sub>1c</sub>)**

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**PREPARATION**

Reconstitute vials with 0.5mL deionized water. Gently mix for 10 minutes. Observe for undissolved material. The reconstituted controls may be dispensed in 0.1 mL aliquots, sealed tightly and frozen at -20°C.

**The lyophilized HbA<sub>1c</sub> controls should be assayed in the same manner as blood specimens including the hemolysate procedure.**

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**STORAGE AND STABILITY**

1. Store controls at 2-8°C. Stable until expiration date if sealed tightly. PROTECT FROM LIGHT AND HEAT.
2. Reconstituted controls retain their assigned values for at least three months if frozen. If not frozen, the reconstituted controls are stable at least one month if stored at 2-8°C and sealed tightly.
3. Do not freeze and thaw more than once.

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**Quantitative determination of glycosylated hemoglobin (HbA<sub>1c</sub>)**
**IVD**

Store 2 - 8°C.

**INTENDED USE**

This product is for the purpose of monitoring accuracy and precision in the quantitative determination of human hemoglobin A1c (HbA<sub>1c</sub>) in blood by automated immunoassay.

**PRODUCT CHARACTERISTICS**

The lyophilized HbA<sub>1c</sub> controls are hemolysates prepared from packed human erythrocytes. The controls provide two levels of HbA<sub>1c</sub>, one level in the normal range and the other level in the elevated range. Stabilizers are added to maintain hemoglobin in the reduced state providing complete control of the HbA<sub>1c</sub> procedure.

**REAGENTS**

<b>HbA<sub>1c</sub> Control Level 1</b>	Human blood. Contains preservatives. The HbA <sub>1c</sub> concentration is indicated on the vial label.
<b>HbA<sub>1c</sub> Control Level 2</b>	

**PREPARATION**

Reconstitute vials with 0,5mL deionized water. Gently mix for 10 minutes. Observe for undissolved material. The reconstituted controls may be dispensed in 0,1 mL aliquots, sealed tightly and frozen at -20°C.

**The lyophilized HbA<sub>1c</sub> controls should be assayed in the same manner as blood specimens including the hemolysate procedure.**
**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

Avoid contact with skin and mucous membranes.

**MATERIAL REQUIRED AND NOT PROVIDED**

Hemoglobin A1c Reagent kit, from Spinreact Ref: 43100 or 43101

Hemoglobin A1c Calibrator kit, from Spinreact Ref:43105

**STORAGE AND STABILITY**

- Store controls at 2-8°C. Stable until expiration date if sealed tightly. PROTECT FROM LIGHT AND HEAT.
- Reconstituted controls retain their assigned values for at least three months if frozen. If not frozen, the reconstituted controls are stable at least one month if stored at 2-8°C and sealed tightly.
- Do not freeze and thaw more than once.

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